# IMMULITE®

## Cortisol

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**SIEMENS**

Siemens Medical Solutions Diagnostics

# IMMULITE®/IMMULITE® 1000 Cortisol

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of cortisol (hydrocortisone, Compound F) in serum, as an aid in the clinical assessment of adrenal status.

Catalog Number: **LKCO1** (100 tests), **LKCO5** (500 tests)

Test Code: **COR** Color: **Dark Pink**

CDC Analyte Identifier Code: 1032

CDC Test System Identifier Code: 10159

CLIA Complexity Category: Moderate

## Summary and Explanation

Cortisol (hydrocortisone, Compound F) is the most abundant circulating steroid and the major glucocorticoid secreted by the adrenal cortex. Physiologically effective in anti-inflammatory activity and blood pressure maintenance, cortisol is also involved in gluconeogenesis, calcium absorption and the secretion of gastric acid and pepsin.

As an indicator of adrenocortical function, measurement of blood cortisol levels is useful in the differential diagnosis of Addison's and Cushing's disease, hypopituitarism, and adrenal hyperplasia and carcinoma.<sup>4</sup>

A variety of stimulation and suppression tests — ACTH Stimulation, ACTH Reserve, Dexamethasone Suppression — can supply supportive information on adrenocortical function.<sup>4,5,6</sup>

Anomalous cortisol concentrations have been shown to exist in patients with acute infections, severe pain, diabetes mellitus or heart failure, and in women either pregnant or on estrogen therapy. In addition, certain virilizing syndromes and iatrogenic conditions raise other naturally occurring steroid levels to high (nonphysiologic) concentrations that may interfere with assaying cortisol directly, unless a highly specific antiserum is used.

## Principle of the Procedure

IMMULITE/IMMULITE 1000 Cortisol is a solid-phase, competitive chemiluminescent enzyme immunoassay.

**Incubation Cycles:** 1 × 30 minutes.

## Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Cortisol has not been tested with all possible variations of tube types.

**Volume Required:** 10 µL serum. (Sample cup must contain at least 100 µL more than the total volume required.)

**Storage:** 7 days at 2–8°C or 3 months at –20°C.<sup>8</sup>

## Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### Cortisol Test Units (LCO1)

Each barcode-labeled unit contains one bead coated with polyclonal rabbit anti-cortisol antibody. Stable at 2–8°C until expiration date.

**LKCO1:** 100 units. **LKCO5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### Cortisol Reagent Wedge (LCO2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to cortisol in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

**LKCO1:** 1 wedge. **LKCO5:** 5 wedges.

### Cortisol Adjustors (LCOL, LCOH)

Two vials (Low and High), 3 mL each, of cortisol in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**LKCO1:** 1 set. **LKCO5:** 2 sets.

## Kit Components Supplied Separately

### Cortisol Sample Diluent (LCOZ)

For the dilution of patient samples. One vial containing 25 mL of processed cortisol-free human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM :** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

**CON6:** Tri-level, multi-constituent control

Also Required

Sample transfer pipets, distilled or deionized water, controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

### Recommended Adjustment Interval:

2 weeks.

**Quality Control Samples:** Use controls or sample pools with at least two levels (low and high) of cortisol.

## Expected Values

The following reference ranges have been reported in the literature for circulating cortisol levels and for the results of stimulation and suppression tests.

Diurnal Variation	<i>a.m.:</i> 5 – 25 µg/dL (138 – 690 nmol/L) <sup>1,8,9</sup> <i>p.m.:</i> Approximately half of <i>a.m.</i> values
ACTH Stimulation	Over twice (usually 3 to 5 times) basal values <sup>4,5</sup>
ACTH Reserve	Below basal (control) values <sup>5,6</sup>
Dexamethasone Suppression	Below basal (control) values for the screening, low-dose and high-dose tests <sup>4,5</sup>

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitations

Circulating cortisol results may be falsely elevated in samples obtained from patients being treated with prednisolone or prednisone (converted to prednisolone *in vivo*). Caution must therefore be exercised with cortisol determinations for patients undergoing therapy with these and structurally related synthetic corticosteroids.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data representative of the assay's performance. Results are expressed in µg/dL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

### Conversion Factor:

µg/dL × 27.59 → nmol/L

**Calibration Range:** 1 – 50 µg/dL  
(28 – 1,380 nmol/L).

**Analytical Sensitivity:** 0.2 µg/dL  
(5.5 nmol/L).

**Precision:** Samples were processed in quadruplicate over the course of 5 days, 4 runs per day, for a total of 20 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three cortisol solutions (109, 215 and 349 µg/dL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for cortisol. (See "Specificity" table.)

Although some steroids exhibit crossreactivity, their normal physiological concentrations are low compared to cortisol, hence they will not significantly interfere in the IMMULITE Cortisol procedure. There is, however, an observed crossreactivity of approximately 49% with prednisolone which must be considered for patients undergoing this therapy. Since prednisone is converted to prednisolone *in vivo*, caution must be exercised with cortisol determinations for patients undergoing either therapy.

**Bilirubin:** Samples spiked with 100 and 200 mg/L of conjugated and unconjugated bilirubin were analyzed. Bilirubin may interfere with the assay, causing elevation of values. (See "Bilirubin" table.)

**Hemolysis:** Presence of packed red blood cells in concentrations up to 30 µL/mL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

**Method Comparison:** The assay was compared to Coat-A-Count Cortisol on 83 patient samples. (Concentration range: approximately 3.3 to 41 µg/dL. See graph.) By linear regression:

$$(IML) = 0.95 (CAC) + 0.32 \mu\text{g/dL} \\ r = 0.988$$

Means:

14.6 µg/dL (IMMULITE)  
15.0 µg/dL (CAC)

## References

- 1) Foster L, Dunn R. Single-antibody technique for radioimmunoassay of cortisol in unextracted serum or plasma. Clin Chem 1974;20:365.
- 2) Farmer R, Pierce C. Plasma cortisol determinations: radioimmunoassay and competitive protein binding compared. Clin Chem 1974; 20:411.
- 3) Ruder H, et al. A radioimmunoassay for cortisol in plasma and urine. J Clin Endo Metab 1972;35:219.
- 4) Rothfeld B, ed. Plasma Cortisol. In: Nuclear medicine in vitro. 1974:120.
- 5) Murphy B, et al. Clinical studies utilizing a new method for the serial determination of plasma corticoids. J Canad Med Assoc 1964;90:775.
- 6) Sparks R. Measurement of serum 11-deoxycortisol and cortisol after metyrapone. Ann Intern Med 1971;75:717.
- 7) Kowalski A, Paul W. A simple

extraction procedure for the determination of free (unconjugated) cortisol in urine by radioimmunoassay. Clin Chem 1979;25:1152. 8) Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: W.B. Saunders, 1994. 9) Wilson JD, Foster DW, editors. Williams textbook of endocrinology. 7th ed. Philadelphia: W.B. Saunders, 1985. 10) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

## Technical Assistance

In the United States, contact Siemens Medical Solutions Diagnostics Technical Services department.  
Tel: 800.372.1782 or 973.927.2828  
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Siemens Medical Solutions Diagnostics is certified to ISO 13485:2003.

## Tables and Graphs

### Precision ( $\mu\text{g/dL}$ )

	Mean <sup>3</sup>	Within-Run <sup>1</sup>			Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV	
1	2.7	0.24	8.8%	0.27	10%	
2	4.4	0.33	7.5%	0.37	8.4%	
3	11.3	0.88	7.8%	0.87	7.7%	
4	19.0	1.1	5.8%	1.2	6.3%	
5	29.1	1.8	6.2%	2.2	7.6%	
6	35.8	2.4	6.7%	2.6	7.3%	

### Linearity ( $\mu\text{g/dL}$ )


	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	16 in 16 <sup>5</sup>	31.9	—	—
	8 in 16	13.3	16.0	83%
	4 in 16	7.5	8.0	94%
	2 in 16	3.8	4.0	95%
	1 in 16	1.9	2.0	95%
2	16 in 16	34.5	—	—
	8 in 16	15.8	17.3	91%
	4 in 16	8.6	8.6	100%
	2 in 16	4.7	4.3	109%
	1 in 16	2.1	2.2	95%
3	16 in 16	45.3	—	—
	8 in 16	20.9	22.7	92%
	4 in 16	11.2	11.3	99%
	2 in 16	5.7	5.7	100%
	1 in 16	3.1	2.8	111%

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EC REP DPC Biemann GmbH