IMMULITE®

Vitamin B12

For use on the IMMULITE®
and IMMULITE® 1000 systems

DPC®

IMMULITE®/IMMULITE® 1000 Vitamin B12

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of vitamin B12 in serum or heparinized plasma, as an aid in clinical diagnosis and treatment of anemia.

Catalog Number: **LKVB1** (100 tests),
LKVB5 (500 tests)

Test Code: **VB** Color: **Violet**

Summary and Explanation

Vitamin B12 (cobalamin) and folate are nutrients essential to hematopoiesis.⁶ Megaloblastic anemia is almost always due to lack of one of these two vitamins.⁶ Vitamin B12 deficiency can also result in severe neurological impairment.^{1,2,7}

Circulating levels of vitamin B12 are usually a good index to tissue stores. That is, vitamin B12 levels as measured in serum or plasma by an optimized assay system are typically low in vitamin B12 deficiency, and normal or elevated otherwise.¹ Exceptions to this rule can occur in those relatively uncommon situations where levels of vitamin B12 transport proteins are abnormal.² Thus, low circulating vitamin B12 levels can occur in the absence of vitamin B12 deficiency where the level of transcobalamin I (a physiologically inactive transport protein) is low.⁸

Conversely, vitamin B12 deficiency can occur in the presence of normal or even elevated plasma vitamin B12 levels, where transcobalamin II levels are low or where levels of inactive vitamin B12 transport proteins are high, as in chronic myelogenous leukemia.^{5,8,11} (Circulating folate levels are usually normal or elevated in vitamin B12 deficiency, but red cell folate levels are frequently low in this condition.¹⁰)

Vitamin B12 deficiency occurs only rarely as a result of dietary lack of this vitamin.^{3,6} More commonly, it results from impaired absorption, as in partial or total gastrectomy, or in pernicious anemia, a condition characterized by absence or

near absence of intrinsic factor.⁶ Since roughly two thirds of all patients with pernicious anemia have blocking antibodies to intrinsic factor (IFbAb), while IFbAb are only very rarely encountered in other situations, IFbAb determinations represent a useful follow-up test for the differential diagnosis of vitamin B12 deficiency.^{4,8,10} (Circulating intrinsic factor antibodies are present in more than half of all pernicious anemia patients. Increased transport protein levels can occur, for example, in chronic myelogenous leukemia.)

Common causes of high vitamin B12 levels include liver disease, myeloproliferative disease (with chronic myelogenous leukemia as a special case) and the use of multivitamin supplements.

Principle of the Procedure

IMMULITE/IMMULITE 1000 Vitamin B12 is a solid-phase, competitive chemiluminescent enzyme immunoassay.

IMMULITE/IMMULITE 1000 Vitamin B12 involves a preliminary heat denaturation step. Vitamin B12 in the patient sample is released from carrier proteins by incubation at 100°C in the presence of dithiothreitol and potassium cyanide to inactivate vitamin B12-binding proteins, even at extreme levels, as well as antibodies to intrinsic factor.

After the heat denaturation step, the treated patient sample and hog intrinsic factor are simultaneously introduced into an IMMULITE/IMMULITE 1000 Test Unit containing a polystyrene bead coated with a B12 analog, and incubated for approximately 30 minutes at 37°C with intermittent agitation. During this incubation, vitamin B12 in the treated sample competes with the B12 analog on the solid phase for a limited number of vitamin B12 binding sites on the purified intrinsic factor. (Endogenous vitamin B12 analogs do not interfere, because the binder is free of R-protein.) Alkaline phosphatase-labeled anti-hog intrinsic factor is introduced, and the Test Unit is incubated for another 30-minute cycle. The unbound enzyme conjugate is removed by a centrifugal wash.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

Note: If sample is to be analyzed for both vitamin B12 and folic acid, the patient must be in a fasting state.¹⁸

Although neither bilirubin nor hemolysis has any clinically significant effects on the vitamin B12 results, the presence of hemoglobin in the sample will increase folic acid values. Additionally, hemolysis may indicate mistreatment of the specimen before receipt by the laboratory, hence the results should be interpreted with caution.

Because EDTA has a significant effect on results, it should not be used as an anticoagulant.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE/IMMULITE 1000 Vitamin B12 has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 200 µL of sample — serum or plasma — is required for the sample pretreatment step. A single determination uses 100 µL of the *treated* sample: the sample cup should contain at least 250 µL more of the treated sample than the total volume required for the number of vitamin B12 determinations to be performed.

Storage: If not assayed within 8 hours, store at –20°C: stable for 6-8 weeks.¹⁴ Avoid excessive exposure to direct light.

Dilution of High Samples: The IMMULITE Vitamin B12 procedure has a

calibration range of 150 to 1,200 pg/mL (111 to 885 pmol/L). All samples expected to have vitamin B12 levels greater than this value should be diluted with Vitamin B12 Sample Diluent before assay.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

The Borate-KCN Buffer Solution contains cyanide. Extreme care must be taken to avoid all bodily contact with this reagent.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chloramphenicol, at concentrations less than 0.1 g/dL has been added as a preservative. Chloramphenicol is known to cause cancer; this disclosure is required by the state of California.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Vitamin B12 Test Units (LVB1)

Each barcode-labeled unit contains one bead coated with vitamin B12 analog. Stable at 2–8°C until expiration date.

LKVB1: 100 units. **LKVB5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

Vitamin B12 Reagent Wedges (LVBA, LVBB)

With barcodes **LVBA**: one wedge (7.5 mL) containing vitamin B12 binding protein (purified hog intrinsic factor), with preservative. **LVBB**: one wedge (7.5 mL) containing alkaline phosphatase (bovine calf intestine) conjugated to murine monoclonal anti-hog intrinsic factor antibody in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKVB1: 1 set. **LKVB5**: 5 sets.

Vitamin B12 Adjustors (LVBL, LVBH)

Two vials (Low and High) of lyophilized vitamin B12 in a human protein-based matrix, with preservative. Reconstitute each vial with **4.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

LKVB1: 1 set. **LKVB5**: 2 sets.

Note: The Adjustors must be subjected to the same heat denaturation step as patient samples. (See Sample Pretreatment section.)

Borate-KCN Buffer Solution (LBCN)

125 mL of Borate-KCN Buffer Solution, with preservative. Stable at 2–8°C for 30 days after opening.

LKVB1: 1 vial. **LKVB5**: 5 vials.

Caution! Contains cyanide. Extreme care must be taken to avoid all bodily contact with this reagent.

Dithiothreitol Solution (LDTT)

3 mL of a dithiothreitol solution. Stable at 2–8°C for 30 days after opening.

LKVB1: 1 vial. **LKVB5**: 5 vials.

Kit Components Supplied Separately**Vitamin B12 Sample Diluent (LVBZ)**

For the manual dilution of patient samples. One vial containing 25 mL of human serum protein-based diluent, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

CON6: Tri-level, multi-constituent control

Also Required

Sample transfer pipets, distilled or deionized water, controls.

Plain 12x75 mm polypropylene tubes, with loose-fitting caps — available from DPC.

Micropipets: 20, 200, and 1,000 µL. For the 1 mL reagent addition, a reliable repeating dispenser is suitable.

Plastic container with lid (for preparing the Working Solution).

Foil — or other means for protecting the hemolysates from light.

Covered boiling waterbath (100°C). Do *not* use an oven or a heat block.

Ambient temperature waterbath — for cooling the tubes after the boiling step.

4.0 mL volumetric pipet (± 0.02 mL) for reconstitution of Adjustors.

Preparation of Working Solution

Note: Samples being tested for both vitamin B12 and folic acid must be pretreated with the Working Solution described in the package insert for IMMULITE Folic Acid, catalog number LKFO.

The amount of each component depends on the number of tests to be performed. The volumes required, in microliters per test, are tabulated below. Be sure to multiply these volumes by a number slightly greater than the number of tests to be run. (The components are supplied in volumes sufficient for making up the Working Solution in approximately a 20% excess.)

	µL/test
Borate-KCN Buffer Solution	1,000
Dithiothreitol Solution	20

Important: The Working Solution should be prepared on a daily basis. If not used immediately, it should be refrigerated at

2–8°C for a period of not more than 24 hours.

Sample Pretreatment

1. Pipet **200 µL** of each Adjustor, control, or patient sample — serum or plasma — into the tubes prepared.

For patient samples expected to have vitamin B12 levels above the assay's calibration range, dilute the serum or plasma with Vitamin B12 Sample Diluent.

2. Add **1,000 µL** of the Working Solution to all tubes. Vortex.
3. Loosely cap all tubes and place them in a covered, boiling waterbath (100°C) for **15 – 20 minutes**.
4. Remove the tubes from the boiling waterbath, and cool them in an ambient waterbath for **5 minutes**.
5. Pipet at least **350 µL** of the treated sample to an IMMULITE/IMMULITE 1000 Sample Cup.

Treated samples, both serum and plasma, are stable at room temperature 15 – 28°C or refrigerated at 2 – 8°C for 1 hour prior to assay. Note that commercial controls may show variable stability after treatment.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Note that both Reagent Wedges A and B must be loaded on the carousel to run this assay.

Each Sample cup holder can be followed by up to four test units.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of vitamin B12.

Expected Values

Serum samples from 245 apparently healthy male and female blood donors who had not donated blood within the last 12 months, and 147 serum samples from apparently healthy volunteers were processed with the IMMULITE Vitamin B12 assays, yielding a reference range of 174 to 878 pg/mL (128 to 648 pmol/L)

Consider these limits as guidelines only. Each laboratory should establish its own reference ranges.

Limitations

EDTA plasma is not recommended for use.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1998;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without anticoagulants, gel barriers or clot-promoting additives.)

Conversion Factor:
 $\text{pg/mL} \times 0.7378 \rightarrow \text{pmol/L}$

Calibration Range: 150 – 1,200 pg/mL (111 – 885 pmol/L).

Analytical Sensitivity: 125 pg/mL (92 pmol/L).

Intraassay Precision (Within-Run): Statistics were calculated for samples from the results of 15 replicates in a single run. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run): Statistics were calculated for samples assayed in 10 different runs. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with four vitamin B12 solutions (2,160, 3,640, 5,460 and 7,640 pg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: As one test of the specificity of the purified intrinsic factor binder, mixtures of vitamin B12 analogs (containing mostly cobinamide) were assayed for vitamin B12 by the IMMULITE Vitamin B12 assay. The results showed no crossreactivity with cobinamide, even at concentrations as high as 1,800,000 pg/mL and 50,000,000 pg/mL.

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 381 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: Samples ($n=22$) collected into plain and heparinized vacutainer tubes were assayed:

(Heparin) = 0.95 (Serum) + 8.5 pg/mL
 $r = 0.987$

Means:
474 pg/mL (Serum)
460 pg/mL (Heparin)

Because EDTA would have a significant effect on results, it should not be used as an anticoagulant.

Protein Effect: Two patient serum samples, spiked with varying amounts of human serum albumin (HSA) and normal

human IgG, were assayed. (See "Protein Effect" table for results.)

Method Comparison 1: The assay was compared to DPC's Dualcount[®] Solid Phase No Boil (SPNB) assay on 83 patient serum samples. (Concentration range: approximately 155 to 1,143 pg/mL. See "Method Comparison 1" graph.) By linear regression:

(IML) = 1.00 (SPNB) + 13.6 pg/mL
 $r = 0.979$

Means:
426 pg/mL (IMMULITE)
414 pg/mL (SPNB)

Method Comparison 2: The IMMULITE Vitamin B12 procedure was also compared to a commercially available radioimmunoassay (Kit A) on 97 patient serum samples. (Concentration range: approximately 153 to 1,218 pg/mL. See "Method Comparison 2" graph.) By linear regression:

(IML) = 0.91 (Kit A) + 18.6 pg/mL
 $r = 0.963$

Means:
406 pg/mL (IMMULITE)
424 pg/mL (Kit A)

References

- 1) Allen RH. Clinical role and current status of serum cobalamin (Vitamin B12) assays. *Ligand Quarterly* 1981 Fall;4(3):37-44, 67.
- 2) Carmel R. Megaloblastic anemia: vitamin B12 and folate. In: Fairbanks VF, editor. *Current hematology*. New York: John Wiley, 1983: 243-80.
- 3) Colman N. The radioisotopic investigation of anemia. *Ligand Quarterly* 1981 Fall;4(3):24-30.
- 4) Dawson DW. Diagnosis of vitamin B12 deficiency. *Br Med J* 1984;289:938-9.
- 5) Grsbeck R. Biochemistry and clinical chemistry of vitamin B12 transport and the related diseases. *Clin Biochem* 1984;17:99-107.
- 6) Herbert V. The nutritional anemias. *Hosp Pract* 1980 Mar;15(3):65-83, 87-9.
- 7) Herbert V. Biology of disease: megaloblastic anemias. *Lab Invest* 1985;52:3-19.
- 8) Kapadia CR, Donaldson RM. Disorders of cobalamin (vitamin B12) absorption and transport. *Ann Rev Med* 1985;36:93-110.
- 9) Kolhouse JF, et al. Cobalamin analogues are present in human plasma and can mask cobalamin deficiency because current radioisotope dilution assays are not specific for true cobalamin. *New Engl J Med* 1978;299:785-92.
- 10) Lindenbaum J. Status of laboratory testing in the diagnosis of megaloblastic anemia. *Blood* 1983;61:624-7.
- 11) Parry TE. The diagnosis of megaloblastic anemia. *Clin Lab Haematol* 1980;2:89-109.
- 12) Mollin DL, Ross GIM. Serum vitamin B12 concentrations in leukemia and in some other

hematological conditions. Brit J Haemat 1995;1:155-72. 13) Gilbert HS, et al. Serum vitamin B12 content and unsaturated vitamin B12 binding capacity in myeloproliferative disease. Ann Intern Med 1969;71:719-29. 14) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:636. 15) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 16) Food and Drug Administration (FDA), Letters from the Bureau of Medical Devices (BMD) to manufacturers of Vitamin B12 radioassay kits. April 1980, January 1981. 17) National Committee for Clinical Laboratory Standards (NCCLS). Guidelines for evaluating a B12 (cobalamin) assay. (Villanova, PA: NCCLS proposed standard PLSA-12, March 1979). Revised as "NCCLS recommendations with respect to vitamin B12 assays" — presented to the Food and Drug Administration (FDA) panel, 9 July 1979. 18) Chanarin I. The megaloblastic anemias. 2nd ed. Oxford: Blackwell Scientific, 1979. 19) Zucker RM, Podell E, Allen RH. Multiple problems with current no-boil assays for serum cobalamin. Ligand Quarterly 1981 Fall;4(3):52-8. See further: letters by Ziering S, et al, and by Zucker RM, et al, *ibid*, 60-3, and by Allen RH, *ibid*, 1982 Spring;5(1):48-9. 20) Kubasik NP, Ricotta M, Sine HE. Commercially-supplied binders for plasma cobalamin (vitamin B12) analysis — "purified" intrinsic factor, "cobinamide"-blocked R-protein binder, and non-purified intrinsic factor-R-protein binder — compared to microbiological assay. Clin Chem 1980;26:598-600.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Intraassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³
1	159	18	11.3%
2	204	21	10.3%
3	401	27	6.7%
4	736	39	5.3%
5	1,308	77	5.9%

Interassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³
1	310	23	7.4%
2	660	40	6.1%
3	1,192	202	17%

Linearity (pg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	558	—	—
	4 in 8	229	279	82%
	2 in 8	140	140	100%
	1 in 8	77	70	110%
2	8 in 8	1,052	—	—
	4 in 8	561	526	107%
	2 in 8	270	263	103%
	1 in 8	147	132	111%


Recovery (pg/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	250	—	—
	A	323	346	93%
	B	408	420	97%
	C	550	511	108%
	D	711	620	115%
2	—	527	—	—
	A	627	609	103%
	B	723	683	106%
	C	815	774	105%
	D	880	883	100%
3	—	816	—	—
	A	830	883	94%
	B	903	957	94%
	C	1,000	1,048	95%
	D	1,122	1,157	97%

Protein Effect (pg/mL)


	Added ¹	Vitamin B12, pg/mL ²
1	—	541
	30 g/L HAS	494
	50 g/L HAS	493
	20 g/L IgG	533
	30 g/L IgG	537
	60 g/L IgG	412
	50 g/L HSA plus	517
	30 g/L IgG	517
2	—	854
	30 g/L HAS	800
	50 g/L HAS	745
	20 g/L IgG	775
	30 g/L IgG	783
	60 g/L IgG	642
	50 g/L HSA plus	741
	30 g/L IgG	741

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
DPC Biermann GmbH