IMMULITE®

Calcitonin

For use on the IMMULITE®
and IMMULITE® 1000 systems

DPC®

IMMULITE® / IMMULITE 1000® Calcitonin

English

Intended Use: For in vitro diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of calcitonin (thyrocalcitonin) in human serum or heparinized plasma, as an aid in the diagnosis and treatment of diseases involving the thyroid and parathyroid glands, including carcinoma and hyperparathyroidism.

Catalog Number: **LKCL1** (100 tests)

Test Code: **CAL** Color: **Dark Pink**

CDC Analyte Identifier Code: 1041

CDC Test System Identifier Code: 10159

CLIA Complexity Category: Moderate

Summary and Explanation

Human calcitonin is a single-chain peptide hormone produced primarily in the thyroid gland, where it is secreted by the parafollicular "C" cells rather than by the thyroxine-secreting cells.¹ The calcitonin chain is 32 amino acids in length, and has a molecular mass of 3,418 daltons.¹ In pathological conditions, the hormone exhibits a certain degree of heterogeneity. Its half-life in plasma is approximately 10 minutes. Metabolized predominantly by the kidney, calcitonin is also degraded to some extent by a plasma factor, making it important to freeze samples soon after collection.^{1,3}

The physiological significance of calcitonin is uncertain, but it does have a modest calcium- and phosphorus-lowering effect. Its secretion is normally modulated by rising and falling calcium levels, with an increase in circulating calcium prompting an increase in the calcitonin level.^{1,18}

A fasting basal range of 3–19 pg/mL for healthy adults, based on radioimmunoassay of extracted plasma, has been reported.² Elevated levels (greater than 100 pg/mL) are encountered in a variety of pathological conditions, notably medullary thyroid carcinoma, a tumor of the calcitonin-secreting cells of the thyroid. The peptide is also frequently

elevated in leukemic and myeloproliferative disorders.⁹

In addition, it may be produced ectopically by tumors, e.g. those of the lung or breast.^{10,11,12,15,19} Its measurement, however, does not appear to be of value in the diagnosis of such disorders.

Elevations may also be seen in connection with hyperparathyroidism, hypergastrinemia, renal failure and chronic inflammatory disease.^{13,18}

Principle of the Procedure

IMMULITE/IMMULITE 1000 Calcitonin is a solid-phase, enzyme-labeled, two-site chemiluminescent immunometric assay.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

A fasting morning specimen is usually obtained to establish baseline values. If the samples are being collected as part of a provocative test, the time of collection should be noted relative to calcium or pentagastrin infusion.

EDTA has an effect on the assay and is not recommended as a sample type.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Calcitonin has not been tested with all possible variations of tube types. Consult the section on

Alternate Sample Types for details on tubes that have been tested.

Volume Required: 75 µL serum or heparinized plasma. (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: Stable 15 days at –20°C or longer at –70°C.²³ Freeze immediately.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Calcitonin Test Units (LCL1)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-calcitonin. Stable at 2–8°C until expiration date.

LKCL1: 100 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

Calcitonin Reagent Wedge (LCL2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-calcitonin, in a protein buffer. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKCL1: 1 wedge.

Calcitonin Adjustors (LCLL, LCLH)

Two vials (Low and High) of lyophilized calcitonin, in a protein/buffer matrix. Reconstitute each vial with **3.0 mL** distilled or deionized water. Aliquot and freeze. Stable at –20°C for 3 months after reconstitution.

LKCL1: 1 set.

Kit Components Supplied Separately

Calcitonin Sample Diluent (LCLZ)

For the manual dilution of patient samples. 25 mL of calcitonin-free protein/buffer matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

LCLCM: Bi-level, protein/buffer-based control module.

Also Required

Sample transfer pipets, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval: 4 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of calcitonin.

Expected Values

Based on its relationship to DPC's IMMULITE 2000 Calcitonin (see Method Comparison 2), the assay can be expected to have essentially the same reference ranges. Serum samples from 120 male and 90 female subjects in apparent good health were analyzed by the IMMULITE 2000 Calcitonin procedure.

Males: Upper 95% range: 8.4 pg/mL
(2.46 pmol/L)
Absolute range: ND – 18.2 pg/mL
(ND – 5.33 pmol/L)

Females: Upper 95% range: 5.0 pg/mL
(1.46 pmol/L)
Absolute range: ND – 11.5 pg/mL
(ND – 3.36 pmol/L)

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Multiple forms of immunoreactive calcitonin are present in the circulation of healthy subjects, patients with medullary thyroid carcinoma, and patients with other malignancies. These different forms of calcitonin (polymers and/or glycosylated forms) can have molecular weights varying from 3.4 to 70 kDa. Therefore, samples may exhibit varied response in an immunoassay and non-linear dilution.²⁴⁻²⁸

Some healthy individuals may be encountered who have elevated immunoreactive calcitonin levels.²

Because calcitonin is highly labile, special attention must be paid to proper specimen collection and preservation.^{1,3,7}

Where several samples are collected from a single patient, e.g. in pentagastrin stimulation testing or serial monitoring, all samples should be assayed by the same procedure.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have

been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

pg/mL × 0.2926 → pmol/L

Reportable Range: 2 – 2,000 pg/mL
(0.6 – 585 pmol/L) (WHO 2nd IRP 89/620).

High-dose Hook Effect:

None up to 25,000 pg/mL.

Analytical Sensitivity: 2 pg/mL
(0.6 pmol/L).

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three calcitonin solutions (600, 6,000 and 18,000 pg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for calcitonin. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 550 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 29 volunteers into

plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of calcitonin, to obtain values throughout the reportable range of the assay, and then assayed by the IMMULITE Calcitonin procedure.

(Heparin) = 1.02 (Serum) + 21 pg/mL
r = 0.998

(EDTA) = 0.72 (Serum) + 10 pg/mL
r = 0.997

(SST) = 0.91 (Plain Serum) + 19 pg/mL
r = 0.996

Means:

756 pg/mL (Plain Serum)

790 pg/mL (Heparin)

555 pg/mL (EDTA)

709 pg/mL (SST)

Method Comparison 1: The assay was compared to a commercially available (Kit A) on 53 samples. (Concentration range: approximately up to 1,200 pg/mL. See graph.) By linear regression:

(IML) = 0.71 (Kit A) + 4.2 pg/mL
r = 0.990

Means:

155 pg/mL (IMMULITE)

212 pg/mL (Kit A)

Method Comparison 2: The assay was compared to DPC's IMMULITE 2000 Calcitonin on 150 samples. (Concentration range: approximately up to 1,400 pg/mL. See graph.) By linear regression:

(IML) = 0.91 (IML 2000) + 13.2 pg/mL
r = 0.991

Means:

253 pg/mL (IMMULITE)

265 pg/mL (IMMULITE 2000)

References

1) Austin L, Heath H. Calcitonin; physiology and pathophysiology. *New Engl J Med* 1981;304:269-78. 2) Body J-J, Heath H. 'Nonspecific' increases in plasma immunoreactive calcitonin in healthy individuals. *Clin Chem* 1984;30:511-4. 3) Catherwood BD, Deftos LJ. General principles, problems, and interpretation in the radioimmunoassay of calcitonin and parathyroid hormone. In: Bikle DD, editor. *Assay of calcium-regulating Hormones*. New York: Springer-Verlag, 1983: 151-68. 4) Charpin C, et al. Calcitonin, somatostatin and ACTH immunoreactive cells in a case of familial bilateral thyroid medullary carcinoma. *Cancer* 1982;50:1806-14. 5) Deftos LJ, et al. Influence of age and sex on plasma calcitonin in human beings. *New Engl J Med*

1980;302:1351-3. 6) Heath H, Sizemore GW. Plasma calcitonin in normal man. *J Clin Invest* 1977;60:1135-40. 7) Heath H, Sizemore GW. Radioimmunoassay for calcitonin. *Clin Chem* 1982;28:1219-26. 8) Hillyard C, et al. Familial medullary thyroid carcinoma. *Lancet* 1978;i:1009-11. 9) Hillyard C, et al. Immunoreactive calcitonin in leukaemia. *Br Med J* 1979;2:1392-93. 10) Hillyard C, Stevenson J. Calcitonin as a tumour marker. *Ligand Rev* 1980 (Spring);2(1):31-3. 11) Milhaud G. Ectopic secretion of calcitonin. In: Pecile A, editor. *Calcitonin* 1980. Amsterdam: Excerpta Medica, 1981: 154-69. 12) Mulder H, Hackeng WHL. Ectopic secretion of calcitonin. *Acta Med Scand* 1978;204:253-6. 13) Mulder H, et al. Hypercalcitoninaemia in patients with chronic inflammatory disease. *Neth J Med* 1980;23:129-31. 14) Parthemore J, Deftos LJ. Calcitonin secretion in normal human subjects. *J Clin Endocrinol Metab* 1978;47:184-8. 15) Schwartz KE, et al. Calcitonin in nonthyroidal cancer. *J Clin Endocrinol Metab* 1979;49:438-44. 16) Sizemore GW, Go VLW. Stimulation tests for diagnosis of medullary thyroid carcinoma. *Mayo Clin Proc* 1975;50:53-6. 17) Stephanas AV, et al. Medullary thyroid carcinoma; importance of serial serum calcitonin measurements. *Cancer* 1979;43:825-37. 18) Toffaletti J. Calcitonin. *Clinical Chemistry News* 1984 (July);10(7):20-1, 27. 19) Wallach S, et al. Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung. *J Clin Endocrinol Metab* 1981;53:602-6. 20) Wells SA. Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland. *Ann Surg* 1978;188:139-41. 21) Wells SA, et al. The importance of early diagnosis in patients with hereditary medullary thyroid carcinoma. *Ann Surg* 1982;195:595-9. 22) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard*. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 23) Tietz NW, editor. *Clinical guide to laboratory tests*. 3rd ed. Philadelphia: W.B. Saunders 1995: 98. 24) Goltzman D, Tischler AS. Characterization of the immunochemical forms of calcitonin released by a medullary thyroid carcinoma in tissue culture. *J Clin Invest* 1978 Feb; 61:449-58. 25) Becker KL, Snider RH, Silva OL, Moore CF. Calcitonin heterogeneity in lung cancer and medullary thyroid cancer. *Acta Endocrinol* 1978 Sep; 89(1):89-99. 26) Zajac JD, Martin TJ, Hudson P, Niall H, Jacobs JW. Biosynthesis of calcitonin by human lung cancer cells. *Endocrinology* 1985 Feb; 116(2):749-55. 27) Becker KL, Gazdar AF, Carney DN, Snider RH, Moore CF, Silva OL. Calcitonin secretion by continuous cultures of small cell carcinoma of the lung: incidence and immunoheterogeneity studies. *Cancer Lett* 1983 Mar; 18(2):179-85. 28) Bertagna XY, Nicholson WE, Pettengill OS, Sorenson GD, Mount CD, Orth DN. Ectopic production of high molecular weight calcitonin and corticotropin by human small cell carcinoma cells in tissue culture:

evidence for separate precursors. J Clin Endocrinol Metab 1978 Dec; 47(6):1390-3.

Technical Assistance

In the United States, Contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (pg/mL)

| | Mean ³ | Within-Run ¹ | | Total ² | |
|---|-------------------|-------------------------|-----------------|--------------------|-------|
| | | SD ⁴ | CV ⁵ | SD | CV |
| 1 | 22 | 0.86 | 3.9% | 2.5 | 11.4% |
| 2 | 44 | 1.3 | 3.0% | 3.3 | 7.5% |
| 3 | 145 | 3.3 | 2.3% | 10.9 | 7.5% |
| 4 | 621 | 13.8 | 2.2% | 65.6 | 10.6% |
| 5 | 1,029 | 24.5 | 2.4% | 87.3 | 8.5% |
| 6 | 1,207 | 27.4 | 2.3% | 115 | 9.5% |


Linearity (pg/mL)

| | Dilution ¹ | Observed ² | Expected ³ | %O/E ⁴ |
|---|-----------------------|-----------------------|-----------------------|-------------------|
| 1 | 8 in 8 ⁵ | 85.3 | — | — |
| | 4 in 8 | 46.6 | 42.7 | 109% |
| | 2 in 8 | 24.7 | 21.3 | 116% |
| | 1 in 8 | 12.8 | 10.7 | 120% |
| 2 | 8 in 8 | 190 | — | — |
| | 4 in 8 | 97.2 | 95.0 | 102% |
| | 2 in 8 | 46.8 | 47.5 | 99% |
| | 1 in 8 | 23.9 | 23.8 | 100% |
| 3 | 8 in 8 | 396 | — | — |
| | 4 in 8 | 188 | 198 | 95% |
| | 2 in 8 | 96 | 99 | 97% |
| | 1 in 8 | 49 | 50 | 98% |
| 4 | 8 in 8 | 918 | — | — |
| | 4 in 8 | 427 | 459 | 93% |
| | 2 in 8 | 205 | 230 | 89% |
| | 1 in 8 | 103 | 115 | 90% |
| 5 | 8 in 8 | 1,642 | — | — |
| | 4 in 8 | 772 | 821 | 94% |
| | 2 in 8 | 361 | 411 | 88% |
| | 1 in 8 | 174 | 205 | 85% |

Recovery (pg/mL)


| | Solution ¹ | Observed ² | Expected ³ | %O/E ⁴ |
|---|-----------------------|-----------------------|-----------------------|-------------------|
| 1 | — | NA | — | — |
| | A | 30.8 | 30 | 103% |
| | B | 302 | 300 | 101% |
| | C | 988 | 900 | 110% |
| 2 | — | 16.9 | — | — |
| | A | 42.8 | 46.1 | 93% |
| | B | 324 | 316 | 103% |
| | C | 960 | 916 | 105% |
| 3 | — | 228 | — | — |
| | A | 256 | 247 | 104% |
| | B | 520 | 517 | 101% |
| | C | 1,256 | 1,117 | 112% |
| 4 | — | 515 | — | — |
| | A | 535 | 519 | 103% |
| | B | 850 | 789 | 108% |
| | C | 1,608 | 1,389 | 116% |
| 5 | — | 885 | — | — |
| | A | 897 | 871 | 103% |
| | B | 1,185 | 1,141 | 104% |
| | C | 2,000 | 1,741 | 115% |

Method Comparison 1

(IML) = 0.71 (Kit A) + 4.2 pg/mL
r = 0.990

Method Comparison 2

(IML) = 0.91 (IML 2000) + 13.2 pg/mL
r = 0.991

Specificity

| Compound ¹ | pg/mL Added ² | Apparent % Cross-pg/mL ³ | reactivity ⁴ |
|-----------------------|--------------------------|-------------------------------------|-------------------------|
| ACTH | 100,000 | ND | ND |
| C-Peptide | 40,000,000 | ND | ND |
| Insulin | 33,000,000 | ND | ND |
| HGH | 910,000 | ND | ND |
| Prolactin | 1,000,000 | ND | ND |
| PTH | 150,000 | ND | ND |
| TSH | 140,000 | ND | ND |
| Salmon Calcitonin | 50,000 | ND | ND |
| Pre-Procalcitonin | 2,000,000 | ND | ND |
| Gastrin (G-17) | 2,000,000 | ND | ND |
| Human Big Gastrin | 1,000,000 | ND | ND |
| CCK | 10,000,000 | ND | ND |

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH