IMMULITE®

Troponin I

For use on the IMMULITE®
and IMMULITE® 1000 systems

SIEMENS

Siemens Medical Solutions Diagnostics

IMMULITE®/IMMULITE® 1000 Troponin I

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of troponin I in serum, heparinized or EDTA plasma, as an aid in the diagnosis of acute myocardial infarction (AMI).

Catalog Number: **LKTI1** (100 tests), **LKTI5** (500 tests)

Test Code: **TPI** Color: **Light Gray**

CDC Test System Identifier Code: 10159
CLIA Complexity Category: Moderate

Summary and Explanation

Acute myocardial infarction (AMI) is usually diagnosed on the basis of chest pain, electrocardiographic changes, and elevations of markers of myocardial injury. The MB isoenzyme of creatine kinase (CK-MB) has been the preferred marker for two decades.² A study by Wu, et al. found excellent clinical sensitivity of the CK-MB assay between 6 and 24 hours following onset of AMI, with decreased sensitivity beginning in the 24 to 48-hour interval.¹³ However, CK-MB levels can also increase in patients with acute or chronic muscle disease who lack apparent cardiac injury. In the same study, myoglobin, a muscle protein considered to be an early marker of AMI, became elevated within 6 hours of onset, achieved peak clinical sensitivity in the interval 6 to 12 hours following onset, and no longer offered diagnostic value by hour 24.¹³ Myoglobin, although valuable for the early information it provides, also lacks specificity for cardiac injury. A marker specific for myocardial injury is therefore highly desirable.

Cummins, et al.^{6,7} reported the release of cardiac troponin I (cTnI) in AMI. Many studies have focused on cTnI as a candidate marker with acceptable sensitivity and specificity for AMI and other cardiac diseases.

Troponin, a molecule that binds to the thin filament (actin) of striated muscle fibers, acts with intracellular calcium to control

the interaction of the thin filament with the thick filament (myosin), thus regulating muscle contraction. Troponin consists of three subunits: T, which connects the troponin complex and tropomyosin (another cardiac muscle regulatory protein); I, which prevents muscle contraction in the absence of calcium; and C, which binds calcium.⁸ Cardiac troponin I (MW 22.5 kDa) and the two skeletal muscle isoforms of troponin I have considerable amino acid sequence homology, but cTnI contains an additional N-terminal sequence¹¹ and is highly specific for myocardium.¹

Clinical studies report several desirable features of cTnI as a marker of myocardial injury. cTnI rises early in AMI patients and attains levels that are clearly separated from baseline values, so that by 7 hours following onset, the cTnI test detects 95 percent of patients in whom AMI will be confirmed.⁹ Plasma values of cTnI remain elevated for several days, providing a long window for detection of cardiac injury.^{3,13} cTnI has also demonstrated value for predicting mortality risk in unstable angina and in non-Q wave myocardial infarction.⁵

cTnI has demonstrated equivalent diagnostic accuracy for AMI when compared with lactate dehydrogenase type 1 and CK-MB,^{3,10} and may clarify diagnosis in contexts where elevated CK-MB cannot be attributed with certainty to cardiac injury alone.³ These include surgery,⁴ traumatic injury, renal failure, seizures, and skeletal muscle myopathies.¹

In addition, a study on patients undergoing coronary artery bypass grafting (CABG) showed cTnI to be a sensitive marker for perioperative myocardial infarction (PMI); the peak concentration and time of peak both served as diagnostic criteria.¹²

Principle of the Procedure

IMMULITE/IMMULITE 1000 Troponin I is a solid-phase, enzyme-labeled chemiluminescent immunometric assay. The solid phase (bead) is coated with monoclonal murine anti-troponin I antibody. The liquid phase consists of alkaline phosphatase (bovine calf

intestine) conjugated to polyclonal goat anti-troponin I antibody.

The patient sample and the reagent are incubated together with the coated bead for 30 minutes. During this time, troponin I in the sample forms the antibody sandwich complex with monoclonal murine anti-troponin I antibody on the bead and the enzyme conjugated polyclonal goat anti-troponin I antibody in the reagent. Unbound patient sample and enzyme conjugate are then removed by centrifugal washes. Finally, chemiluminescent substrate is added to the test unit containing the bead and the signal is generated in proportion to the bound enzyme.

Incubation Cycles: 1 × 30 minutes.
Time to First Result: 42 minutes.

Specimen Collection

EDTA plasma yields lower assay results. When monitoring, the baseline sample and all subsequent samples should be collected using the same type of tubes, either plain, or EDTA, or heparinized. Do *not* interchange the type of tubes. Please see the Alternate Sample Type section for further information.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Icteric or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Troponin I has not been tested with all possible variations of tube types. Consult the

section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 µL serum or plasma. (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: 5 days at 2–8°C¹⁷ or 1 month at –20°C.¹⁸

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Troponin I Test Units (LT11)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-troponin I antibody. Stable at 2–8°C until expiration date.

LKT11: 100 units. **LKT15:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

Troponin I Reagent Wedge (LT12)

With barcode. 7.5 mL of alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-troponin I antibody in buffer, with

preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKTI1: 1 wedge. **LKTI5:** 5 wedges.

Troponin I Adjustors (LTIL, LTIH)

Two vials (Low and High) of lyophilized troponin I in a nonhuman serum matrix, with preservative Reconstitute each vial with **3.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Let stand for 30 minutes and **use immediately**, or aliquot and freeze. Stable at –20°C for 2 months (aliquotted) after reconstitution.

LKTI1: 1 set. **LKTI5:** 2 sets.

Kit Components Supplied Separately

Troponin I Sample Diluent (LTIZ)

For the manual dilution of patient samples. One vial (25 mL) containing a troponin I-free nonhuman serum matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

CCCM: A bi-level, nonhuman serum-based Cardiac Marker Control Module, containing troponin I.

Also Required

Sample transfer pipets, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval: 2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of troponin I.

Expected Values

Two hundred and fifty-five serum samples from healthy laboratory volunteers and from hospitalized patients that had been shown to be negative for troponin I by another immunometric method were analyzed using the IMMULITE Troponin I assay. The median value for these samples was nondetectable; 98% of the values were below 1.0 ng/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

The assay yields lower values when used with EDTA plasma. (See the "Alternate Sample Type" section.)

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: Up to 180 ng/mL

Analytical Sensitivity: 0.1 ng/mL

High-dose Hook Effect:

None up to 11,000 ng/mL.

Precision: Using a study design comparable to NCCLS EP5-A, a representative spectrum of samples was processed in quadruplicate, one run per day (alternating between morning and afternoon shifts), for a total of at least 20 runs per sample. Results were analyzed using one-way analysis of variance. (See "Precision" table.)

Specificity: The antibodies are highly specific for troponin I. (See "Specificity" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three troponin I solutions (160, 850 and 2,600 ng/mL) were assayed. (See "Recovery" table for representative data.)

Bilirubin: May cause a depression of values. (See "Bilirubin" table.)

Hemolysis: Presence of hemoglobin in concentrations up to 570 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 5,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: Samples were collected into plain, heparinized and EDTA vacutainer tubes.

(Heparin) = 0.94 (Serum) + 0.13 ng/mL
r = 0.99

(EDTA) = 0.71 (Serum) + 1.6 ng/mL
r = 0.97

Method Comparison A:¹⁴ The assay was compared to a commercially available immunofluorescent assay for troponin I (Kit A) on 144 patient serum samples (Concentration range: nondetectable to approximately 26 ng/mL. See "Method Comparison A" graph.) By linear regression:

(IML) = 1.58 (Kit A) – 0.03 ng/mL
r = 0.969

Means:
2.13 ng/mL (IMMULITE)
1.37 ng/mL (Kit A)

Method Comparison B:¹⁴ The assay was compared to a commercially available

immunofluorescent assay for troponin I (Kit B) on 144 patient serum samples (Concentration range: nondetectable to approximately 26 ng/mL. See "Method Comparison B" graph.) By linear regression:

(IML) = 1.04 (Kit B) + 0.22 ng/mL
r = 0.982

Means:
2.13 ng/mL (IMMULITE)
1.83 ng/mL (Kit B)

References

- 1) Adams J, Bodor G, Davila-Roman V, Delmez JA, Apple FS, Ladenson JH, et al. Cardiac troponin I: a marker with high specificity for cardiac injury. *Circulation* 1993;88:101-6. 2) Adams JE III, Abendschein DR, Jaffe AS. Biochemical markers of myocardial injury. Is MB creatine kinase the choice for the 1990's? *Circulation* 1993;88:750-63. 3) Adams JE III, Schechtman KB, Landt Y, Ladenson JH, Jaffe AS. Comparable detection of acute myocardial infarction by creatine kinase MB isoenzyme and cardiac troponin I. *Clin Chem* 1994;40:1291-5. 4) Adams JE, Sicard G, Allan BT, Bridwell KH, Lenke LG, Davila-Roman VG, et al. More accurate diagnosis of perioperative myocardial infarction with measurement of cardiac troponin I. *N Engl J Med* 1994;330:670-4. 5) Antman EM, Tanasijevic MJ, Thompson B, Schactman M, McCabe CH, Cannon CP, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. *N Engl J Med* 1996;335:1342-9. 6) Cummins B, Cummins P. Cardiac specific troponin I release in canine experimental myocardial infarction: development of a sensitive enzyme-linked immunoassay. *J Mol Cell Cardiol* 1987;19:999-1010. 7) Cummins B, Auckland ML, Cummins P. Cardiac-specific troponin I radioimmunoassay in the diagnosis of acute myocardial infarction. *Am Heart J* 1987;113:1333-44. 8) Darnell J, Lodish H, Baltimore D. *Molecular cell biology*. New York: Scientific American Books, 1986: 827-8. 9) Fonarow GC. *UCLA Clinical Practice Guideline* 1996 July;2(3) on Internet at <http://www.cost-quality.com/2,3art.html>. 10) Jaffe AS, Landt Y, Parvin CA, Abendschein DR, Geltman EM, Ladenson JH. Comparative sensitivity of cardiac troponin I and lactate dehydrogenase isoenzymes for diagnosing acute myocardial infarction. *Clin Chem* 1996;42:1770-6. 11) Larue C, Defacque-Lacquement H, Calzolari C, Nguyen DL, Pau B. New monoclonal antibodies as probes for human cardiac troponin I: epitopic analysis with synthetic peptides. *Molec Immunol* 1992;29:271-8. 12) Mair J, Larue C, Mair P, Balogh D, Calzolari C, Puschendorf B. Use of cardiac troponin I to diagnose perioperative myocardial infarction in coronary artery bypass grafting. *Clin Chem* 1994;40:2066-70. 13) Wu AHB, Fen YJ, Contois JH, Pervaiz S.

Comparison of myoglobin, creatine kinase-MB, and cardiac troponin I for diagnosis of acute myocardial infarction. *Ann Clin Lab Sci* 1996;26:291-300. 14) Peetz D, Hafner G. Institute of Clinical and Laboratory Medicine, University of Mainz. Unpublished results. 15) Wong SS. Strategic utilization of cardiac markers for the diagnosis of acute myocardial infarction. *Ann Clin Lab Sci* 1996;26(4):301-12. 16) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 17) Data on file. 18) Tietz NW, editor. *Clinical guide to laboratory tests*. 3rd ed. Philadelphia: W.B. Saunders, 1995:614.

Technical Assistance

In the United States, contact Siemens Medical Solutions Diagnostics Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Siemens Medical Solutions Diagnostics is certified to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	0.80	0.046	5.8%	0.067	8.4%
2	1.56	0.071	4.6%	0.117	7.5%
3	3.9	0.133	3.4%	0.27	6.9%
4	8.0	0.217	2.7%	0.61	7.6%
5	18.9	0.67	3.5%	1.16	6.1%
6	34	1.01	3.0%	2.08	6.1%
7	86	3.1	3.6%	6.5	7.6%

Recovery (ng/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	0	—	—
	A	8.3	8.0	104%
	B	42.7	42.5	100%
	C	128	130	98%
2	—	19.2	—	—
	A	26.6	26.2	102%
	B	63.0	61.0	103%
	C	175	148	118%
3	—	36.6	—	—
	A	41.6	42.8	98%
	B	79.0	77.8	102%
	C	180	165	109%

Specificity

Compound ¹	ng/mL Added ²	Apparent ng/mL ³	% Cross-reactivity ⁴
Human Skeletal Troponin I	200	0.34	0.17%
Human Cardiac Troponin C	200	ND	ND
Human Cardiac Troponin T	45	0.32	0.71%

ND: not detectable⁵


Linearity


	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	5.33	—	—
	4 in 8	2.84	2.67	106%
	2 in 8	1.30	1.33	98%
	1 in 8	0.77	0.67	115%
2	8 in 8	46.2	—	—
	4 in 8	21.0	23.1	91%
	2 in 8	10.4	11.6	90%
	1 in 8	5.75	5.78	99%
3	8 in 8	95.9	—	—
	4 in 8	48.5	48.0	101%
	2 in 8	24.7	24.0	103%
	1 in 8	12.6	12.0	105%
4	8 in 8	120	—	—
	4 in 8	51.0	60.0	85%
	2 in 8	24.3	30.0	81%
	1 in 8	12.4	15.0	83%

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EC REP DPC Biemann GmbH