# IMMULITE®

## IL-6

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**DPC®**

# IMMULITE®/IMMULITE® 1000 IL-6

## English

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of interleukin 6 (IL-6) in serum, EDTA or heparinized plasma, as an aid in the study of inflammatory diseases.

Catalog Number: **LK6PZ** (50 tests), **LK6P1** (100 tests)

Test Code: **I6** Color: **Dark Green**

### Summary and Explanation

Interleukin-6 (IL-6) is a mediator of the immune system, which has a wide variety of biological actions. It is also known as B cell stimulatory factor (BCSF), B cell stimulation factor 2 (BSF-2), hybridoma growth factor (HGF), hepatocyte stimulatory factor (HSF), cytolytic differentiation factor for T lymphocytes (CDF) and many others.

The IL-6 cDNA is coding for a polypeptide consisting of 212 amino acids. This protein is cleaved to a mature protein of 184 amino acids. Due to different degrees of glycosylation — at positions 73 and 172, respectively — and phosphorylation, IL-6 has a molecular weight ranging between 21.5 and 28 kDa.

Many different cells are capable of IL-6 synthesis including monocytes/macrophages, fibroblasts, endothelial cells, keratinocytes, mast cells, T cells and many tumor cell lines. *In vivo* and *in vitro*, IL-6 acts as a differentiation factor for B cells and an activation factor for T cells. In the presence of IL-2, it drives the differentiation of T cells to cytotoxic T cells and induces the proliferation of thymocytes. After stimulation with IL-4, an IL-6 stimulus is necessary for B cells to differentiate into Ig-secreting plasma cells. IL-6 is a potent growth factor of different human myelomas and is active in concentrations less than 10 pg/mL. IL-3 and IL-6 show *in vitro* synergistic effects in the differentiation of hematopoietic progenitor cells.

Elevated IL-6 serum or plasma levels may occur in different diseases including sepsis, autoimmune diseases, lymphomas, AIDS, alcoholic liver disease and in patients with infections, or transplant rejection.

### Principle of the Procedure

IMMULITE/IMMULITE 1000 IL-6 is a solid-phase, enzyme-labeled, chemiluminescent sequential immunometric assay.

**Incubation Cycles:** 2 × 30 minutes.

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on tube materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 IL-6 has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 100 µL serum, EDTA or heparinized plasma. (Sample cup must contain at least 250 µL more than the total volume required.)

**Storage:** 1 day at 2–8°C or 6 months at –20°C.

### Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### IL-6 Test Units (L6P1)

Each barcode-labeled unit contains one bead coated with a monoclonal murine anti-IL-6 antibody. Stable at 2–8°C until expiration date.

**LK6PZ:** 50 units. **LK6P1:** 100 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### IL-6 Reagent Wedges (L6PA, L6PB)

With barcodes. **L6PA:** 7.5 mL of a protein/buffer matrix, with preservative. **L6PB:** 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal sheep anti-IL-6 antibody in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

**LK6PZ:** 1 set. **LK6P1:** 2 sets.

### IL-6 Adjustors (L6PL, L6PH)

Two vials (Low and High) of lyophilized IL-6 in a protein/buffer matrix, with preservative. Reconstitute each vial with **3.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 3 days after

reconstitution, or for 6 months (aliquotted) at –20°C.

**LK6PZ:** 1 set. **LK6P1:** 1 set.

## Kit Components Supplied Separately

### IL-6 Sample Diluent (L6PZ)

For the manual dilution of patient samples. One vial containing 25 mL of IL-6-free protein/buffer matrix, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM :** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

**LILCM:** A bi-level, human serum-based IMMULITE Cytokine Control Module containing IL-6.

Also Required

Sample transfer pipets, distilled or deionized water, controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

The sample cup must contain at least 250 µL more than the total volume required. Each sample cup holder can be followed by up to four test units.

Note that both Reagent Wedges A and B must be loaded on the carousel to run this assay.

**Recommended Adjustment Interval:**

2 weeks.

**Quality Control Samples:** Use controls or sample pools with at least two levels (low and high) of IL-6.

## Expected Values

Based on its relationship to DPC's IMMULITE 2000 IL-6 (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

A reference range study for IMMULITE 2000 IL-6 was performed on serum samples collected from 60 healthy laboratory volunteers into plain glass serum tubes. The non-parametric lower 95% range was from nondetectable to 3.4 pg/mL. The absolute range was nondetectable to 5.9 pg/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Calibration Range:** Up to 1,000 pg/mL  
NIBSC 1st IS 89/548 (IU/mL):  
IU/mL = 0.1 × pg/mL.

**Analytical Sensitivity:** 2 pg/mL

**High-dose Hook Effect:**  
None up to 60,000 pg/mL.

**Precision:** Samples were repeatedly assayed in duplicate over the course of several days for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three IL-6 solutions (1,000, 4,000, and 8,000 pg/mL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for IL-6. (See "Specificity" table.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Presence of hemoglobin in concentrations up to 550 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from 12 volunteers into plain glass and plastic serum tubes, plastic heparinized and EDTA tubes, and plastic gel barrier tubes (SST<sup>®</sup>). All tubes were from Becton Dickinson. Equal volumes of the matched samples were spiked with various concentrations of IL-6, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE/IMMULITE 1000 IL-6 procedure.

(EDTA Plastic) = 0.99 (Serum Glass) + 4.8 pg/mL

(Heparin Plastic) = 0.98 (Serum Glass) + 4.6 pg/mL

(SST Plastic) = 0.71 (Serum Glass) + 8.4 pg/mL

(Serum Plastic) = 0.76 (Serum Glass) + 6.6 pg/mL

Means:

167 pg/mL (Serum Glass)

170 pg/mL (EDTA Plastic)

168 pg/mL (Heparin Plastic)

127 pg/mL (SST Plastic)

134 pg/mL (Serum Plastic)

**Method Comparison 1:** The assay was compared to a commercially available assay for IL-6 (Kit A) on 78 patient samples. (Concentration range: approximately 3 to 300 pg/mL. See graph 1.) By linear regression:

(IML) = 0.75 (Kit A) + 5.2 pg/mL  
r = 0.940

Means:  
104 pg/mL (IMMULITE)  
132 pg/mL (Kit A)

**Method Comparison 2:** The assay was also compared to DPC's IMMULITE 2000 IL-6 (L2K6P) on 150 patient samples. (Concentration range: approximately 3 to 915 pg/mL. See graph 2.) By linear regression:

(IML) = 0.92 (IML 2000) + 6.1 pg/mL  
r = 0.994

Means:  
203 pg/mL (IMMULITE)  
215 pg/mL (IMMULITE 2000)

**Method Comparison 3:** A total of 110 specimens were tested by the IMMULITE/IMMULITE 1000 IL-6 assay on the two systems. (Concentration range: approximately 3 to 900 pg/mL. See graph 3.) By linear regression:

(IML 1000) = 0.98 (IML) – 3.0 pg/mL  
r = 0.994

Means:  
206 pg/mL (IMMULITE 1000)  
213 pg/mL (IMMULITE)

## References

1) Hirano T, Akira S, Taga T, Kishimoto T. Biological and clinical aspects of interleukin-6. *Immunol Today* 1990;11:443-9. 2) Saito S, Kasahara T, Kato Y, Ishihara Y, Ichijo M. Elevation of amniotic fluid interleukin 6 (IL-6), IL-8 and granulocyte colony stimulating factor (G-CSF) in term and preterm parturition. *Cytokine* 1993;5:81-8. 3) Waage A, Brandtzaeg P, Halstensen A, Kierulf P, Espevik, T. The complex pattern of cytokines in serum from patients with meningococcal septic shock. *J Exp Med* 1989; 169:333-8. 4) Hack CE, De Groot ER, Felt-Bersma RJF, Nuijens JH, Strack Van Schijndel RJM, Eerenberg-Belmer AJM, Thijs LG, Aarden LA. Increased plasma levels of interleukin-6 in sepsis. *Blood* 1989;74/5:1704-10. 5) Steinmetz HT, Herberitz A, Bertram M, Diehl V. Increase in interleukin-6 serum level preceding fever in granulocytopenia and correlation with death from sepsis. *J Infect Dis* 1995;171:225-8. 6) Buck C, Bundschu J, Gallati H, Bartmann P, Pohlandt F. Interleukin-6: A sensitive parameter for the early diagnosis of neonatal bacterial infection. *Pediatrics* 1994;93:54-8. 7) Hummel M, Czerlinski S, Friedel N, Liebenthal C, Hapser D, von Baehr R, Hetzer R, Volk HD. Interleukin-6 and interleukin-8 concentrations as predictors of outcome in ventricular assist device patients before heart transplantation. *Critical Care Med* 1994;22:448-54. 8) National Committee for Clinical

Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

## Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

## Tables and Graphs

### Precision (pg/mL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	88	3.1	3.5%	4.5	5.1%
2	112	5.8	5.2%	8.3	7.4%
3	230	9.1	4.0%	12.2	5.3%
4	388	17.1	4.4%	22.8	5.9%
5	638	29.6	4.6%	46.6	7.3%
6	1001	62.0	6.2%	74.9	7.5%

### Specificity

Compound <sup>1</sup>	pg/mL Added <sup>2</sup>	Apparent pg/mL <sup>3</sup>	% Cross-reactivity <sup>4</sup>
IL-1 $\alpha$	10,000	ND	ND
IL-1 $\beta$	10,000	ND	ND
IL-2	10,000	ND	ND
IL-4	10,000	ND	ND
IL-8	10,000	ND	ND
TNF $\alpha$	10,000	5.3	0.053%
Interferon Gamma	10,000	6.5	0.065%

ND: not detectable.<sup>5</sup>


### Linearity (pg/mL)



	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	103	—	—
	4 in 8	53.6	51.5	104%
	2 in 8	27.1	25.8	105%
	1 in 8	14.2	12.9	110%
2	8 in 8	418	—	—
	4 in 8	202	209	97%
	2 in 8	108	105	103%
	1 in 8	56.7	52.3	109%
3	8 in 8	516	—	—
	4 in 8	279	258	108%
	2 in 8	148	129	115%
	1 in 8	76.7	64.5	119%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
**EURO/DPC LTD**
Diagnostic Products Corporation
EC REP DPC Biemann GmbH