# IMMULITE®

## IL-8

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**DPC®**

# IMMULITE®/IMMULITE® 1000 IL-8

## English

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of interleukin 8 (IL-8) in serum or EDTA plasma, as an aid in the study of inflammatory diseases.

Catalog Number: **LK8PZ** (50 tests),  
**LK8P1** (100 tests)

Test Code: **I8** Color: **Dark Pink**

### Summary and Explanation

Interleukin 8 (IL-8) is a non-glycosylated protein with a molecular weight of 8 kDa. One of its most important biochemical functions is to act as a chemoattractant for neutrophils. Due to this property IL-8 was previously called neutrophil activating protein (NAP-1),<sup>1</sup> neutrophil activating factor (NAF)<sup>2</sup> and monocyte-derived neutrophil chemotactic factor (MDNCF).<sup>3</sup>

Initially, IL-8 was purified from monocytes which are believed to be the major source for the protein.<sup>2,3,4,5</sup> But it was shown that many other cell types, such as endothelial cells,<sup>6</sup> epithelial cells,<sup>7</sup> hepatocytes,<sup>7</sup> fibroblasts<sup>8</sup> and chondrocytes, can also produce IL-8. The release of the protein by these cells can be initiated by a wide variety of stimuli including LPS, IL-1, TNF, viruses and urate crystals.<sup>9</sup>

IL-8 is a member of the superfamily of chemokines, small 8–10 kDa proteins, which are referred to as pro-inflammatory cytokines. Besides their biological function as chemoattractants, members of these proteins show a 20–50% homology based on their amino acid sequence and a conserved four-cysteine motif. Other proteins which belong to this group are GRO $\alpha$ , RANTES, MCP-1, MIP-1 $\alpha$ , MIP-1 $\beta$ , IP-10, ENA-78, etc.<sup>10</sup>

Due to its biological properties, IL-8 plays a crucial role in many inflammatory events. As such, elevated levels of IL-8 are reported in various diseases including psoriasis,<sup>11</sup> cystic fibrosis,<sup>12</sup> idiopathic pulmonary fibrosis,<sup>12</sup> pleural disease<sup>12</sup> and rheumatoid arthritis.<sup>13</sup>

### Principle of the Procedure

IMMULITE/IMMULITE 1000 IL-8 is a solid-phase, two-site chemiluminescent immunometric assay.

**Incubation Cycles:** 1 × 30 minutes.

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 IL-8 has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 50  $\mu$ L serum or plasma. (Sample cup must contain at least 100  $\mu$ L more than the total volume required.) Avoid hemolysis; and separate the serum or plasma from the cells within two hours.<sup>15</sup>

**Storage:** 2 days at 2–8°C or for longer periods at –20°C.

### Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were

tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### IL-8 Test Units (L8P1)

Each barcode-labeled unit contains one bead coated with a murine monoclonal anti-IL-8 antibody. Stable at 2–8°C until expiration date.

**LK8PZ:** 50 units. **LK8P1:** 100 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### IL-8 Reagent Wedge (L8P2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to a rabbit polyclonal anti-IL-8 antibody in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

**LK8PZ:** 1 wedge. **LK8P1:** 1 wedge.

### IL-8 Adjustors (L8PL, L8PH)

Two vials (Low and High) of lyophilized IL-8 in a nonhuman serum matrix, with preservative. Reconstitute each vial with **3.0 mL** distilled or deionized water. Mix by gentle, intermittent swirling. Aliquot and freeze any unused material at –20°C within 3 hours of reconstitution. Stable at –20°C for 6 months after reconstitution.

**LK8PZ:** 1 set. **LK8P1:** 1 set.

## Kit Components Supplied Separately

### IL-8 Sample Diluent (L8PZ)

For the manual dilution of patient samples. One vial (25 mL) of an IL-8 and IL-10-free nonhuman serum matrix, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM :** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

**LILCM:** IMMULITE Cytokine Control Module (Bi-level, human serum-based).

Also Required

Sample transfer pipets, distilled or deionized water, controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

**Recommended Adjustment Interval:** 2 weeks.

**Quality Control Samples:** Use controls or sample pools with at least two levels (low and high) of IL-8.

## Expected Values

Samples from a total of 50 healthy individuals were assayed by the IMMULITE IL-8 procedure. The absolute range of these samples was from nondetectable to 62 pg/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without anticoagulants, gel barriers or clot-promoting additives.)

**Calibration Range:** Up to 7,500 pg/mL (NIBSC 89/520)

**Analytical Sensitivity:** 2 pg/mL

**Intraassay Precision (Within-Run):** Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

**Interassay Precision (Run-to-Run):** Statistics were calculated for samples assayed in 20 different runs. (See "Interassay Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three IL-8 solutions (158, 386 and 8,000 pg/mL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for interleukin 8. (See "Specificity" table.)

**Alternate Sample Type:** Samples ( $n = 5$ ) were collected into plain, heparinized and EDTA vacutainer tubes.

(EDTA) = 0.97 (Serum) + 13 pg/mL  
 $r = 1.00$

(Heparin) = 0.96 (Serum) – 17 pg/mL  
 $r = 1.00$

Means:  
3,403 pg/mL (Serum)  
3,339 pg/mL (EDTA)  
3,256 pg/mL (Heparin)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Method Comparison:** The assay was compared to the Milenia IL-8 assay on 32 patient samples. (Concentration range: approximately 19 to 6,700 pg/mL. See graph 1.) By linear regression:

(IML) = 1.04 (Milenia) – 62 pg/mL  
 $r = 0.982$

Means:  
2,086 pg/mL (IMMULITE)  
2,070 pg/mL (Milenia)

The IMMULITE IL-8 procedure (mono/poly, LK8P) was also compared to the IMMULITE IL-8 (mono/mono, LKL8) assay on 30 patient samples. (Concentration range: approximately 11 to 7,258 pg/mL. See graph 2.) By linear regression:

(LK8P) = 0.92 (LKL8) – 8.2 pg/mL  
 $r = 0.996$

Means:  
1,756 pg/mL (LK8P)  
1,908 pg/mL (LKL8)

## References

- 1) Larsen CG, et al. The neutrophil-activating protein (NAP-1) is also chemotactic for T lymphocytes. Science 1989;243:1464–6. 2) Walz AP, et al. Purification and amino acid sequencing of NAF, a novel neutrophil-activating factor produced by monocytes. Biochem Biophys Res Commun 1987;149:755–61. 3) Yoshimura TK, et al. Purification of a monocyte-derived neutrophil chemotactic factor that has peptide sequence similarity to other host defense cytokines. Proc Natl Acad Sci USA 84 1987;9233–7. 4) Peveri PA, et al. A novel neutrophil-activating factor produced by human mononuclear phagocytes. J Exp Med 1988;167:1547–59. 5) Damme van JJ, et al. A novel NH2-terminal sequence-characterized human monokine possessing neutrophil chemotactic, skin-reactive and granulocytosis-promoting activity. J Exp Med 1988;167:1364–76. 6) Schroeder JM and Christophers E. Secretion of novel and homologous neutrophil-activating peptides by LPS-stimulated human

endothelial cells. *J Immunol* 1989;142:244–51. 7) Kunkel SL, et al. Regulation of tumor necrosis factor- $\alpha$  and neutrophil activating protein-1 gene expression: potential role of cytokine directed cell communication during multiple organ injury. In: *Respiratory Distress Syndrome: Molecular to Man*, 1989. 8) Strieter RM, et al. Monokine-induced neutrophil chemotactic factor gene expression in human fibroblasts. *J Biol Chem* 1989;264:10621–6. 9) Baggiolini M, et al. Neutrophil-activating peptide-1/interleukin 8, a novel cytokine that activates neutrophils. *J Clin Invest* 1989;84:1045–7. 10) Sager R, et al. In: *J Westwick, et al, editors. Chemotactic Cytokines*. New York: Plenum Press, 1991:73. 11) Sticherling M, et al. Localization of neutrophil-activating peptide-1/interleukin 8-immunoreactivity in normal and psoriatic skin. *J Invest Dermatol* 1991;96:26–30. 12) Strieter RM, et al. The immunopathology of chemotactic cytokines: the role of interleukin-8 and monocyte chemoattractant protein-1. *J Lab Clin Med* 1994;123:183–97. 13) Peichl P, et al. Presence of NAP-1/IL-8 in synovial fluids indicates a possible pathogenic role in rheumatoid arthritis. *Scand J Immunol* 1991;34:333–9. 14) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard*. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 15) Data on file.

### Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

## Tables and Graphs

### Intraassay Precision (pg/mL)

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	89	3.2	3.6%
2	1,067	41	3.8%
3	4,187	156	3.7%

### Interassay Precision (pg/mL)

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	40	2.6	6.5%
2	120	8.9	7.4%
3	465	30	6.5%
4	1,018	53	5.2%
5	5,177	342	6.6%

### Linearity (pg/mL)

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	67	—	—
	4 in 8	34	34	100%
	2 in 8	19	17	112%
	1 in 8	10.5	8.4	125%
2	8 in 8	174	—	—
	4 in 8	84	87	97%
	2 in 8	38	44	86%
	1 in 8	19	22	86%
3	8 in 8	363	—	—
	4 in 8	188	182	103%
	2 in 8	91	91	100%
	1 in 8	45	45	100%

### Recovery (pg/mL)


	Spiking Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	64	—	—
	A	68	69	99%
	B	77	80	96%
	C	439	461	95%
2	—	147	—	—
	A	146	148	99%
	B	166	159	104%
	C	545	540	101%
3	—	290	—	—
	A	286	283	101%
	B	285	295	97%
	C	666	676	99%

## Specificity


Compound <sup>1</sup>	pg/mL Added <sup>2</sup>	Apparent pg/mL <sup>3</sup>	% Cross-reactivity <sup>4</sup>
IL-1 $\beta$	10,000	ND	ND
IL-2	10,000	ND	ND
IL2R	10,000	ND	ND
IL-4	10,000	ND	ND
IL-6	10,000	ND	ND
IFN	10,000	ND	ND
TNF $\alpha$	10,000	ND	ND
Pure Horse Serum	—	ND	ND

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
**EURO/DPC** LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH