IMMULITE®

High Sensitivity CRP

For use on the IMMULITE®
and IMMULITE® 1000 systems

DPC®

IMMULITE®/ IMMULITE® 1000

High Sensitivity CRP

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of C-reactive protein (CRP) in serum or plasma, as an aid in the detection and evaluation of the infection, tissue injury, inflammatory disorders and associated diseases.

Catalog Number: **LKCR1** (100 tests),
LKCR5 (500 tests)

Test Code: **CRP** Color: **Dark Blue**

CDC Analyte Identifier Code: 1001

CDC Test System Identifier Code: 10159

CLIA Complexity Category: Moderate

Summary and Explanation

CRP is an acute phase β -globulin with a molecular mass of approximately 118,000 daltons.¹ CRP is highly conserved, composed of five identical cyclic globular subunits, and is classified as a member of the pentraxin superfamily of proteins.^{1,2}

CRP levels in serum or plasma may rise during general, nonspecific response to infectious and non-infectious inflammatory conditions such as rheumatoid arthritis, cardiovascular disease and peripheral vascular disease.¹⁷ CRP is synthesized in the liver and is normally present as a trace constituent of serum or plasma.

While the precise *in vivo* functions of CRP during the inflammation state are not known, there is considerable evidence indicating a role in the recognition and elimination of foreign pathogens, as well as potentially toxic endogenous substances released from damaged tissue, by assisting humoral and cellular immunity.^{1,3} In normal individuals, hepatic cells produce CRP at constitutive levels of less than 1 mg/dL.³

Within 24 to 48 hours following acute tissue damage, production dramatically rises to approximately 1,000 times the constitutive level. CRP levels may remain

elevated for several days before returning to normal.^{1,3,4}

As elevated CRP values are always associated with pathological changes, the CRP assay provides useful information for the diagnosis, therapy, and monitoring of inflammatory conditions and associated diseases.⁷⁻⁹ Additionally, measurement of CRP by high sensitivity CRP assays may add to the predictive value of other markers used to assess the risk of cardiovascular and peripheral vascular diseases.¹¹⁻¹⁶

Principle of the Procedure

IMMULITE/IMMULITE 1000 High Sensitivity CRP is a solid-phase, chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE [IMMULITE 1000] High Sensitivity CRP has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 10 μ L prediluted sample: serum or heparinized or EDTA

plasma. (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: 3 days at 2–8°C or 2 months at –20°C.

Patient Sample Dilution: All samples must be prediluted 1-in-101 in CRP Sample Diluent (LCRZ4), e.g. by adding 10 µL sample to 1,000 µL Diluent.

Note that printed results have been corrected for the 1-in-101 dilution. Samples expected to have CRP levels outside of the calibration range should be further diluted. Correction for further dilution must be entered manually.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

CRP Test Units (LCR1)

Each barcode-labeled unit contains one bead coated with anti-ligand. Stable at 2–8°C until expiration date.

LKCR1: 100 units. **LKCR5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

CRP Reagent Wedge (LCR2)

With barcode. 7.5 mL ligand-labeled anti-CRP murine monoclonal antibody and alkaline phosphatase (bovine calf intestine) conjugated to rabbit polyclonal anti-CRP in buffer. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKCR1: 1 wedge. **LKCR5:** 5 wedges.

CRP Adjustors (LCRL, LCRH)

Two vials (Low and High), 2 mL each, of CRP in a protein/buffer matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C. **Do not predilute the Adjustors.**

LKCR1: 1 set. **LKCR5:** 2 sets.

CRP Sample Diluent (LCRZ4)

For the manual dilution of samples.

100 mL of a CRP-free protein/buffer matrix, with preservative, Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKCR1: 1 vial. **LKCR5:** 5 vials.

Kit Components Supplied Separately

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

LCRCM: Tri-level CRP Control Module

Also Required

Sample transfer pipets, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval: 2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of CRP.

Expected Values

A study performed on 100 apparently healthy laboratory volunteers yielded a median of 0.14 mg/dL and an upper 97.5th percentile of 1.1 mg/dL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

For diagnostic purposes, CRP results should be used in conjunction with other test results, the overall clinical presentation to the physician, and all other appropriate information.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in mg/dL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor: mg/dL \times 10 \rightarrow mg/L

Calibration Range: Up to 15 mg/dL (150 mg/L). WHO 1st IS 85/506.

Analytical Sensitivity: 0.01 mg/dL (0.1 mg/L).

Functional Sensitivity: <0.02 mg/dL (<0.2 mg/L).

High-dose Hook Effect: None up to 890 mg/dL.

Precision: Samples were repeatedly assayed in quadruplicate over the course of several days, for a total of 20 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples (prediluted 1-in-101) were assayed under various further dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three CRP solutions (60, 120 and 225 mg/dL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for CRP. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 570 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: Samples ($n = 16$) collected into plain and heparinized vacutainer tubes and spiked with CRP were assayed:

(Heparin) = 0.97 (Serum) + 0.29 mg/dL
 $r = 0.993$

Means:
13.1 mg/dL (Serum)
13.0 mg/dL (Heparin)

In a second study, samples ($n = 17$) collected into plain and EDTA vacutainer tubes and spiked with CRP were assayed.

(EDTA) = 1.03 (Serum) - 0.3 mg/dL
 $r = 0.995$

Means:
13.7 mg/dL (Serum)
13.8 mg/dL (EDTA)

Method Comparison 1: The IMMULITE High Sensitivity CRP procedure was compared to a commercially available immunoradiometric assay (Kit A) on 97 patient samples, (Concentration range: nondetectable to approximately

8.2 mg/dL. See graph.) By linear regression:

(IMMULITE) = 1.04 (Kit A) – 0.01 mg/dL
r = 0.989

Means:
0.89 mg/dL (IMMULITE)
0.87 mg/dL (Kit A)

Method Comparison 2: The assay was compared to another commercially available assay (Kit B) on 199 samples. (Concentration range: nondetectable to approximately 2.5 mg/dL. See graph.) By linear regression:

(IML) = 0.953 (Kit B) + 0.002 mg/dL
r = 0.996

Means:
0.23 mg/dL (IMMULITE)
0.24 mg/dL (Kit B)

References

1) Schultz DR, Arnold PI. Properties of four acute phase proteins: CRP, serum amyloid A protein, α_1 -acid glycoprotein, and fibrinogen. *Semin Arthritis Rheum* 1990;20:129-47. 2) Osmand AP, Friedenson B, Gewurz H, et al. Characterization of CRP and the complement subcomponent C1t as homologous proteins displaying cyclic pentameric symmetry (pentraxins). *Proc Nat Acad Sci USA* 1977;74:739-43. 3) Tietz NW, editor. *Clinical guide to laboratory tests*. 3rd ed. Philadelphia: WB Saunders, 1995:358. 4) Kushner I. The phenomenon of acute phase response. *Ann N Y Acad Sci* 1982;389:39-48. 5) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture: approved standard*. 4th ed. NCCLS Document H3-A4. Wayne, PA: NCCLS, 1998. 6) Kushner I, Rzewnicki DL. The acute phase response: general aspects. *Baillieres Clin Rheumatol* 1994;8:513-30. 7) Dixon JS, Bird HA, Sitton NG, et al. CRP in the serial assessment of disease activity in rheumatoid arthritis. *Scand J Rheumatol* 1984;13:39-44. 8) Hind CR, Pepys MB. The role of CRP measurement in clinical practice. *Int Med* 1984;5:112-51. 9) van Leeuwen MA, van Rijswijk MH. Acute phase proteins in the monitoring of inflammatory disorders. *Baillieres Clin Rheumatol* 1994;8:531-52. 10) Ridker PM, Cushman M, Stampfer MJ, et al. Plasma concentration of CRP and risk of developing peripheral vascular disease. *Circulation* 1998;97:425-8. 11) Ridker PM, Buring JE, Shih J, et al. Prospective study of CRP and the risk of future cardiovascular events among apparently healthy women. *Circulation* 1998;98:731-3. 12) Ridker PM, Glynn RJ, Hennekens CH. CRP adds to the predictive value of total and HDL cholesterol in determining risk of first myocardial infarction. *Circulation* 1998;97:2007-11. 13)

Ridker PM, Rifai N, Pfeffer MA, et al. Inflammation, pravastatin, and the risk of coronary events after myocardial infarction in patients with average cholesterol levels. *Cholesterol and Recurrent Events (CARE) Investigators. Circulation* 1998;98:839-44. 14) Tracy RP, Lemaitre RN, Psaty BM, et al. Relationship of CRP to risk of cardiovascular disease in the elderly. Results from the Cardiovascular Health Study and the Rural Health Promotion Project. *Arterioscler Thromb Vasc Biol* 1997;17:1121-7. 15) Macy EM, Hayes TE, Tracy RP. Variability in the measurement of CRP in healthy subjects: implications for reference intervals and epidemiological applications. *Clin Chem* 1997;43:52-8. 16) Ridker PM, Cushman M, Stampfer MJ, et al. Inflammation, aspirin, and the risk of cardiovascular disease in apparently healthy men. *N Engl J Med* 1997;336:973-9. Erratum: *ibid* 1997;337:356. 17) Du Clos TW. Function of CRP. *Ann Med* 2000; 32:274-8.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Precision (mg/dL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	0.02	0.001	5.0%	0.002	10.0%
2	0.04	0.002	5.0%	0.004	10.0%
3	0.08	0.004	5.0%	0.006	7.5%
4	0.15	0.008	5.3%	0.009	6.0%
5	0.31	0.013	4.2%	0.015	4.8%
6	1.5	0.063	4.2%	0.073	4.9%
7	7.8	0.465	6.0%	0.570	7.3%
8	14.8	0.951	6.4%	1.14	7.7%

Linearity (mg/dL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	16 in 16 ⁵	1.17	—	—
	8 in 16	0.58	0.59	98%
	4 in 16	0.28	0.29	97%
	2 in 16	0.13	0.15	87%
	1 in 16	0.07	0.07	100%
2	16 in 16	4.48	—	—
	8 in 16	2.08	2.24	93%
	4 in 16	1.07	1.12	96%
	2 in 16	0.51	0.56	91%
	1 in 16	0.28	0.28	100%
3	16 in 16	8.67	—	—
	8 in 16	4.32	4.34	100%
	4 in 16	2.09	2.17	96%
	2 in 16	1.05	1.08	97%
	1 in 16	0.51	0.54	94%
4	16 in 16	11.4	—	—
	8 in 16	5.70	5.70	100%
	4 in 16	2.55	2.85	89%
	2 in 16	1.20	1.43	84%
	1 in 16	0.65	0.71	92%

Recovery (mg/dL)


	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	0.46	—	—
	A	3.52	3.44	102%
	B	6.06	6.44	94%
	C	12.2	11.7	104%
2	—	0.97	—	—
	A	3.80	3.92	97%
	B	6.69	6.92	97%
	C	11.4	12.2	93%
3	—	1.66	—	—
	A	4.34	4.58	95%
	B	7.01	7.58	93%
	C	12.6	12.8	98%

Specificity


Compound ¹	Amount Added mg/dL ²	% Cross-reactivity ³
Serum Amyloid A	1,000	ND
Human Serum Albumin	1,000	ND
Human IgG	1,000	ND
Transferrin	1,000	ND

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH