**IMMULITE®**

**IL2R**

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**DPC®**

# IMMULITE/IMMULITE 1000® IL2R

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of soluble interleukin 2 receptor (IL2R) in serum and EDTA plasma, as an aid in the study of inflammatory diseases.

Catalog Number: **LKIPZ** (50 tests), **LKIP1** (100 tests), **LKIP5** (500 tests)

Test Code: **I2R** Color: **Light Blue**

### Summary and Explanation

The receptor of the cytokine interleukin 2 (IL-2) plays a crucial role in the regulation of the immune response. Binding of IL-2 to its receptor (IL2R) on the surface of T-lymphocytes triggers a series of intracellular signaling events that result in the activation and proliferation of resting T cells and, ultimately, in the generation of helper, suppressor and cytotoxic T cells, which mediate immune reactions.

The IL-2 receptor is unique among growth factor receptors in that it is made up of at least three distinct membrane components: the  $\alpha$  chain (IL2R $\alpha$ ) — the so-called Tac antigen (MW 55kDa), a  $\beta$ -chain (IL2R $\beta$ ; MW 70-75 kDa), and the  $\gamma$ -chain (IL2R $\gamma$ ; MW 64 kDa). Different combinations of these three components give rise to the generation of various forms of the IL-2 receptor, each of which manifests different binding affinities to IL-2.

Most resting T cells, B cells, large granular lymphocytes (LGLs) and monocytes do not express significant numbers of this receptor on their surfaces. Upon activation, receptor molecules are expressed on the surface of the cells, and a soluble form (sIL2R) is released, which is about 10 kDa smaller than the membrane-bound protein.

It has been found that sIL2R is present at low levels in serum of healthy individuals and at significantly elevated levels in a broad range of disorders such as neoplastic diseases, autoimmune

diseases, organ allograft rejection and different infections. Thus, it appears that sIL2R can serve as a marker for diagnosis, therapeutic evaluation and management of cancer, as well as an indicator of a wide spectrum of disorders involving immune activation.

### Principle of the Procedure

IMMULITE/IMMULITE 1000 IL2R is a solid-phase, two-site chemiluminescent immunometric assay.

**Incubation Cycles:** 1 × 30 minutes.

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 IL2R has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 50  $\mu$ L serum or EDTA plasma. (Sample cup must contain at least 100  $\mu$ L more than the total volume required.)

**Storage:** 2 days at 2–8°C or for longer at –20°C.

### Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### IL2R Test Units (LIP1)

Each barcode-labeled unit contains one bead coated with a murine monoclonal anti-IL2R antibody. Stable at 2–8°C until expiration date.

**LKIPZ:** 50 units. **LKIP1:** 100 units.

**LKIP5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### IL2R Reagent Wedge (LIP2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to rabbit polyclonal anti-IL2R antibody in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

**LKIPZ:** 1 wedge. **LKIP1:** 1 wedge.

**LKIP5:** 5 wedges.

### IL2R Adjustors (LIPL, LIPH)

Two vials (Low and High) of lyophilized IL2R in a protein-based buffer, with preservative. Reconstitute each vial with **3.0 mL** distilled or deionized water. Mix by

gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 3 months after reconstitution. Do not freeze. Following reconstitution, the Adjustors should be aliquoted, to avoid contamination from repeated usage.

**LKIPZ:** 2 sets. **LKIP1:** 2 sets.

**LKIP5:** 3 sets.

## Kit Components Supplied Separately

### IL2R Sample Diluent (LIPZ)

For the manual dilution of patient samples. One vial containing 25 mL of an IL2R-free protein-based buffer, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM :** Probe Cleaning Kit

**LCHX-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

**LILCM:** IMMULITE Cytokine Control Module (Bi-level, human serum-based).

Also Required

Sample transfer pipets, distilled or deionized water, controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manuals.

See the IMMULITE or IMMULITE 1000 Operator's Manuals for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

**Recommended Adjustment Interval:** 4 weeks.

**Quality Control Samples:** Use controls or sample pools with at least two levels (low and high) of IL2R.

## Expected Values

A study performed on 87 apparently healthy adults yielded a median of 391 U/mL and a 95% reference range of

223 to 710 U/mL

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

### Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

### Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in U/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Calibration Range:** Up to 7,500 U/mL

**Analytical Sensitivity:** 5 U/mL

**High-dose Hook Effect:** None up to 225,000 U/mL.

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three IL2R solutions (7,280, 13,780 and 27,440 U/mL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for IL2R. (See "Specificity" table.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no

effect on results, within the precision of the assay.

**Hemolysis:** Presence of hemoglobin in concentrations up to 384 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from 9 volunteers into plain, heparinized, EDTA and Becton Dickinson SST<sup>®</sup> vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of IL2R, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE IL2R procedure.

(EDTA) = 0.97 (Serum) + 120 U/mL  
r = 0.988

(Heparin) = 1.04 (Serum) – 4.8 U/mL  
r = 0.999

(SST) = 0.98 (PlainTubes) + 25.5 U/mL  
r = 1.00

Means:

2,234 U/mL (Serum)  
2,186 U/mL (EDTA)  
2,321 U/mL (Heparin)  
2,293 U/mL (SST)

**Method Comparison:** The assay (LKIP, mono/poly) was compared to DPC's IMMULITE IL2R (LK2R, mono/mono) on 66 samples. (Concentration range: approximately 200 to 2,500 U/mL. See graph.) By linear regression:

(LKIP) = 0.69 (LK2R) + 54 U/mL  
r = 0.986

Means:

832 U/mL (LKIP)  
1,125 U/mL (LK2R)

### References

- 1) Rubin LA, et al. Identification and characterization of a released form of the Interleukin-2 Receptor. Leukocytes and host defense 1986;95-102. 2) Rubin LA, et al. Soluble interleukin-2 receptors are released from activated human lymphoid cells *in vitro*. J Immunol 1985;135:3172-7. 3) Rubin LA, et al. The released interleukin-2 receptor binds interleukin-2 efficiently. J Immunol 1986;137:3841-4. 4) Pizzolo G. The soluble Interleukin-2 receptor as a new biological marker in diseases. Immunol Clin 1988;VII:1. 5) Pui Ch-

H. Serum Interleukin-2 Receptor: Clinical and biological Implications. *Leukemia* 1989;3:323-7. 6) Rubin LA. The soluble Interleukin-2 Receptor: biology function, and clinical application. *Annals Internal Med* 1990;113:619-27. 7) Zerler B. The soluble Interleukin-2 Receptor as a marker for Human Neoplasia and Immune status. *Cancer Cells* 1991;3:471-9. 8) Taniguchi T, Minami Y. The IL-2/IL-2 Receptor System: A current overview. *Cell* 1993;73:5-8. 9) Waldmann ThA. The IL-2/IL-2 receptor system: a target for rational immune intervention. *Immunol Today* 1993;14: 264-70. 10) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

### Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

## Tables and Graphs

### Precision (U/mL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	213	6.7	3.2%	13	6.1%
2	411	15	3.7%	26	6.3%
3	752	25	3.3%	49	6.5%
4	1,167	41	3.5%	85	7.3%
5	1,795	57	3.2%	135	7.5%
6	2,463	88	3.6%	189	7.7%
7	3,015	107	3.6%	241	8.0%
8	3,763	109	2.9%	303	8.1%

### Linearity (U/mL)

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	4 in 4 <sup>5</sup>	345	—	—
	2 in 4	185	173	107%
	1 in 4	88	86	102%
2	8 in 8	842	—	—
	4 in 8	455	421	108%
	2 in 8	224	211	106%
3	8 in 8	1,267	—	—
	4 in 8	655	634	103%
	2 in 8	340	317	107%
4	8 in 8	1,492	—	—
	4 in 8	816	746	109%
	2 in 8	423	373	113%
5	8 in 8	208	187	111%
	8 in 8	3,196	—	—
	4 in 8	1,540	1,598	96%
6	2 in 8	831	799	104%
	1 in 8	388	400	97%
	8 in 8	3,327	—	—
7	4 in 8	1,703	1,664	102%
	2 in 8	840	832	101%
	1 in 8	426	416	102%

## Recovery (U/mL)


	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	332	—	—
	A	586	679	86%
	B	874	1,004	87%
	C	1,372	1,687	81%
2	—	1,267	—	—
	A	1,450	1,568	92%
	B	1,765	1,893	93%
	C	2,418	2,576	94%
3	—	1,492	—	—
	A	1,689	1,781	95%
	B	2,055	2,106	98%
	C	2,495	2,789	89%
4	—	2,348	—	—
	A	2,579	2,595	99%
	B	2,835	2,920	97%
	C	3,441	3,603	96%
5	—	3,196	—	—
	A	3,184	3,400	94%
	B	3,398	3,725	91%
	C	3,934	4,408	89%
6	—	3,327	—	—
	A	3,237	3,524	92%
	B	3,650	3,850	95%
	C	4,237	4,533	93%

## Specificity

Compound <sup>1</sup>	pg/mL Added <sup>2</sup>	Apparent U/mL <sup>3</sup>	% Cross-reactivity <sup>4</sup>
IL-1 $\alpha$	10,000	ND	ND
IL-2	10,000	ND	ND
IL-4	10,000	ND	ND
IL-6	10,000	ND	ND
IL-8	10,000	ND	ND
IL-1 $\beta$	10,000	ND	ND
Interferony	10,000	ND	ND
TNF $\alpha$	10,000	ND	ND

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
**EURO/DPC** LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH