IMMULITE®

IL-1 β

For use on the IMMULITE®
and IMMULITE® 1000 systems

DPC®

IMMULITE®/IMMULITE® 1000 IL-1 β

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of interleukin-1 β (IL-1 β) in serum or heparinized plasma, as an aid in the clinical assessment of various immunological diseases.

Catalog Number: **LKL1Z** (50 tests), **LKL11** (100 tests)

Test Code: **I1B** Color: **Aqua**

Summary and Explanation

Interleukin-1 is an inflammatory cytokine, mainly produced by activated macrophages and peripheral neutrophils, but other cells like smooth muscle cells, astrocytes, fibroblasts, keratinocytes, T cells, B cells and NK cells are capable of IL-1 synthesis, too. The biological activity of IL-1 is exhibited in two forms: IL-1 α and IL-1 β .

The gene sequences of IL-1 α and IL-1 β were published in 1984,^{1,2} and it was determined that both genes code for proteins which show only 27% homology in their amino acid sequence. Both forms are produced as precursor proteins, showing a molecular weight of about 35 kD. In contrast to proIL-1 β , proIL-1 α shows biological activity. ProIL-1 α , as well as proIL-1 β , is proteolytically cleaved to form the 17 kD mature proteins. Because both forms lack a leader sequence, most of the IL-1 remains in the cytosol. Only in mononuclear phagocytes, between 40 and 60% of the IL-1 β is released from the cell.

Both forms of IL-1, IL-1 α and IL-1 β , bind to the same receptor^{3,4} on the cell surface of their target cells and therefore mediate similar biological actions. Different types of receptors occur on different cells. An 80 kD receptor, also called Type I receptor or p80, is primarily located on T cells, whereas the 60 kD Type II receptor can be found on B cells, granulocytes and monocytes.

The synthesis of IL-1 can be induced by a combination of other cytokines,

endotoxins, viruses, mitogens and antigens. In contrast, prostaglandin E2 and corticosteroids, lipoproteins, lipids, α 2-macroglobulin and a naturally occurring antagonist, called IL-1 receptor antagonist, inhibit IL-1 synthesis.

IL-1 mediates a wide variety of biological actions. It stimulates the production and secretion of IL-2 and the expression of IL-2 receptors by helper cells. IL-1 acts synergistically with other factors in the activation and differentiation of B cells to immunoglobulin secreting cells. Furthermore, it stimulates the activation and differentiation of NK cells, fibroblasts and thymocytes. IL-1 acts antiproliferatively on many tumor cell types, increases the tumor cytotoxicity of macrophages and induces tumor regression.

In synergy with TNF α , IL-1 activates osteoclasts and therefore plays an important role in the regulation of bone metabolism.⁵ IL-1 has different effects on the central nervous system. It is an endogenous pyrogen and causes fever in humans at doses less than 1 ng/kg.⁶ Additionally, it induces the synthesis of ACTH, endorphins, vasopressin and somatostatin.⁷

Due to the biological properties of interleukin-1 described above, the measurement of this mediator is of interest for the monitoring and diagnosis of bone diseases, inflammatory diseases, various types of cancer and many other immunological diseases.

Principle of the Procedure

IMMULITE/IMMULITE 1000 IL-1 β is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 \times 60 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 IL-1 β has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 75 μ L serum or heparinized plasma. (Sample cup must contain at least 100 μ L more than the total volume required.)

Storage: 2 days at 2–8°C or 6 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

IL-1 β Test Units (LL11)

Each barcode-labeled unit contains one bead coated with a murine monoclonal anti-IL-1 β antibody. Stable at 2–8°C until expiration date.

LKL1Z: 50 units. **LKL11:** 100 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

IL-1 β Reagent Wedge (LL12)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to a rabbit polyclonal anti-IL-1 β antibody in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKL1Z: 1 wedge. **LKL11:** 2 wedges.

IL-1 β Adjustors (LL1L, LL1H)

Two vials (Low and High) containing lyophilized IL-1 β in a nonhuman serum/buffer matrix, with preservative. Reconstitute each vial with 3.0 mL distilled or deionized water. Mix by gentle, intermittent swirling. Aliquot and freeze immediately. Stable at –20°C for 2 months after reconstitution.

LKL1Z: 1 set. **LKL11:** 1 set.

Kit Components Supplied Separately

IL-1 β Sample Diluent (LL1Z)

For the manual dilution of patient samples. One vial 25 mL containing an IL-1 β -free nonhuman serum/buffer matrix. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM : Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

LILCM: IMMULITE Cytokine Control Module (Bi-level, human serum-based).

Also Required

Sample transfer pipets, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of IL-1 β .

Expected Values

Samples from a total of 47 healthy individuals were assayed by the IMMULITE IL-1 β procedure. All results were below 5 pg/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless

otherwise specified, all results were generated on serum samples collected in tubes without anticoagulants, gel barriers, or clot-promoting additives.)

Calibration Range: up to 1,000 pg/mL (1st IS 86/680)

Analytical Sensitivity: 1.5 pg/mL

High-dose Hook Effect: None up to 100,000 pg/mL.

Intraassay Precision (Within-Run): Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run): Statistics were calculated for samples assayed in 10 different runs. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three IL-1 β solutions (424, 1,058 and 1,923 pg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for IL-1 β . (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 381 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: Samples ($n = 10$) were collected into plain and heparinized vacutainer tubes.

(Heparin) = 1.02 (Serum) + 2.64 pg/mL
 $r = 0.996$

Means:
298 pg/mL (Serum)
320 pg/mL (Heparin)

EDTA and citrate plasma samples are not recommended in the IMMULITE IL-1 β procedure since they lead to erroneous results.

Method Comparison: The IMMULITE IL-1 β procedure was compared to Milenia

IL-1 β assay on 38 patient samples.
(Concentration range: approximately 13 to 604 pg/mL. See graph.) By linear regression:

$$(IML) = 1.1 (\text{Milenia}) + 6.4 \text{ pg/mL}$$

$$r = 0.995$$

Means:
150 pg/mL (IMMULITE)
136 pg/mL (Milenia)

References

1) Lomedico PT, Gubler R, Hellmann CP, Dukovich M, Giri JG, Pan JE, Collier K, Semionow R, Chua AO, Mizel SB. Cloning and expression of murine interleukin-1 cDNA in *Escherichia coli*. *Nature* 1984;312:458. 2) Auron PE, Webb AC, Rosenwasser LJ, Mucci SF, Rich A, Wolff SM, Dinarello CA. Nucleotide sequence of human monocyte interleukin-1 precursor cDNA. *Proc Natl Acad Sci USA* 1984;81:7907. 3) Sims JE, March CJ, Cosman D, Widmer MB, MacDonald HR, McMahan CJ, Grubin CE, Wignall JM, Jackson JL, Call SM, Friend B, Alpert AR, Gillis S, Urdal DL, Dower SK. cDNA expression cloning of the IL-1 receptor, a member of the immunoglobulin superfamily. *Science* 1988;241:585. 4) Chizzonite R, Truitt T, Kilian P, Kern A, Nunes P, Parker KP, Kaffka KL, Chua AO, Lugg DK, Gubler U. Two high affinity interleukin-1 receptors represent separate gene products. *Proc Natl Acad Sci USA* 1989;86:8029. 5) Konig A, Muhlbauer RC, Fleisch H. Tumor necrosis factor and interleukin-1 stimulate bone resorption *in vivo* as measured by ³H-tetracycline excretion from prelabeled mice. *J Bone Miner Res* 1988;3:621. 6) Tewari A, Buhles WC Jr, Starnes HF Jr. Preliminary report: effects of interleukin-1 on platelet counts. *Lancet* 1990;336: 712. 7) Rerhaug A, Michie HR, Manson JM, Watters JM, Dinarello CA, Wolff SM, Wilmore DW. Inhibition of cyclooxygenase attenuates the metabolic response to endotoxin in humans. *Arch Surg* 1988;123:162. 8) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

Available outside the United States only.
For technical assistance, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Intraassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³
1	39	1.1	2.8%
2	190	9.3	4.9%
3	669	32	4.8%

Interassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³
1	13	1.0	7.7%
2	33	3.0	9.1%
3	57	3.7	6.5%
4	111	5.4	4.9%
5	206	10.3	5.0%
6	397	23	5.8%
7	785	38	4.8%
8	973	56	5.8%

Recovery (pg/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	10.8	—	—
	A	32	31	103%
	B	69	63	110%
	C	121	106	114%
2	—	126	—	—
	A	141	141	100%
	B	163	173	94%
	C	213	216	99%
3	—	369	—	—
	A	372	372	100%
	B	387	403	96%
	C	437	447	98%

Specificity


Compound ¹	pg/mL Added ²	Apparent pg/mL ³	% Cross-reactivity ⁴
IL-1 α	10,000	ND	ND
IL-2	10,000	ND	ND
IL2R	10,000	ND	ND
IL-4	10,000	ND	ND
IL-5	10,000	ND	ND
IL-6	10,000	ND	ND
IL-8	10,000	ND	ND
IL-10	10,000	ND	ND
IFN γ	10,000	ND	ND
TNF α	10,000	ND	ND
Pure Horse Serum	—	ND	ND

ND: not detectable.⁵

Linearity (pg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	42.2	—	—
	4 in 8	20.1	21.1	95%
	2 in 8	10.2	10.6	96%
	1 in 8	4.8	5.3	91%
2	8 in 8	184	—	—
	4 in 8	84	92	91%
	2 in 8	43	46	93%
	1 in 8	22	23	96%
3	8 in 8	190	—	—
	4 in 8	98	95	103%
	2 in 8	45	48	94%
	1 in 8	23	24	96%
4	8 in 8	980	—	—
	4 in 8	500	490	102%
	2 in 8	269	245	110%
	1 in 8	142	123	115%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation