IMMULITE®

Androstenedione

**For use on the IMMULITE®
and IMMULITE® 1000 systems**

DPC®

IMMULITE/IMMULITE 1000[®]

Androstenedione

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of Δ^4 -androstenedione in human serum.

Catalog Numbers: **LKA01** (100 tests), **LKA05** (500 tests)

Test Code: **AND** Color: **Brown**

Summary and Explanation

Δ^4 -Androstenedione is a steroid which serves as a major precursor for testosterone and estrone.^{1,2,3} Its clinical interest derives from the fact that it is often elevated in cases of abnormal hair growth (hirsutism) and virilization.^{1,6,8-12}

Unlike the adrenal androgens dehydroepiandrosterone and its sulfate, circulating androstenedione originates both from the adrenals and from the ovaries.^{1,12,13} Plasma levels increase steadily from about the seventh year of life, then gradually decline after the third decade.^{2,4,13} Androstenedione exhibits a diurnal variation, being highest in the morning, and also a cyclical variation during the menstrual cycle, being highest near midcycle.^{1,5,14,15} During pregnancy, there is an increase in the plasma level.^{3,11}

Principle of the Procedure

IMMULITE/IMMULITE 1000 Androstenedione is a solid-phase, competitive chemiluminescent enzyme immunoassay.

Incubation Cycles: 1 × 60 minutes.

Specimen Collection

Neither EDTA nor heparinized plasma is recommended for use.

Lipemic, hemolyzed, icteric or grossly contaminated samples may give erroneous results. The use of an ultracentrifuge is recommended to clear lipemic samples.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on tube materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Androstenedione has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 25 μ L serum. (Sample cup must contain at least 100 μ L more than the total volume required.)

Storage: 24 hours at 2–8°C or 2 months (aliquotted) at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Androstenedione Test Units (LAO1)

Each barcode-labeled unit contains one bead coated with polyclonal rabbit anti-androstenedione antibody. Stable at 2–8°C until expiration date.
LKA01: 100 units. **LKA05:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

Androstenedione Reagent Wedge (LAO2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to androstenedione in buffer. Store capped and refrigerated: stable at 2–8°C until expiration date.

Recommended usage is within 30 days after opening when stored as indicated.
LKA01: 1 wedge. **LKA05:** 5 wedges.

Androstenedione Adjustors (LAOL, LAOH)

Two vials (Low and High), 2 mL each, of androstenedione in a processed human serum matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKA01: 1 set. **LKA05:** 2 sets.

Kit Components Supplied Separately

LSUBX: Chemiluminescent Substrate
LPWS2: Probe Wash Module
LKPM: Probe Cleaning Kit
LCHx-y: Sample Cup Holders (barcoded)
LSCP: Sample Cups (disposable)
LSCC: Sample Cup Caps (optional)

CON6: Tri-level, multi-constituent control

Also Required

Sample transfer pipets, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of androstenedione.

Expected Values

A reference range study was performed on serum samples from 48 male and 58 female adult laboratory volunteers using the IMMULITE/IMMULITE 1000 Androstenedione procedure.

Males: Median: 1.8 ng/mL (6.3 nmol/L)
Central 95% range: 0.7 – 3.6 ng/mL (2.4 – 12.6 nmol/L)

Females: Median: 1.9 ng/mL (6.6 nmol/L)
Central 95% range: 0.3 – 3.5 ng/mL (1.0 – 12.2 nmol/L)

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

The growing graafian follicle secretes androstenedione, and a two-fold increase in its production occurs near the midcycle. Specimens drawn during this time may therefore yield androstenedione results which fall above the upper central 95% limit.^{1,2}

Diurnal variation has been observed. For example, in 11 healthy men 18 to 30 years of age, the mean 24-hour androstenedione concentration was 0.93 ng/mL, with a circadian amplitude equaling 42% of the mean. The major acrophase occurred in the morning at 07:10 hours; the major nadir occurred in the evening at 23:00; and there was a secondary acrophase at 16:35. The 95% confidence interval for acrophase/nadir was 68 minutes.¹³ Accordingly, the time of collection should be noted.

It should be remembered that exercise can have a significant impact on circulating androstenedione levels. Thus, in a study of 18 professional soccer players, the mean androstenedione concentration increased from a baseline of approximately 0.9 ng/mL to 1.8 ng/mL at end of game, declining to 0.7 ng/mL ninety minutes afterwards.¹⁴ In another study, the

mean androstenedione level for 20 moderately trained men increased significantly, from a baseline mean of 1.11 ng/mL, to 1.41 ng/mL five minutes after a one-hour swim.¹⁵

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data representative of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

ng/mL \times 3.4916 \rightarrow nmol/L

Calibration Range: 0.3 – 10 ng/mL
(1.04 to 35 nmol/L).

Analytical Sensitivity: 0.3 ng/mL
(1.0 nmol/L).

Precision: Samples were assayed in duplicate over the course of multiple days, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with 3 androstenedione solutions (10, 40 and 100 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for androstenedione, with an extremely low crossreactivity to other naturally occurring steroids or therapeutic drugs

that may be present in patient samples. (See "Specificity" table.)

Bilirubin: Samples spiked with 100 and 200 mg/L of conjugated and unconjugated bilirubin were analyzed. Bilirubin may interfere with the assay, causing elevations. (See "Bilirubin" table.)

Hemolysis: Samples spiked with 135, 270 and 540 mg/dL of hemoglobin were analyzed. Hemolysis may interfere with the assay, causing elevations. (See "Hemolysis" table.)

Lipemia: May cause a depression of values. (See "Lipemia" table.)

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 26 volunteers into plastic plain serum, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of androstenedione, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE Androstenedione procedure.

(Heparin) = 0.88 (Serum) – 0.15 ng/mL
 $r = 0.981$

(SST) = 1.04 (Serum) – 0.10 ng/mL
 $r = 0.989$

(EDTA) = All samples were >10 ng/mL.

Means:

3.83 ng/mL (Serum Plastic)
3.20 ng/mL (Heparin Plastic)
3.87 ng/mL (SST Plastic)
>10 ng/mL (EDTA Plastic)

Neither EDTA nor heparinized plasma are recommended for use.

Method Comparison 1: The assay was compared to a commercially available solid phase radioimmunoassay for androstenedione (Kit A) on 78 endogenous patient samples. (Concentration range: approximately 0.5 to 6 ng/mL. See graph 1.) By linear regression:

(IML) = 1.17 (Kit A) – 0.37 ng/mL
 $r = 0.924$

Means:

1.83 ng/mL (IMMULITE)
1.88 ng/mL (Kit A)

Method Comparison 2: The assay was also compared to DPC's IMMULITE 2000 Androstenedione on 334 endogenous and

spiked samples. (Concentration range: approximately 0.3 to 9.5 ng/mL. See graph 2.) By linear regression:

$$(IML) = 1.04 (IML\ 2000) + 0.09\ ng/mL$$

$$r = 0.992$$

Means:

2.69 ng/mL (IMMULITE)
2.49 ng/mL (IMMULITE 2000)

Method Comparison 3: A total of 95 patient samples were tested by the IMMULITE/IMMULITE 1000 Androstenedione assay on the two systems. (Concentration range: approximately 0.4 to 10 ng/mL. See "Method Comparison 3" graph). By linear regression:

$$(IML\ 1000) = 0.99 (IML) + 0.04\ ng/mL$$

$$r = 0.991$$

Means:

3.53 ng/mL (IMMULITE 1000)
3.52 ng/mL (IMMULITE)

References

- 1) Abraham GE. Ovarian and adrenal contribution to peripheral androgens during the menstrual cycle. *J Clin Endocrinol Metab* 1974;39:340.
- 2) Baird DT. Ovarian steroid secretion and metabolism in women. In: James VHT, Serio M, Giusti G, editors. *The Endocrine Function of the Human Ovary*. New York: Academic Press, 1976: 125-33.
- 3) Fietz J et al. Simultaneous radioimmunoassay of androstenedione, dehydroepiandrosterone and 11 β -hydroxyandrostenedione in plasma. *Horm Res* 1980;13: 133-49.
- 4) Hummer L et al. An easy and reliable radioimmunoassay of serum androstenedione: age-related normal values in 252 females aged 2 to 70 years. *Scand J Clin Lab Invest* 1983;43:301-6.
- 5) Judd HL, Yen SSC. Serum androstenedione and testosterone levels during the menstrual cycle. *J Clin Endocrinol Metab* 1973;39:475.
- 6) Keenan BS et al. Plasma androgens in congenital adrenal hyperplasia; androstenedione concentration as an index of adrenal androgen suppression. *J Lab Clin Med* 1979;94:799-808.
- 7) Putz Z. A selective radioimmunoassay of androstenedione in plasma and saliva. *J Clin Chem Clin Biochem* 1982;20:761-4.
- 8) Raj S et al. Diagnostic value of androgen measurements in polycystic ovary syndrome. *Obstet Gynecol* 1978;52:169-71.
- 9) Redi B et al. Diagnostic strategy in hyperandrogenic syndrome. *Horm Res* 1983;18:117-24.
- 10) Sciarra F. Diagnosis of virilizing syndromes: endocrinological parameters. In: Milinatti G et al, editors. *Androgenization in Women*. New York: Raven Press 1983;85-115.
- 11) Tourniaire J, Pugeat M. Strategic approach of hyperandrogenism in women. *Horm Res* 1983;18:125-34.
- 12) Vermeulen A. Androgen secretion by adrenals

and gonads. In: Mahesh V, Greenblatt RB, editors. *Hirsutism and Virilism*. Boston: John Wright PSG 1983;17-34.
- 13) Lejeune-Lenain C, Van Cauter E, Desir D, et al. Control of circadian and episodic variations of adrenal androgens secretion in man. *J Endocrinol Invest* 1987;10:267-76.
- 14) Lupo C, Baldi L, Bonifazi M, et al. Androgen levels following a football match. *Eur J Appl Physiol* 1985;54:494-6.
- 15) Velardo A, Pantaleoni M, Valerio L, et al. Influence of exercise on dehydroepiandrosterone sulphate and d₄-androstenedione plasma levels in man. *Exp Clin Endocrinol* 1991;97:99-101.
- 16) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard*. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	0.66	0.06	9.1%	0.10	15.2%
2	1.57	0.08	5.1%	0.13	8.3%
3	3.69	0.23	6.2%	0.27	7.3%
4	4.35	0.23	5.3%	0.28	6.4%
5	9.11	0.58	6.4%	0.87	9.5%

Linearity (ng/mL)




	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8	4.03	—	—
	4 in 8	1.93	2.02	96%
	2 in 8	0.93	1.01	92%
	1 in 8	0.45	0.50	89%
2	8 in 8	4.33	—	—
	4 in 8	2.15	2.17	99%
	2 in 8	0.99	1.08	91%
	1 in 8	0.50	0.54	92%
3	8 in 8	4.88	—	—
	4 in 8	2.75	2.44	113%
	2 in 8	1.10	1.22	90%
	1 in 8	0.51	0.61	84%
4	8 in 8	6.62	—	—
	4 in 8	3.73	3.31	113%
	2 in 8	1.61	1.66	97%
	1 in 8	0.75	0.83	91%
5	8 in 8	7.70	—	—
	4 in 8	4.18	3.85	109%
	2 in 8	1.97	1.93	102%
	1 in 8	1.02	0.96	106%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH