Vitamin B12

DPC®

IMMULITE[®] 2000 Vitamin B12

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of vitamin B12 in serum or heparinized plasma, as an aid in clinical diagnosis and treatment of anemia.

Catalog Number: **L2KVB2** (200 tests),
L2KVB6 (600 tests)

Test Code: **VB** Color: **Violet**

Summary and Explanation

Vitamin B12 (cobalamin) and folate are nutrients essential to hematopoiesis.⁶ Megaloblastic anemia is almost always due to lack of one of these two vitamins.⁶ Vitamin B12 deficiency can also result in severe neurological impairment.^{1,2,7}

Circulating levels of vitamin B12 are usually a good index to tissue stores. That is, vitamin B12 levels as measured in serum or plasma by an optimized assay system are typically low in vitamin B12 deficiency, and normal or elevated otherwise.¹ Exceptions to this rule can occur in those relatively uncommon situations where levels of vitamin B12 transport proteins are abnormal.² Thus, low circulating vitamin B12 levels can occur in the absence of vitamin B12 deficiency where the level of transcobalamin I (a physiologically inactive transport protein) is low.⁸

Conversely, vitamin B12 deficiency can occur in the presence of normal or even elevated plasma vitamin B12 levels, where transcobalamin II levels are low or where levels of inactive vitamin B12 transport proteins are high, as in chronic myelogenous leukemia.^{5,8,11} (Circulating folate levels are usually normal or elevated in vitamin B12 deficiency, but red cell folate levels are frequently low in this condition.¹⁰)

Vitamin B12 deficiency occurs only rarely as a result of dietary lack of this vitamin.^{3,6} More commonly, it results from impaired absorption, as in partial or total gastrectomy, or in pernicious anemia, a

condition characterized by absence or near absence of intrinsic factor.⁶ Since roughly two thirds of all patients with pernicious anemia have blocking antibodies to intrinsic factor (IFbAb), while IFbAb are only very rarely encountered in other situations, IFbAb determinations represent a useful follow-up test for the differential diagnosis of vitamin B12 deficiency.^{4,8,10} (Circulating intrinsic factor antibodies are present in more than half of all pernicious anemia patients. Increased transport protein levels can occur, for example, in chronic myelogenous leukemia.)

Common causes of high vitamin B12 levels include liver disease, myeloproliferative disease (with chronic myelogenous leukemia as a special case) and the use of multivitamin supplements.

Principle of the Procedure

IMMULITE 2000 Vitamin B12 is a solid-phase, competitive chemiluminescent enzyme immunoassay involving an automated alkaline denaturation procedure.

The IMMULITE 2000 performs a one-cycle, sample treatment of patient serum or plasma with dithiothreitol (DTT) and a sodium hydroxide/potassium cyanide solution (NaOH/KCN) in a reaction tube containing no bead. After a 30-minute incubation, the treated sample is transferred to a second reaction tube containing a vitamin B12-coated polystyrene bead and hog intrinsic factor (HIF). During this subsequent 30-minute incubation, the vitamin B12 released from endogenous binding proteins during sample treatment competes with immobilized vitamin B12 for binding with HIF. Alkaline phosphatase-labeled anti-hog intrinsic factor is then introduced which binds to any HIF that is immobilized on the B12-coated bead in the final 30 minute incubation. Unbound enzyme conjugate is removed by centrifugal wash. Substrate is added and the procedure continues as described for typical immunoassays in the Operator's Manual.

Incubation Cycles: 3 × 30 minutes. 2 test positions per assay: 1 sample treatment cup; 1 immunoreaction cup.

Cycle 1: Alkaline denaturation of endogenous binding proteins.

Cycles 2 & 3: Immunoreaction.

Folic Acid Panel / Vitamin B12: To achieve maximum assay throughput, patient samples for folic acid and vitamin B12 should be run alternately, beginning with a folic acid sample. See Operator's Manual on creating a panel.

Specimen Collection

Note: If sample is to be analyzed for both vitamin B12 and folic acid, the patient must be in a fasting state.¹⁷

Although neither bilirubin nor hemolysis has any clinically significant effects on the vitamin B12 results, the presence of hemoglobin in the sample will increase folic acid values. Additionally, hemolysis may indicate mistreatment of the specimen before receipt by the laboratory, hence the results should be interpreted with caution.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Because EDTA has a significant effect on results, it should not be used as an anticoagulant.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Vitamin B12 has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 75 µL serum or heparinized plasma.

Storage: If not assayed within 8 hours, store at –20°C: stable for 6–8 weeks.¹⁴ Avoid excessive exposure to direct light.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chloramphenicol, at concentrations less than 0.1 g/dL has been added as a preservative. Chloramphenicol is known to cause cancer; this disclosure is required by the state of California.

NaOH/KCN: Solution contains cyanide. Avoid all bodily contact.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Vitamin B12 Bead Pack (L2VB12)

With barcode. 200 beads, coated with vitamin B12 analog. Stable at 2–8°C until the expiration date.

L2KVB2: 1 pack. **L2KVB6:** 3 packs.

Vitamin B12 "A" Reagent Wedge (L2VBA2)

With barcodes. Wedges (labeled "A"), each containing 15 mL dithiothreitol (DTT) in buffer, and 15 mL NaOH/KCN (NH₄CN) solution in buffer, with preservative. Prior to opening, stable at 2–8°C until expiration date. *After opening*, stable at 2–8°C for 30 days.

L2KVB2: 3 wedges. **L2KVB6:** 6 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Vitamin B12 “D” Reagent Wedge (L2VBD2)

With barcodes. Wedges (labeled “D”), each containing 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to murine monoclonal anti-hog intrinsic factor in buffer, and 11.5 mL vitamin B12 binding protein (purified hog intrinsic factor) in buffer, with preservative. Stable at 2–8°C until expiration date.
L2KVB2: 1 wedge. **L2KVB6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Vitamin B12 Adjustors (LVBL, LVBH)

Two vials (Low and High), lyophilized vitamin B12 in a human protein-based matrix, with preservative. Reconstitute each vial with **4.0 mL** sterile distilled water, and mix by gentle inversion. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.
L2KVB2: 2 sets. **L2KVB6:** 3 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Vitamin B12/Folic Acid Sample Diluent (L2FVZ)

For on-board dilution of high samples. 25 mL of concentrated (ready-to-use), vitamin B12/folic acid-free human protein-based matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2FVZ: 3 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

CON6: Tri-level, multi-constituent control

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Both Reagent Wedge A (L2VBA2) and Reagent Wedge D (L2VBD2) must be loaded onto the instrument prior to assay.

Recommended Adjustment Interval:

2 weeks.

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of vitamin B12.

Expected Values

Serum samples from 147 apparently healthy male and female volunteers were processed by the IMMULITE 2000 Vitamin B12 assay. Nonparametric analysis yielded a median value of 432 pg/mL (319 pmol/L) and a central 95% range of 193 – 982 pg/mL (142 – 725 pmol/L).

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

EDTA plasma is not recommended for use.

Blocking antibodies specific for intrinsic factor are present in more than half of all pernicious anemia patients. In order to accurately measure vitamin B12 levels in such samples it is necessary to inactivate these antibodies. While the alkaline hydrolysis has been shown to be effective in inactivating high titers of intrinsic factor blocking antibodies that may be present, the possibility that these antibodies are

not fully inactivated in a rare sample, especially samples with extremely high titers of these antibodies, can not be ruled out. When results obtained are in conflict with the clinical examination, patient medical history and other findings, the sample should be tested for intrinsic factor blocking antibodies. A sample pool selected for a high titer of intrinsic factor blocking antibodies could also be incorporated into the routine quality control of this assay to specifically monitor the inactivation of such antibodies during the pretreatment step.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:
pg/mL \times 0.7378 \rightarrow pmol/L

Reportable Range: 150 – 1,000 pg/mL
(111 – 738 pmol/L)

Analytical Sensitivity: 125 pg/mL
(92 pmol/L)

Intraassay Precision: Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

Interassay Precision: Statistics were calculated for samples assayed in at least

10 different runs. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three vitamin B12 solutions (2,000, 4,500 and 7,500 pg/mL), were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for vitamin B12. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 417 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 4,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: Samples were collected using plain, heparinized, and Becton Dickinson SST vacutainer tubes. EDTA plasma is not recommended for use. By linear regression:

(Heparin) = 0.96 (Serum) + 6 pg/mL
r = 0.995

(SST) = 0.92 (Plain Tubes) + 46 pg/mL
r = 0.988

Means:
448 pg/mL (Heparin)
467 pg/mL (SST)
459 pg/mL (Serum)

Method Comparison: The assay was compared to DPC's IMMULITE Vitamin B12 on 153 samples (Concentration range: approximately 152 to 982 pg/mL. See graph.) By linear regression:

(IML 2000) = 0.91 (IML) – 16.4 pg/mL
r = 0.913

Means:
413 pg/mL (IMMULITE 2000)
471 pg/mL (IMMULITE)

References

- 1) Allen RH. Clinical role and current status of serum cobalamin (Vitamin B12) assays. Ligand Quarterly 1981 Fall;4(3):37-44, 67. 2) Carmel R. Megaloblastic anemia: vitamin B12 and folate. In: Fairbanks VF, editor. Current hematology. New York: John Wiley, 1983: 243-80. 3) Colman

N. The radioisotopic investigation of anemia. *Ligand Quarterly* 1981 Fall;4(3):24-30. 4) Dawson DW. Diagnosis of vitamin B12 deficiency. *Br Med J* 1984;289:938-9. 5) Gräsbeck R. Biochemistry and clinical chemistry of vitamin B12 transport and the related diseases. *Clin Biochem* 1984;17:99-107. 6) Herbert V. The nutritional anemias. *Hosp Pract* 1980 Mar;15(3):65-83, 87-9. 7) Herbert V. Biology of disease: megaloblastic anemias. *Lab Invest* 1985;52:3-19. 8) Kapadia CR, Donaldson RM. Disorders of cobalamin (vitamin B12) absorption and transport. *Ann Rev Med* 1985;36:93-110. 9) Kolhouse JF, et al. Cobalamin analogues are present in human plasma and can mask cobalamin deficiency because current radioisotope dilution assays are not specific for true cobalamin. *New Engl J Med* 1978;299:785-92. 10) Lindenbaum J. Status of laboratory testing in the diagnosis of megaloblastic anemia. *Blood* 1983;61:624-7. 11) Parry TE. The diagnosis of megaloblastic anaemia. *Clin Lab Haematol* 1980;2:89-109. 12) Mollin DL, Ross GIM. Serum vitamin B12 concentrations in leukaemia and in some other haematological conditions. *Brit J Haemat* 1995;1:155-72. 13) Gilbert HS, et al. Serum vitamin B12 content and unsaturated vitamin B12 binding capacity in myeloproliferative disease. *Ann Intern Med* 1969;71:719-29. 14) Tietz NW, editor. *Clinical guide to laboratory tests*. 3rd ed. Philadelphia: WB Saunders, 1995:636. 15) Food and Drug Administration (FDA), Letters from the Bureau of Medical Devices (BMD) to manufacturers of Vitamin B12 radioassay kits. April 1980, January 1981. 16) National Committee for Clinical Laboratory Standards (NCCLS). Guidelines for evaluating a B12 (cobalamin) assay. (Villanova, PA: NCCLS proposed standard PLSA-12, March 1979). Revised as "NCCLS recommendations with respect to vitamin B12 assays" — presented to the Food and Drug Administration (FDA) panel, 9 July 1979. 17) Chanarin I. *The megaloblastic anaemias*. 2nd ed. Oxford: Blackwell Scientific, 1979. 18) Zucker RM, Podell E, Allen RH. Multiple problems with current no-boil assays for serum cobalamin. *Ligand Quarterly* 1981 Fall;4(3):52-8. See further: letters by Ziering S, et al, and by Zucker RM, et al, *ibid*, 60-3, and by Allen RH, *ibid*, 1982 Spring;5(1):48-9. 19) Kubasik NP, Ricotta M, Sine HE. Commercially-supplied binders for plasma cobalamin (vitamin B12) analysis — "purified" intrinsic factor, "cobinamide"-blocked R-protein binder, and non-purified intrinsic factor-R-protein binder — compared to microbiological assay. *Clin Chem* 1980;26:598-600.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828

Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Intraassay Precision (pg/mL)

| | Mean ¹ | SD ² | CV ³ |
|---|-------------------|-----------------|-----------------|
| 1 | 190 | 25 | 13% |
| 2 | 343 | 24 | 7.0% |
| 3 | 744 | 50 | 6.7% |

Interassay Precision (pg/mL)

| | Mean ¹ | SD ² | CV ³ |
|---|-------------------|-----------------|-----------------|
| 1 | 191 | 28 | 15% |
| 2 | 369 | 22 | 6.0% |
| 3 | 698 | 55 | 7.9% |


Specificity

| Compound ¹ | Amount Added ² | % Cross reactivity ³ |
|-----------------------|---------------------------|---------------------------------|
| Cobinamide | 50 µg/mL | ND |

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determinabile. **Method Comparison:** Vitamin B12: Vitamina B12.

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Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH