 IMMULITE[®] 2000

CK-MB

DPC[®]

IMMULITE® 2000 CK-MB

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of creatine kinase isoenzyme MB (CK-MB) in heparinized plasma or serum, as an aid in patient management and the assessment of prognosis of myocardial infarction.

Catalog Number: **L2KMB2** (200 tests),
L2KMB6 (600 tests)

Test Code: **CMB** Color: **Dark Blue**

Summary and Explanation

Creatine kinase (CK) is an enzyme, found primarily in muscle and brain tissue, which exists as three dimeric isoenzymes — CK-MM (CK-3), CK-MB (CK-2), and CK-BB (CK-1) — built from subunits designated M and B. The CK-MB isoenzyme, which has a molecular mass of approximately 87,000 daltons, accounts for 5 to 50% of total CK activity in myocardium. In skeletal muscle, by contrast, it normally accounts for just 1% or less, CK-MM being the dominant form, though the percentage can be as high as 10% in conditions reflecting skeletal muscle injury and regeneration (e.g. severe exercise, muscular dystrophy, polymyositis).²

CK-MB is one of the most important myocardial markers (in spite of not being altogether cardiac-specific), with well-established roles in confirming acute myocardial infarction (AMI) and in monitoring reperfusion during thrombolytic therapy following AMI.²

In AMI, plasma CK-MB typically rises some 3 to 8 hours after the onset of chest pains, peaks within 9 to 30 hours, and returns to baseline levels within 48 to 72 hours.¹⁵ The pattern of serial CK-MB determinations is more informative than a single determination: one CK-MB measurement, even when taken at an appropriate time, cannot definitively confirm or rule out the occurrence of AMI. High levels might reflect skeletal injury rather than myocardial damage. A value within the reference range might be significant if it represents an

increase from the patient's baseline level. (Low baseline levels are sometimes encountered in the elderly.) Accordingly, it has been recommended that CK-MB be measured on admission to the emergency room, and at intervals thereafter: for example, at 3-hour intervals over a 6- to 9-hour period in patients with nonspecific electrocardiogram changes;^{2,5} or at 6- to 8-hour intervals over a 24-hour period and more frequently if thrombolytic therapy has been instituted.²

Thrombolytic therapy, if successful, leads to an "enzyme washout" evidenced by a sharp increase in circulating CK-MB levels as early as 90 minutes after the initiation of therapy.² Accordingly, serial CK-MB levels have been used in this context to assess reperfusion.

Principle of the Procedure

IMMULITE 2000 CK-MB is a chemiluminescent enzyme-labeled immunometric assay, based on ligand-labeled anti-CK-MB monoclonal antibody and separation by anti-ligand-coated solid phase. The bound CK-MB complex — and thus also the photon output, as measured by the luminometer — is proportional to the concentration of CK-MB in the sample.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

An ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Icteric or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 CK-MB has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 µL heparinized plasma or serum.

Storage: 24 hours at 2–8°C, or 30 days at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Note that the antibodies to CK-MB used in the procedure are licensed from Washington University.

CK-MB Bead Pack (L2MB12)

With barcode. 200 beads, coated with anti-ligand. Stable at 2–8°C until expiration date.

L2KMB2: 1 pack. **L2KMB6:** 3 packs.

CK-MB Reagent Wedge (L2MBA2)

With barcode. 11.5 mL ligand-labeled anti-CK-MB murine monoclonal antibody and alkaline phosphatase (bovine calf intestine) conjugated to murine monoclonal anti-CK-BB antibody in buffer. Stable at 2–8°C until expiration date.

L2KMB2: 1 wedge. **L2KMB6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

CK-MB Adjustors (LMBL, LMBH)

Two vials (Low and High) of lyophilized CK-MB in a nonhuman serum matrix, with preservative. At least 30 minutes before use, reconstitute each vial by adding 4.0 mL distilled or deionized water. Mix by *gentle* swirling or inversion. Aliquot and freeze after reconstitution: stable for 2 months (aliquotted) at –20°C.

L2KMB2: 1 set. **L2KMB6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 2 (L2M2Z, L2M2Z4)

For the on-board dilution of patient samples. One vial of concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels **L2M2Z4:** 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

CCCM: Cardiac Marker Control Module containing CK-MB.

Also Required
Distilled or deionized water; test tubes;
controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of CK-MB.

Expected Values

A reference range study was conducted on serum samples from 171 apparently healthy laboratory volunteers. This population yielded a median value of 0.9 ng/mL, an upper 97.5th percentile of 2.6 ng/mL, and an upper 99 percentile of 4.8 ng/mL.

A second study was conducted on serum samples from the following population categories:

- 68 patients, hospitalized for non-cardiac related conditions
- 11 non-AMI emergency department patients with chest pain
- 12 non-trauma emergency department patients, with no evidence of cardiovascular disease.

This population of 91 samples yielded a median value of 2.3 ng/mL, with an upper 97.5th percentile of 5.7 ng/mL and an upper 99 percentile of 7.2 ng/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-

33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on samples collected in heparinized tubes.)

Conversion Factor: ng/mL \times 1 \rightarrow μ g/L

Calibration Range: Up to 500 ng/mL.

Analytical Sensitivity: 0.6 ng/mL.

High-dose Hook Effect: None up to 200,000 ng/mL.

Precision: Samples were repeatedly assayed in duplicate over the course of several days, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three CK-MB solutions (500, 2,000 and 7,000 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for CK-MB, with particularly low crossreactivity to other naturally occurring forms of creatine kinase isoenzymes. (See "Specificity" table.)

Bilirubin: May cause a decrease in values. (See "Bilirubin" table.)

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To determine whether serum can be used in the IMMULITE 2000 CK-MB procedure, blood was collected from 36 laboratory volunteers into heparinized and plain vacutainer tubes. All samples were spiked with CK-MB and assayed by the IMMULITE 2000 CK-MB procedure, with the following results.

(Serum) = 0.96 (Heparin) + 2.9 ng/mL
r = 0.993

Means:

122 ng/mL (Heparin)
120 ng/mL (Serum)

A separate study on 44 endogenous samples from healthy volunteers was also conducted.

(Serum) = 1.02 (Heparin) + 0.03 ng/mL
r = 0.995

Means:

1.44 ng/mL (Heparin)
1.49 ng/mL (Serum)

Method Comparison 1: The assay was compared to DPC's IMMULITE CK-MB (LKMB) on 190 samples. (Concentration range: up to approximately 500 ng/mL. See graph 1.) By linear regression:

(IML 2000) = 1.04 (LKMB) – 0.82 ng/mL
r = 0.998

Means:

39.5 ng/mL (IMMULITE 2000)
38.7 ng/mL (IMMULITE – LKMB)

Method Comparison 2: One hundred twenty-five of these samples, with values less than 10 ng/mL, were analyzed by linear regression separately.

(See graph 2.)

(IML 2000) = 0.88 (LKMB) + 0.009 ng/mL
r = 0.961

Means:

3.12 ng/mL (IMMULITE 2000)
3.53 ng/mL (IMMULITE – LKMB)

References

- 1) Adams JE, Schechtman KB, Landt Y, et al. Comparable detection of acute myocardial infarction by creatine kinase MB isoenzyme and cardiac troponin I. Clin Chem 1994;40:1291-5.
- 2) Apple FS, Preese LM. Creatine kinase-MB: detection of myocardial infarction and monitoring reperfusion. J Clin Immunoassay 1994;17:24-9.
- 3) Bhayana V, Cohoe S, Leung FY, et al. Diagnostic evaluation of creatine kinase-2 mass and creatine kinase-3 and -2 isoform ratios in early diagnosis of acute myocardial infarction. Clin Chem 1993;39:488-95.
- 4) Bruns DE. Diagnosis of acute myocardial infarction when

skeletal muscle damage is present: a caveat regarding use of creatine kinase isoenzymes. Clin Chem 1989;35:705.

- 5) Gibler WB, Lewis LM, Erb RE, et al. Early detection of acute myocardial infarction in patients presenting with chest pain and nondiagnostic ECGs; serial CK-MB sampling in the emergency department. Ann Emerg Med 1990;19:1359-66.
- 6) Henderson AR, Stark JA, McQueen MJ, et al. Is determination of creatine kinase-2 after electrophoretic separation accurate? Clin Chem 1994;40:177-83.
- 7) Kallner A, Sylven C, Brodin U, et al. Early diagnosis of acute myocardial infarction; a comparison between chemical predictors. Scand J Clin Lab Invest 1989;49:633-9.
- 8) Kiyasu Y, John. Current status of detecting CK-MB for patient management. Am Clin Prod Rev 1985;4:29-31.
- 9) Lang H, Wuerzburg U. Creatine kinase, an enzyme of many forms. Clin Chem 1982;28:1439-47.
- 10) Lee KN, Csako G, Bernhardt P, Elin RJ. Relevance of macro creatine kinase type 1 and type 2 isoenzymes to laboratory and clinical data. Clin Chem 1994;40:1278-83.
- 11) Lee TH, Rouan GW, Weisberg MC, et al. Sensitivity of routine clinical criteria for diagnosing myocardial infarction within 24 hours of hospitalization. Ann Intern Med 1987;106:181-6.
- 12) Panteghini M. Creatine kinase MB isoforms. J Clin Immunoassay 1994;17:30-4.
- 13) Thompson WG, Mahr RG, Yohannan WS, Pincus MR. Use of creatine kinase MB isoenzyme for diagnosing myocardial infarction when total creatine kinase activity is high. Clin Chem 1988;34:2208-10.
- 14) Vaidya HC. Creatine kinase-MB. Clin Chem News 1988;14:11-2.
- 15) Wong SS. Strategic utilization of cardiac markers for the diagnosis of acute myocardial infarction. Ann Clin Lab Sci 1996;26:301-12.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	13.9	0.81	5.8%	0.90	6.5%
2	54.3	2.97	5.5%	3.68	6.8%
3	95.4	5.83	6.1%	7.10	7.4%
4	216	10.5	4.9%	12.0	5.6%
5	331	24.2	7.3%	30.5	9.2%

Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	20	—	—
	4 in 8	9.5	10	95%
	2 in 8	5.1	5.0	102%
	1 in 8	2.6	2.5	104%
2	8 in 8	30	—	—
	4 in 8	15	15	100%
	2 in 8	7.7	7.5	103%
	1 in 8	3.9	3.8	103%
3	8 in 8	145	—	—
	4 in 8	67	73	92%
	2 in 8	33	36	92%
	1 in 8	16	18	89%
4	8 in 8	149	—	—
	4 in 8	70	75	93%
	2 in 8	35	37	95%
	1 in 8	18	19	95%

Recovery (ng/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	5.9	—	—
	A	29	31	94%
	B	105	106	99%
	C	341	356	96%
2	—	18	—	—
	A	41	42	98%
	B	118	117	101%
	C	339	367	92%
3	—	54	—	—
	A	75	76	99%
	B	155	151	103%
	C	391	401	98%
4	—	71	—	—
	A	92	92	100%
	B	172	167	103%
	C	420	417	101%
5	—	126	—	—
	A	150	145	103%
	B	236	220	107%
	C	488	470	104%

Specificity



Compound ¹	Amount Added ² (ng/mL)	% Cross reactivity ³
CK-MM	5000	ND
CK-MM	100,000	0.27%
CK-BB	2500	ND
CK-BB	5000	ND

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH