Myoglobin

DPC®

IMMULITE[®] 2000 Myoglobin

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of myoglobin in serum and heparinized plasma, as an aid in the diagnosis of acute myocardial infarction (AMI).

Catalog Number: **L2KMY2** (200 tests),
L2KMY6 (600 tests)

Test Code: **MYO** Color: **Brown**

Summary and Explanation

Acute myocardial infarction (AMI) disrupts cardiac cell membranes, releasing intracellular cardiac proteins into the vascular system. Some of the proteins, including myoglobin, creatine kinase-MB (CK-MB), lactate dehydrogenase type 1 (LD1), and cardiac troponin subunits I and T, have proven useful in diagnosing AMI.^{1,2} The optimal clinical utility of each marker depends on specific protein characteristics. Myoglobin, the smallest of the markers, diffuses rapidly throughout the vascular system and provides the earliest indication of AMI. Myoglobin levels become elevated 0.5–2 hours after chest pain onset and peak within 5–12 hours.¹ The kidneys rapidly eliminate the 17.8 kDa protein from the system, restoring normal circulating concentrations within 16–36 hours.¹ Since the protein rapidly clears from the system, myoglobin concentrations can reliably indicate reinfarction.³ Additionally, myoglobin measurements can preclude AMI: two consecutive low measurements, the first upon patient admission and the second 1–2 hours later, negatively predict AMI in nearly all cases.^{3–5} Myoglobin measurements provide early detection of reperfusion after thrombolytic treatment as well.⁶ Some unrelated patient conditions produce high myoglobin levels, decreasing test specificity. Myoglobin participates in aerobic metabolism in both skeletal and cardiac muscle cells, and high levels accompany various muscle traumas. Renal failure and other kidney problems produce high myoglobin levels as well. Most of the complications have

distinct clinical symptoms that allow reliable differential diagnosis.^{1,2}

Principle of the Procedure

IMMULITE 2000 Myoglobin is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Myoglobin has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 10 µL serum or heparinized plasma.

Storage: 10 days at 2–8°C or 2 months at –20°C.⁸

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis;

for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Myoglobin Bead Pack ((L2MY12)

With barcode. 200 beads, coated with monoclonal murine anti-myoglobin antibody. Stable at 2–8°C until expiration date.

L2KMY2: 1 pack. **L2KMY6:** 3 packs.

Myoglobin Reagent Wedge (L2MYA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-myoglobin antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KMY2: 1 wedge. **L2KMY6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Myoglobin Adjustors (LMYL, LMYH)

Two vials (Low and High) of lyophilized myoglobin in a nonhuman serum matrix, with preservative. At least 30 minutes before use, reconstitute each vial with **2.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KMY2: 1 set. **L2KMY6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 2 (L2M2Z, L2M2Z4)

For the on-board dilution of high samples. One vial of concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2M2Z: 25 mL. **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels. **L2M2Z4:** 5 labels.

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

CCCM: A bi-level, nonhuman serum-based Cardiac Marker Control Module, containing myoglobin as one of several different constituents.

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of myoglobin.

Expected Values

Based on its relationship to DPC's IMMULITE Myoglobin (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

Samples from healthy laboratory volunteers and from hospitalized patients ($n = 258$), that had been shown to yield myoglobin values within the normal reference range of another immunometric method, were analyzed using the IMMULITE Myoglobin assay. The median value for these samples was 25 ng/mL; 97.5% of the values were below 70 ng/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL.

Calibration Range: Up to 1,000 ng/mL

Analytical Sensitivity: 0.5 ng/mL

High-dose Hook Effect:

None up to 20,000 ng/mL.

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: : Samples spiked 1 to 19 with three myoglobin solutions (800, 6,000 and

14,000 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for myoglobin.

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 417 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 24 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of myoglobin, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Myoglobin procedure.

(Heparin) = 1.08 (Serum) – 5.7 ng/mL
 $r = 0.999$

(EDTA) = 1.07 (Serum) – 9.6 ng/mL
 $r = 0.997$

(SST) = 1.04 (Plain Tubes) – 10.7 ng/mL
 $r = 0.345$

Means:

342 ng/mL (Serum)
364 ng/mL (Heparin)
355 ng/mL (EDTA)
345 ng/mL (SST)

Method Comparison 1: The assay was compared to DPC's IMMULITE Myoglobin on 154 samples. (Concentration range: approximately 7 to 950 ng/mL. See graph 1.) By linear regression:

(IML 2000) = 0.96 (IML) + 7.9 ng/mL
 $r = 0.997$

Means:

209 ng/mL (IMMULITE 2000)
209 ng/mL (IMMULITE)

Method Comparison 2: The assay was compared to DPC's IMMULITE *Turbo* Myoglobin on 148 samples. (Concentration range: approximately 7 to 750 ng/mL. See graph 2.) By linear regression:

(IML 2000) = 1.05 (IML Turbo) + 3.6 ng/mL
 r = 0.996

Means:

184 ng/mL (IMMULITE 2000)

172 ng/mL (IMMULITE Turbo)

References

1) Wong SS. Strategic Utilization of cardiac markers for the diagnosis of acute myocardial infarction. *Ann Clin Lab Sci* 1996;26(4):301-12.
 2) Wu AH. Use of cardiac markers as assessed by outcomes analysis. *Clin Biochem* 1997;30(4):339-50. 3) Gornall DA, Roth SN. Serial myoglobin quantitation in the early assessment of myocardial damage: a clinical study. *Clin Biochem* 1996;29(4):379084. 4) Wu AH, Feng YJ, Contois JH, Pervaiz S. Comparison of myoglobin, creatine kinase-MB, and cardiac troponin I for diagnosis of acute myocardial infarction. *Ann Clin Lab Sci* 1996;26(4):291-300. 5) Woo J, Lacbawan FL, Sunheimer R, LeFever D, McCabe JB. Is myoglobin useful in the diagnosis of acute myocardial infarction in the emergency department setting? *Am J Clin Pathol* 1995;103(6):725-9. 6) Lavin F, Kane M, Forde A, Gannon F, Daly K. Comparison of five cardiac markers in the detection of reperfusion after thrombolysis in acute myocardial infarction. *British Heart Journal* 1995;73(5):422-7. 7) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 8) Tietz NW, editor. *Clinical guide to laboratory tests*. 3rd ed. Philadelphia: W.B. Saunders, 1995:442.

Technical Assistance

In the United States, contact DPC's Technical Services department.
 Tel: 800.372.1782 or 973.927.2828
 Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹			Total ²	
		SD ⁴	CV ⁵	SD	CV	
1	24.8	1.2	4.8%	1.4	5.6%	
2	93	3.8	4.1%	5.0	5.4%	
3	210	9.3	4.4%	11.3	5.4%	
4	350	18.6	5.3%	18.9	5.4%	
5	621	38.1	6.1%	43.4	7.0%	

Linearity (ng/mL)



	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	232	—	—
	4 in 8	119	116	103%
	2 in 8	57.9	58.0	100%
	1 in 8	29.6	29.0	102%
2	8 in 8	325	—	—
	4 in 8	176	163	108%
	2 in 8	79.0	81.3	97%
	1 in 8	40.7	40.6	100%
3	8 in 8	517	—	—
	4 in 8	246	259	95%
	2 in 8	125	129	97%
	1 in 8	61.6	64.6	95%
4	8 in 8	673	—	—
	4 in 8	360	337	107%
	2 in 8	175	168	104%
	1 in 8	88.1	84.1	105%
5	8 in 8	967	—	—
	4 in 8	479	484	99%
	2 in 8	225	242	93%
	1 in 8	116	121	96%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH