 IMMULITE[®]
2000

Troponin I

DPC[®]

IMMULITE® 2000 Troponin I

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of troponin I in serum, heparinized or EDTA plasma, as an aid in the diagnosis of acute myocardial infarction (AMI).

Catalog Number: **L2KTI2** (200 tests),
L2KTI6 (600 tests)

Test Code: **TPI** Color: **Light Gray**

Summary and Explanation

Acute myocardial infarction (AMI) is usually diagnosed on the basis of chest pain, electrocardiographic changes, and elevations of markers of myocardial injury. The MB isoenzyme of creatine kinase (CK-MB) has been the preferred marker for two decades.² A study by Wu, et al. found excellent clinical sensitivity of the CK-MB assay between 6 and 24 hours following onset of AMI, with decreased sensitivity beginning in the 24 to 48-hour interval.¹³ However, CK-MB levels can also increase in patients with acute or chronic muscle disease who lack apparent cardiac injury. In the same study, myoglobin, a muscle protein considered to be an early marker of AMI, became elevated within 6 hours of onset, achieved peak clinical sensitivity in the interval 6 to 12 hours following onset, and no longer offered diagnostic value by hour 24.¹³ Myoglobin, although valuable for the early information it provides, also lacks specificity for cardiac injury. A marker specific for myocardial injury is therefore highly desirable.

Cummins, et al.^{6,7} reported the release of cardiac troponin I (cTnI) in AMI. Many studies have focused on cTnI as a candidate marker with acceptable sensitivity and specificity for AMI and other cardiac diseases.

Troponin, a molecule that binds to the thin filament (actin) of striated muscle fibers, acts with intracellular calcium to control the interaction of the thin filament with the thick filament (myosin), thus regulating muscle contraction. Troponin consists of

three subunits: T, which connects the troponin complex and tropomyosin (another cardiac muscle regulatory protein); I, which prevents muscle contraction in the absence of calcium; and C, which binds calcium.⁸ Cardiac troponin I (MW 22.5 kDa) and the two skeletal muscle isoforms of troponin I have considerable amino acid sequence homology, but cTnI contains an additional N-terminal sequence¹¹ and is highly specific for myocardium.¹

Clinical studies report several desirable features of cTnI as a marker of myocardial injury. cTnI rises early in AMI patients and attains levels that are clearly separated from baseline values, so that by 7 hours following onset, the cTnI test detects 95 percent of patients in whom AMI will be confirmed.⁹ Plasma values of cTnI remain elevated for several days, providing a long window for detection of cardiac injury.^{3,13} cTnI has also demonstrated value for predicting mortality risk in unstable angina and in non-Q wave myocardial infarction.⁵

cTnI has demonstrated equivalent diagnostic accuracy for AMI when compared with lactate dehydrogenase type 1 and CK-MB,^{3,10} and may clarify diagnosis in contexts where elevated CK-MB cannot be attributed with certainty to cardiac injury alone.³ These include surgery,⁴ traumatic injury, renal failure, seizures, and skeletal muscle myopathies.¹ In addition, a study on patients undergoing coronary artery bypass grafting (CABG) showed cTnI to be a sensitive marker for perioperative myocardial infarction (PMI); the peak concentration and time of peak both served as diagnostic criteria.¹²

Principle of the Procedure

IMMULITE 2000 Troponin I is a solid-phase, two-site chemiluminescent immunometric assay. The bound complex – and thus also the photon output, as measured by the luminometer – is proportional to the concentration of troponin I in the sample.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

EDTA plasma yields lower assay results. When monitoring, the baseline sample and all subsequent samples should be collected using the same type of tubes, either plain, or EDTA, or heparinized. Do *not* interchange the type of tubes. Please see the Alternate Sample Type section for further information.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Icteric or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Troponin I has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 100 μ L serum or plasma.

Storage: 5 days at 2–8°C,¹⁶ or 1 month at –20°C.¹⁷

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis

B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Troponin I Bead Pack (L2TI12)

With barcode. 200 beads, coated with monoclonal murine anti-troponin I. Stable at 2–8°C until expiration date.

L2KT12: 1 pack. **L2KT16:** 3 packs.

Troponin I Reagent Wedge (L2TIA2)

With barcode. 18.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-troponin I antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KT12: 1 wedge. **L2KT16:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Troponin I Adjustors (LTIL, LTIH)

Two vials (Low and High) of troponin I in a lyophilized nonhuman serum matrix, with preservative. Reconstitute each vial with **3.0 mL** distilled or deionized water. Let stand for 30 minutes, then mix by *gentle* swirling or inversion. Stable for 2 months after reconstitution (aliquotted) at –20°C.

L2KT12: 1 set. **L2KT16:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 1

For the on-board dilution of high samples. One vial of concentrated (ready-to-use) processed, normal human serum, containing undetectable to low levels of troponin I, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2M1Z: 25 mL.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M1Z: 3 labels.

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

CCCM: Bi-level cardiac marker control containing troponin I.

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of troponin I.

Expected Values

Based on its relationship to DPC's IMMULITE Troponin I (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

Two hundred and fifty-five serum samples from healthy laboratory volunteers and from hospitalized patients that had been shown to be negative for troponin I by another immunometric method were analyzed using the IMMULITE Troponin I assay. The median value for these samples was nondetectable; 98% of the values were below 1.0 ng/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

The assay yields lower values when used with EDTA plasma. Please refer to the Alternate Sample Type section for more information.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: Up to 180 ng/mL.

Analytical Sensitivity: 0.2 ng/mL

High Dose Hook Effect: None up to 90,000 ng/mL.

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three troponin I solutions (200, 400 and 900 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for troponin I. (See "Specificity" table.)

Alternate Sample Type: Samples ($n = 12$) were collected using plain and heparinized vacutainer tubes. Matched samples were spiked with troponin I. By linear regression:

(Heparin) = 0.9 (Serum) + 1.84 ng/mL
 $r = 0.980$

Means:
40 ng/mL (Serum)
37 ng/mL (Heparin)

In another study, samples ($n = 13$) were collected using plain and EDTA vacutainer tubes. Matched samples were spiked with troponin I. By linear regression:

(EDTA) = 0.8 (Serum) – 3.1 ng/mL
 $r = 0.988$

Means:
72 ng/mL (Serum)
55 ng/mL (EDTA)

The results show that blood samples collected in EDTA vacutainer tubes yield lower troponin I results than blood samples collected in plain and heparinized tubes.

Bilirubin: Would have an effect on the assay, causing a depression of values. (See "Bilirubin" table.)

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison 1: The assay was compared to DPC's IMMULITE Troponin I on 162 patient samples. (Concentration range: approximately up to 120 ng/mL, as measured by IMMULITE Troponin I. See graph 1.) By linear regression:

(IML 2000) = 0.91 (IML) + 0.73 ng/mL
 $r = 0.996$

Means:
12.4 ng/mL (IMMULITE 2000)
12.9 ng/mL (IMMULITE)

Method Comparison 2: The assay was compared to DPC's IMMULITE Troponin I on 128 samples out of 162 patient samples used for Method Comparison 1. (Concentration range: approximately up to 15 ng/mL, as measured by IMMULITE Troponin I. See graph 2.) By linear regression:

(IML 2000) = 1.15 (IML) – 0.06 ng/mL
 $r = 0.995$

Means:
3.0 ng/mL (IMMULITE 2000)
2.7 ng/mL (IMMULITE)

References

- 1) Adams J, Bodor G, Davila-Roman V, Delmez JA, Apple FS, Ladenson JH, et al. Cardiac troponin I: a marker with high specificity for cardiac injury. *Circulation* 1993;88:101-6.
- 2) Adams JE III, Abendschein DR, Jaffe AS. Biochemical markers of myocardial injury. Is MB creatine kinase the choice for the 1990's? *Circulation* 1993;88:750-63.
- 3) Adams JE III, Schechtman KB, Landt Y, Ladenson JH, Jaffe AS. Comparable detection of acute myocardial infarction by creatine kinase MB isoenzyme and cardiac troponin I. *Clin Chem* 1994;40:1291-5.
- 4) Adams JE, Sicard G, Allan BT, Bridwell KH, Lenke LG, Davila-Roman VG, et al. More accurate diagnosis of perioperative myocardial infarction with measurement of cardiac troponin I. *N Engl J Med* 1994;330:670-4.
- 5) Antman EM, Tanasijevic MJ, Thompson B, Schactman M, McCabe CH, Cannon CP, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. *N Engl J Med* 1996;335:1342-9.
- 6) Cummins B, Cummins P. Cardiac specific troponin I release in canine experimental myocardial infarction: development of a sensitive enzyme-linked immunoassay. *J Mol Cell Cardiol* 1987;19:999-1010.
- 7) Cummins B, Auckland ML, Cummins P. Cardiac-specific troponin I radioimmunoassay in the diagnosis of acute myocardial infarction. *Am Heart J* 1987;113:1333-44.
- 8) Darnell J, Lodish H, Baltimore D. *Molecular cell biology*. New York: Scientific American Books, 1986: 827-8.
- 9) Fonarow GC. UCLA Clinical Practice Guideline 1996 July;2(3) on Internet at http://www.cost-quality.com/2_3art.html.
- 10) Jaffe AS, Landt Y, Parvin CA, Abendschein DR, Geltman EM, Ladenson JH. Comparative sensitivity of cardiac troponin I and lactate dehydrogenase isoenzymes for diagnosing acute myocardial infarction. *Clin Chem* 1996;42:1770-6.
- 11) Larue C, Defacque-Lacquement H, Calzolari C, Nguyen DL, Pau B. New monoclonal antibodies as probes for human cardiac troponin I: epitopic analysis with synthetic peptides. *Molec Immunol* 1992;29:271-8.
- 12) Mair J, Larue C, Mair P,

Balogh D, Calzolari C, Puschendorf B. Use of cardiac troponin I to diagnose perioperative myocardial infarction in coronary artery bypass grafting. Clin Chem 1994;40:2066-70. 13) Wu AHB, Fen YJ, Contois JH, Pervaiz S. Comparison of myoglobin, creatine kinase-MB, and cardiac troponin I for diagnosis of acute myocardial infarction. Ann Clin Lab Sci 1996;26:291-300. 14) Peetz D, Hafner G. Institute of Clinical and Laboratory Medicine, University of Mainz. Unpublished results. 15) Wong SS. Strategic utilization of cardiac markers for the diagnosis of acute myocardial infarction. Ann Clin Lab Sci 1996;26(4):301-12. 16) Data on file. 17) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: W.B. Saunders, 1995:614.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.
The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	0.50	0.05	10%	0.08	16%
2	1.33	0.05	3.8%	0.18	13.5%
3	7.6	0.27	3.6%	0.37	4.9%
4	16.2	0.53	3.3%	1.23	7.6%
5	164	6.45	3.9%	15.2	9.3%

Specificity (ng/mL)

Compound ¹	Amount Added ² ng/mL	% Cross reactivity ³
Human Skeletal Troponin I	10	ND
Human Skeletal Troponin I	100	0.2%
Human Skeletal Troponin I	1,000	0.1%
Human Skeletal Troponin T	10	ND
Human Skeletal Troponin T	100	ND
Human Skeletal Troponin T	1,000	ND
Human Cardiac Troponin T	10	ND
Human Cardiac Troponin T	100	ND
Human Cardiac Troponin T	1,000	0.1%

ND: Not detectable.⁴

Linearity (ng/mL)



	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	5.0	—	—
	4 in 8	2.4	2.5	96%
	2 in 8	1.2	1.3	92%
	1 in 8	0.63	0.63	100%
2	8 in 8	11	—	—
	4 in 8	5.6	5.5	102%
	2 in 8	2.8	2.8	100%
	1 in 8	1.6	1.4	114%
3	8 in 8	20	—	—
	4 in 8	10	10.0	100%
	2 in 8	4.5	5.0	90%
	1 in 8	2.3	2.5	92%
4	8 in 8	51	—	—
	4 in 8	24	26	92%
	2 in 8	13	13	100%
	1 in 8	6.0	6.4	94%
5	8 in 8	60	—	—
	4 in 8	28	30	93%
	2 in 8	13	15	87%
	1 in 8	6.2	7.5	83%

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH