 IMMULITE<sup>®</sup>  
2000

**HCG**

**DPC<sup>®</sup>**

# IMMULITE<sup>®</sup> 2000 HCG

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of human chorionic gonadotropin (HCG) in serum, and for strictly qualitative determinations in urine, as an aid in the detection of pregnancy.

Catalog Number: **L2KCG2** (200 tests), **L2KCG6** (600 tests)

Test Code: **HCG** Color: **Orange**

CDC Analyte Identifier Code: 2502

CDC Test System Identifier Code: 10418

CLIA Complexity Category: Moderate

### Summary and Explanation

Human chorionic gonadotropin (HCG) is a two-chain glycoprotein hormone (MW ~37,000) normally found in blood and urine only during pregnancy. It is secreted by placental tissue, beginning with the primitive trophoblast, almost from the time of implantation, and serves to support the corpus luteum during the early weeks of pregnancy.

According to the literature, circulating HCG typically reaches levels of approximately 2,000 mIU/mL one month after conception. A peak level on the order of 100,000 mIU/mL is attained in the third month, after which a gradual decline sets in. Following delivery, the HCG level normally undergoes rapid descent, reaching nonpregnant concentrations (usually less than 5 mIU/mL) some two weeks later.

Ectopic pregnancies and pregnancies terminating in spontaneous abortion tend to have lower than normal circulating HCG levels, while somewhat higher levels are often seen in multiple pregnancies.

### Principle of the Procedure

IMMULITE 2000 HCG is a solid-phase, two-site chemiluminescent immunometric assay.

**Incubation Cycles:** 1 × 30 minutes.

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 HCG has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 5 µL serum or urine.

**Urine:** A freshly voided, first-morning specimen is preferred.<sup>22</sup> Centrifuging or filtering (Whatman #1) is recommended to clear cloudy samples.

#### Storage

**Serum:** 7 days at 2–8°C, or 2 months at –20°C.<sup>18</sup>

**Urine:** 48 hours at 2–8°C,<sup>22</sup> or at –20°C until assayed.<sup>13</sup>

### Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

### HCG Bead Pack (L2CG12)

With barcode. 200 beads, coated with monoclonal murine anti-HCG. Stable at 2–8°C until expiration date. **L2KCG2:** 1 pack. **L2KCG6:** 3 packs.

### HCG Reagent Wedge (L2CGA2)

With barcode. 11.5 mL of alkaline phosphatase (bovine calf intestine) conjugated to polyclonal ovine anti-HCG, in buffer, with preservative. Stable at 2–8°C until expiration date.

**L2KCG2:** 1 wedge. **L2KCG6:** 3 wedges.

Before use, tear off the top of the label at the perforation, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### HCG Adjustors (LCGL, LCGH)

Two vials (Low and High) of lyophilized HCG in an HCG-free human serum matrix, with preservative. *At least 30 minutes before use*, reconstitute each vial with **3.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

**L2KCG2:** 1 set. **L2KCG6:** 2 sets.

Before running adjustors, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

## Kit Components Supplied Separately

### HCG Sample Diluent (L2CGZ)

For the on-board dilution of high samples. 50 mL of concentrated (ready-to-use) HCG-free human serum, with barcode labels. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2CGZ:** 5 labels

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**L2ZT:** 250 Sample Diluent Test Tubes (16 × 100 mm)

**L2ZC:** 250 Sample Diluent Tube Caps

**CON6:** Tri-level, multi-constituent control

**LCGCM:** HCG High Level Control Module

Also Required

Distilled or deionized water; test tubes; controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

**Recommended Adjustment Interval:** 4 weeks.

**Quality Control Samples:** Use controls or serum pools with at least two levels (low and high) of HCG.

## Qualitative Procedure for Urine

The qualitative procedure allows for determining whether a urine sample has an HCG concentration above or below 30 mIU/mL, within the precision of the assay. Concentrations above this cutoff may be considered positive for pregnancy providing other causes of high HCG levels have been ruled out, viz. trophoblastic and nontrophoblastic neoplasia.<sup>4,9,10,15</sup>

Concentrations below this cutoff are to be expected for samples from individuals who are not pregnant, as well as for women in

such an early stage of pregnancy that the HCG concentration has not yet reached the decision level set for the procedure. Assaying another sample collected two or more days later should settle the matter, since in early pregnancy HCG levels can be expected to double in approximately this time.<sup>11</sup>

### Interpretation of Results

A result of greater than or equal to 30 mIU/mL indicates that the patient sample is **positive** for pregnancy.

A result of less than 30 mIU/mL indicates that the patient sample is **negative** for pregnancy.

### Expected Values

Based on its relationship to DPC's IMMULITE HCG (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

#### Males and Nonpregnant Females

In a study performed with the IMMULITE HCG assay and involving serum samples from 797 adult volunteers, including 428 males and 369 nonpregnant females, 95% of the results for males were below that assay's detection limit (1.1 mIU/mL) and all were below 2.5 mIU/mL, while 83% of the results for nonpregnant females were below the detection limit, 95% were below 2.7 mIU/mL, and all were below 5.3 mIU/mL. The subjects were in apparent good health based on a questionnaire.

#### Pregnant Females

A total of 593 serum samples from apparently healthy pregnant women were processed by the IMMULITE HCG assay. The results are summarized below (in mIU/mL) by gestational age (Gest.) — equivalently, by weeks since last menstrual period (LMP). Also see "Expected Values During Pregnancy" graph.

### Weeks of Pregnancy

Gest.	LMP	n	Median	Central 95%
1.3–2	3.3–4	30	71	16–156
2–3	4–5	54	607	101–4,870
3–4	5–6	34	5,243	1,110–31,500
4–5	6–7	34	26,983	2,560–82,300
5–6	7–8	36	52,090	23,100–151,000
6–7	8–9	33	93,598	27,300–233,000
7–11	9–13	116	117,678	20,900–291,000
11–16	13–18	72	40,989	6,140–103,000
16–21	18–23	80	20,868	4,720–80,100
21–39	23–41	104	15,352	2,700–78,100

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

### Pathological Conditions

Ectopic pregnancies or abnormal intrauterine pregnancies have circulating HCG levels below normal for that gestational period; after an initial rise these levels plateau or fall.<sup>3,17,18</sup>

### Limitations

The procedure for urine samples is designed solely as a screen for pregnancy, using a 30 mIU/mL cutoff. Hence urinary results should be reported in a strictly qualitative manner, e.g. as "positive" or "negative" for pregnancy — concentration values should *not* be reported for this sample type.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See tables and graphs for data *representative* of the assay's performance. Results are expressed as mIU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Calibration Range:** Up to 5,000 mIU/mL, in terms of the WHO 3rd IS 75/537 for quantitative determinations in serum.

**Analytical Sensitivity:**  
0.4 mIU/mL.

**High-Dose Hook Effect:** None up to 2,000,000 mIU/mL.

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three solutions (1,313, 5,720 and 30,340 mIU/mL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The assay is highly specific for HCG, with low crossreactivity to other glycoprotein hormones present in patient samples. (See "Specificity" table.)

For concentrations as high as 10,000 ng/mL, the HCG  $\alpha$ -subunit (WHO 75/569) showed less than 1% crossreactivity. HCG  $\beta$ -subunit (WHO 75/551), spiked into HCG Sample Diluent, yielded the tabulated results. (See "Specificity" table.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Presence of hemoglobin in concentrations up to 384 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of lipemia in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from 14 volunteers into

plain, heparinized, EDTA vacutainer tubes, and SST. Equal volumes of the matched samples were spiked with various concentrations of HCG, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 HCG procedure.

(Heparin) = 1.06 (Serum) – 7.06 mIU/mL  
 $r = 0.99$

(EDTA) = 1.08 (Serum) – 6.7 mIU/mL  
 $r = 0.99$

(SST) = 0.98 (PlainTubes) – 3.25 mIU/mL  
 $r = 0.99$

Means:

204 mIU/mL (Serum)  
211 mIU/mL (Heparin)  
197 mIU/mL (EDTA)  
214 mIU/mL (SST)

**Method Comparison (Serum):** The assay was compared to DPC's IMMULITE HCG on 115 patient samples.

(Concentration range: approximately 1 to over 3,000 mIU/mL. See graph.) By linear regression:

(IML 2000) = 1.0 (IML) – 21 mIU/mL  
 $r = 0.994$

Means:

668 mIU/mL (IMMULITE 2000)  
689 mIU/mL (IMMULITE)

For the 63 data pairs (out of 116) with results by both methods from 1 to 500 mIU/mL:

(IML 2000) = 0.91 (IML) + 0.9 mIU/mL  
 $r = 0.987$

For the 22 data pairs (out of 116) with results by both methods from 1 to 50 mIU/mL:

(IML 2000) = 0.96 (IML) + 0.3 mIU/mL  
 $r = 0.953$

**Method Comparison (Urine):** A total of 123 urine samples from pregnant women were processed by the IMMULITE 2000 HCG assay and by DPC's Double Antibody HCG. The samples, some of which were diluted, had HCG values ranging up to approximately 300 mIU/mL. By linear regression:

IML 2000 = 1.04 (DAb) – 10 mIU/mL  
 $r = 0.988$

In another study with the IMMULITE 2000 HCG assay, as a check on the appropriateness of the 30 mIU/mL cutoff used in the qualitative procedure, urinary results for nonpregnant women showed a

clean separation from the distribution of repeated measurements made on a 30 mIU/mL sample. (All results for the nonpregnant women: less than 3 mIU/mL. Mean – 6SD for the 30 mIU/mL sample: 23 mIU/mL.)

## References

1) Bangham DR, Storing PL. Standardisation of human chorionic gonadotropin, HCG subunits, and pregnancy tests. *Lancet* 1982;1:390. 2) Batzer F. Hormonal evaluation of early pregnancy. *Fertil Steril* 1980;34:1–12. 3) Braunstein GD, et al. First-trimester chorionic gonadotropin measurements as an aid in the diagnosis of early pregnancy disorders. *Amer J Obstet Gynecol* 1978;131:25–32. 4) Clayton L, et al. Endocrine aspects of trophoblastic neoplasia. *J Reprod Med* 1981;26:192–9. 5) Corson SL, Horwitz CA, Lau HL, Soderstrom R. Early diagnosis of pregnancy: a symposium. *J Reprod Med* 1981;26:149–78. 6) DeCherney AH, Minkin MJ, Spangler S. Contemporary management of ectopic pregnancy. *J Reprod Med* 1981;26:519–23. 7) Goldstein DP, Berkowitz RS, Kosasa TS. Measurement of human chorionic gonadotropin. Major Problems in Obstet Gynecol 1982;14:64–97. 8) Hager H, Van Weemen BK. Standardisation of human chorionic gonadotropin, HCG subunits, and pregnancy tests. *Lancet* 1982;1:629. 9) Kohorn EI. Hydatidiform mole and gestational trophoblastic disease in southern Connecticut. *Obstet and Gynecol* 1982;59:78–84. *Erratum:* The first diagram belongs with the caption to Figure 3, and vice versa. 10) Kosasa TS. Measurement of human chorionic gonadotropin. *J Reprod Med* 1981;26:201–6. 11) Lenton E, Neal L, Sulaiman R. Plasma concentrations of human gonadotropin from the time of implantation until the second week of pregnancy. *Fertil Steril* 1982;37:773–8. 12) Lindstedt G, Janson PO, Thorburn J. Sensitivity of serum chorionic gonadotropin assay for ectopic pregnancy. *Lancet* 1981;1:781–2. 13) McCready J, Braunstein GD, et al. Modification of the choriogonadotropin beta-subunit radioimmunoassay for determination of urinary choriogonadotropin. *Clin Chem* 1978;24:1958–61. 14) Storing PL, Gaines-Das R, Bangham DR. International Reference Preparation of Human Chorionic Gonadotrophin for Immunoassay: potency estimates in various bioassay and protein binding assay systems. *J Endocrinol* 1980;84:295–310. 15) Vaitukaitis J. Human chorionic gonadotropin – a hormone secreted for many reasons. *N Eng J Med* 1979;301:324–6. 16) Danzer H, Braunstein GD, et al. Maternal serum human chorionic gonadotropin concentrations and fetal sex predictions. *Fertil Steril* 1980;34:336–40. 17) Husa RO. *The Clinical Marker HCG*. New York: Praeger Publishers, 1987. 18) Greene MF, de Montserrat MF, Tulchinsky D. In: Tietz NW, editor. *Textbook of Clinical Chemistry*.

Philadelphia: Saunders, 1986: 1745. 19) Braunstein GD, et al. Serum human chorionic gonadotropin levels throughout normal pregnancy. *Am J Obstet Gynecol* 1976;126:678. 20) Babson AL. The IMMULITE Automated Immunoassay System. *J Clin Immunoassay* 1991;14:83–8. 21) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 22) Tietz NW, ed. *Clinical Guide to Laboratory Tests*. 2nd Ed. Philadelphia: W.B. Saunders, 1990.

## Technical Assistance

In the United States, contact DPC's Technical Services department.  
Tel: 800.372.1782 or 973.927.2828.  
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

## Tables and Graphs

### Precision (mIU/mL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>			Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV	
1	6.5	0.43	6.6%	0.48	7.4%	
2	103	3.3	3.2%	4.6	4.5%	
3	239	7.6	3.2%	11.1	4.6%	
4	438	13	3.0%	21	4.8%	
5	3,120	79	2.5%	150	4.8%	

### Linearity (mIU/mL)

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	16 in 16 <sup>5</sup>	204	—	—
	8 in 16	100	102	98%
	4 in 16	50	51	98%
	2 in 16	24	26	92%
	1 in 16	12	13	92%
2	16 in 16	552	—	—
	8 in 16	273	276	99%
	4 in 16	134	138	97%
	2 in 16	68	69	99%
	1 in 16	33	35	94%
3	16 in 16	1,144	—	—
	8 in 16	556	572	97%
	4 in 16	267	286	93%
	2 in 16	140	143	98%
	1 in 16	70	72	97%
4	16 in 16	2,088	—	—
	8 in 16	1,020	1,044	98%
	4 in 16	503	522	96%
	2 in 16	240	261	92%
	1 in 16	119	131	91%

### Recovery (mIU/mL)

	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	270	—	—
	A	339	322	105%
	B	597	543	110%
	C	1,770	1,774	100%
2	—	554	—	—
	A	617	592	104%
	B	913	812	112%
	C	2,114	2,043	103%
3	—	949	—	—
	A	984	967	102%
	B	1,290	1,188	109%
	C	2,610	2,419	108%

### Specificity

	Amount Added (ng/mL) <sup>1</sup>	Apparent Concentration (ng/mL) <sup>2</sup>	% Cross-reactivity <sup>3</sup>
FSH	26.8	ND	ND
LH	16.5	ND	ND
TSH	860	ND	ND


ND: not detectable.<sup>4</sup>


β-HCG WHO 75/551 ng/mL	Apparent HCG Concentration			Cross-reactivity
	HCG mIU/mL	HCG ng/mL*	β-HCG ng/mL†	
42.2	825	88.8	53.7	127%
84.4	1,579	170	103	122%
169	2,789	300	182	108%
338	4,803	517	313	93%

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH