 IMMULITE[®]
2000

Free Beta HCG

DPC[®]

IMMULITE[®] 2000 Free Beta HCG

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of free β -HCG subunit in serum.

Catalog Numbers: **L2KFB2** (200 tests), **L2KFB6** (600 tests)

Test Code: **FBC** Color: **Violet**

CDC Analyte Identifier Code: 2502
CDC Test System Identifier Code: 10418

Summary and Explanation

Human chorionic gonadotropin (HCG) is a glycoprotein hormone normally found in blood and urine only during pregnancy. It is secreted by placental tissue, beginning with the primitive trophoblast, almost from the time of implantation, and serves to support the corpus luteum during the early weeks of pregnancy. HCG or HCG-like material is also produced by a variety of trophoblastic and nontrophoblastic neoplasia.^{1,3,4,10}

Intact HCG is a 39,500 dalton molecule composed of two nonidentical subunits, α and β , that are bound to each other in a noncovalent manner. These subunits can also occur in a free, or unbound, form.^{6,8} Only intact HCG has biological activity. The HCG α subunit is structurally identical to the α subunit of the homologous pituitary glycoprotein hormones, luteinizing hormone, follicle stimulating hormone and thyroid stimulating hormone. The β subunit is specific for each of these hormones, and confers upon them their differing biological activities. The β subunits of HCG and luteinizing hormone, however, are structurally quite similar, accounting for the essentially identical biological activities of these two hormones.⁶

Measurements of free β -HCG levels have been found to aid in the diagnosis and monitoring of trophoblastic diseases (moles, choriocarcinoma) and certain testicular tumors, where the ratios of the free β subunit to intact HCG can be quite high. Occasionally, tumors are found to secrete only free β -HCG subunit, and

virtually no intact HCG can be detected.^{5,7,9} Application in diagnosis and monitoring of tumor recurrence requires high sensitivity in the assay's low end, and careful interpretation is required (see Limitation section). Measurements of free β -HCG levels may also have value in the early monitoring of the outcome of *in vitro* fertilization procedures.²

Maternal serum free β -HCG assessment is reported to have significant utility in first- and second-trimester prenatal screening for Down syndrome and other chromosomal anomalies. A first-trimester combination of maternal age, serum PAPP-A and fetal nuchal translucency (NT) measurements might substantially increase the efficiency of prenatal screening compared to second-trimester screening; using this approach, various investigators have reported detection rates for Down syndrome of 85 to 90 percent at a 5 percent false-positive rate.^{12,13}

Principle of the Procedure

IMMULITE 2000 Free Beta HCG is a solid-phase, two-site sequential chemiluminescent immunometric assay.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE 2000 Free Beta HCG has not been tested with all possible variations of tube types.

Volume Required: 5 μ L serum.

Storage: 7 days at 2–8°C or longer at –20°C. Avoid repeated freeze/thaw cycles.

Dilution of High Samples: Samples expected to have free β -HCG levels above the assay's calibration range should be diluted with HCG Sample Diluent (L2CGZ) before assay.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Free Beta HCG Bead Pack (L2FB12)

With barcode. 200 beads, coated with monoclonal murine anti-free β -HCG antibody. Stable at 2–8°C until expiration date.

L2KFB2: 1 pack. **L2KFB6:** 3 packs.

Free Beta HCG Reagent Wedge (L2FBA2)

With barcode. 11.5 mL of protein buffer matrix. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-free β -HCG antibody

in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KFB2: 1 wedge. **L2KFB6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Free Beta HCG Adjustors (LFBL, LFBH)

Two vials (Low and High) of lyophilized free β -HCG in a human serum-based matrix, with preservative. Reconstitute each vial with 3.0 mL distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KFB2: 1 set. **L2KFB6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

HCG Sample Diluent (L2CGZ)

For the on-board dilution of high samples. 50 mL of concentrated (ready-to-use) HCG-free human serum. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 x 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2CGZ: 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 x 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LFBCM: Tri-level, human serum based control module containing free β -HCG

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of free beta HCG.

Expected Values

The reference range for healthy males, nonpregnant females and post-menopausal females can be expected to be less than 0.1 ng/mL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Notice that the analytical sensitivity of the assay is 0.1 ng/mL and it is generally recognized that a free beta HCG assay should have a lower detection limit of < 0.1 ng/mL for an application in monitoring and detection of trophoblastic and testicular tumors. Thus, appropriate interpretation is needed when the assay is used for this application.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL, where 1 ng/mL of free β -HCG is equivalent to 1 mIU/mL of free β -HCG in terms of the World Health Organization's First International Reference Preparation of Chorionic Gonadotrophin Beta Subunit Human, number 75/551 (1st IRP 75/551).

Unless otherwise noted, all results were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.

Calibration Range: Up to 200 ng/mL

Analytical Sensitivity: 0.1 ng/mL

High-dose Hook Effect: None up to at least 190,000 ng/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table.)

Recovery: Samples spiked 1 to 19 with three free β -HCG solutions (70, 700 and 2,000 ng/mL) were assayed. (See "Recovery" table.)

Specificity: The antibody is highly specific for the free β subunit of HCG. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 540 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to a commercially available Free Beta HCG assay (Kit A) on 145 samples. (Concentration range: approximately 1 to 150 ng/mL. See graph.) By linear regression:

$$(IML\ 2000) = 0.97 (\text{Kit A}) + 2.7\ \text{ng/mL}$$

$$r = 0.991$$

Means:
48.1 ng/mL (IMMULITE 2000)
47.0 ng/mL (Kit A)

References

1) Clayton L, Tyrey L, Weed JC Jr, Hammond CB. Endocrine aspects of trophoblastic neoplasia. *J Reprod Med* 1981;26:192-9. 2) Cole LA, Restrepo-Candelo H, Lavy G, Decherney A. HCG free β -subunit an early marker of outcome of *in vitro* fertilization clinical pregnancies. *J Clin Endocrinol Metab* 1987;64:1328-30. 3) Kohorn EI. Hydatidiform mole and gestational trophoblastic disease in southern Connecticut. *Obstet Gynecol* 1982;59:78-84. Erratum: the first diagram belongs with the caption to Figure 3, and vice versa. 4) Kosasa TS. Measurement of human chorionic gonadotropin. *J Reprod Med* 1981;26:201-6. 5) Light PA, Forster JP, Felton T, Eckert H, Toverly KC. Molecular heterogeneity in chorionic gonadotropin in some testicular cancer patients. *Lancet* 1983;1:1284. 6) Norman RJ, Buck RH, De Medeiros SF. Measurement of human chorionic gonadotropin (hCG): indications and techniques for the clinical laboratory. *Ann Clin Biochem* 1990;27:183-94. 7) Ozturk M, Berkowitz R, Goldstein D, Bellet D, Wands JR. Differential production of human chorionic gonadotropin and free subunits in gestational trophoblastic disease. *Am J Obstet Gynecol* 1988;158:193-8. 8) Pierce JG, Parsons TF. Glycoprotein hormones: structure and function [review]. *Ann Rev Biochem* 1981;50:465-95. 9) Saller B, Clara R, Spottl G, Siddle K, Mann K. Testicular cancer secretes intact human chorionic gonadotropin (hCG) and its free β -subunit: evidence that hCG (+hCG- β) assays are the most reliable in diagnosis and follow-up. *Clin Chem* 1990;36:234-9. 10) Vaitukaitis JL. Human chorionic gonadotropin – a hormone secreted for many reasons. *N Eng J Med* 1979;301:324-6. 11) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 12) Bindra R, *et al.* One-stop clinic for assessment of risk for trisomy 21 at 11-14 weeks: a prospective study of 15,030 pregnancies. *Ultrasound Obstet Gynecol* 2002;3:219-25. 13) Wald NJ, Hackshaw AK. Advances in antenatal screening for Down syndrome. *Baillieres Best Pract Res Clin Obstet Gynaecol* 2000;14(4):563-80.

Technical Assistance

Available outside the United States only.
For technical assistance, contact your
National Distributor.

Manufactured by EURO/DPC Ltd. under a
Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹			Total ²	
		SD ⁴	CV ⁵	SD	CV	
1	0.91	0.04	4.4%	0.06	6.6%	
2	2.66	0.08	3.0%	0.11	4.1%	
3	17.0	0.34	2.0%	0.55	3.2%	
4	27.8	0.58	2.1%	0.86	3.1%	
5	60.2	1.05	1.7%	2.10	3.5%	
6	111	3.13	2.8%	4.37	3.9%	
7	164	4.62	2.8%	8.16	5.0%	

Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	16 in 16 ⁵	14.1	—	—
	8 in 16	7.4	7.1	104%
	4 in 16	3.6	3.5	103%
	2 in 16	1.7	1.8	94%
	1 in 16	0.88	0.88	100%
2	16 in 16	17.3	—	—
	8 in 16	8.4	8.7	97%
	4 in 16	4.1	4.3	95%
	2 in 16	2.1	2.2	95%
	1 in 16	1.0	1.1	91%
3	16 in 16	55.7	—	—
	8 in 16	29.8	27.9	107%
	4 in 16	15.0	13.9	108%
	2 in 16	7.2	7.0	103%
	1 in 16	3.7	3.5	106%
4	16 in 16	100	—	—
	8 in 16	50.6	50.0	101%
	4 in 16	25.9	25.0	104%
	2 in 16	12.5	12.5	100%
	1 in 16	5.9	6.3	94%

Recovery (ng/mL)


	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	ND	—	—
	A	3.7	3.5	106%
	B	36.0	35.0	103%
	C	107	100	107%
2	—	3.3	—	—
	A	7.0	6.6	106%
	B	37.5	38.1	98%
	C	105	103	102%
3	—	7.1	—	—
	A	10.6	10.2	104%
	B	42.7	41.7	102%
	C	111	107	104%
4	—	11.5	—	—
	A	15.0	14.4	104%
	B	46.4	45.9	101%
	C	112	111	101%
5	—	29.3	—	—
	A	32.3	31.3	103%
	B	64.6	62.8	103%
	C	129	128	101%

Specificity

Compound ¹	ng/mL Added ²	Apparent ng/mL ³	% Cross-reactivity ⁴
HCG α subunit	3,500	0.52	ND
Intact HCG	10,200	0.24	ND
FSH	100,000	ND	ND
LH	100,000	0.17	ND
TSH	150,000	ND	ND

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH