 IMMULITE[®]
2000

FSH

DPC[®]

IMMULITE® 2000 FSH

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of FSH in serum, as an aid in the diagnosis and treatment of pituitary and gonadal disorders.

Catalog Number: **L2KFS2** (200 tests),
L2KFS6 (600 tests)
Test Code: **FSH** Color: **Light Gray**

CDC Analyte Identifier Code: 1908
CDC Test System Identifier Code: 10418
CLIA Complexity Category: Moderate

Summary and Explanation

Follicle stimulating hormone (folliotropin, FSH) is secreted by the β -cells of the anterior pituitary under the control of the gonadotropin releasing hormone produced in the hypothalamus. FSH facilitates the development and maintenance of gonadal tissues, which synthesize and secrete steroid hormones. Circulating levels of FSH are controlled by a negative feedback mechanism on the hypothalamus by steroidal hormones. Although FSH and LH are required for normal sexual function in both males and females, the secretory patterns are very different for the two sexes.

In mature females, FSH initiates the growth and development of ovarian follicles. During ovulation, when the follicle is ruptured, the follicle, now called the corpus luteum, secretes estradiol and progesterone, which control the circulating levels of FSH by a negative feedback effect on the hypothalamus. In menopause, with diminished ovarian function, there is a resulting decrease in estradiol secretion. Due to the lack of a negative feedback effect, with diminished estradiol, the circulating FSH levels become significantly increased.

In the mature male, FSH is associated with the stimulation and maintenance of spermatogenesis. Testosterone and estradiol have the role of providing the negative feedback signal to the hypothalamus for controlling the release of

FSH. Infertility in males may be due to hypogonadism as a result of primary testicular failure. Testicular failure may be functional failure to mature, or a result of germ cell damage. Whatever the etiology, the conditions of hypogonadism have the net result of dramatically raising the circulating FSH levels, due to the lack of a negative feedback effect.

Principle of the Procedure

IMMULITE 2000 FSH is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

EDTA plasma should not be used as a sample type.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 FSH has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 μ L serum.

Storage: 7 days at 2–8°C, or 2 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

FSH Bead Pack (L2FS12)

With barcode. 200 beads, coated with monoclonal murine anti-FSH. Stable at 2–8°C until expiration date.

L2KFS2: 1 pack. **L2KFS6:** 3 packs.

FSH Reagent Wedge (L2FSA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to murine monoclonal anti-FSH antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KFS2: 1 wedge. **L2KFS6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

FSH Adjustors (LFSL, LFSH)

Two vials (Low and High), 3.0 mL each, of FSH in a nonhuman serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C

L2KFS2: 1 set. **L2KFS6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

FSH Sample Diluent (L2FSZ)

For the on-board dilution of high samples. 25 mL of concentrated (ready-to-use) FSH-free nonhuman serum matrix, with preservative. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C. Dispose of in accordance with applicable laws.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2FSZ: 3 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

CON6: Tri-level, multi-constituent control

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 4 weeks.

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of FSH.

Expected Values

Based on its relationship to DPC's IMMULITE FSH (mono/poly) assay (see Method Comparison 2), the assay can be expected to have essentially the same reference ranges.

Reference ranges were generated using IMMULITE FSH (mono/poly) in a multi-national study involving women in apparent good health (age: 16–44 years), who volunteered to have blood samples drawn, on a daily basis, throughout one complete ovulatory cycle. (Also see "Menstrual Cycle" Graph in Tables and Graphs section.)

Ovulatory Cycles	n*	FSH, mIU/mL	
		Median	Central 95%
Follicular Phase	54 (762)	6.2	2.8 – 11.3
Follicular Phase, Days 2 to 3	54 (108)	6.6	3.0 – 14.4
Midcycle	54 (54)	13.6	5.8 – 21
Luteal Phase	54 (604)	3.4	1.2 – 9.0

*Number of subjects (total number of results)

The table below summarizes results for adult males and (untreated) postmenopausal females from studies with IMMULITE FSH (mono/poly). The ranges for women on oral contraceptives and for postmenopausal women on ERT are based on the IMMULITE FSH assay's relationship to DPC's Coat-A-Count FSH IRMA.

Group	n	FSH, mIU/mL	
		Median	Central 95%
Adult Males	135	3.8	0.7 – 11.1
Adult Females:			
Postmenopausal*	76	90.5	21.7 – 153
Postmenopausal (ERT)	16	27	9.7 – 111
Oral Contraceptives	12	1.7	ND – 4.9

*Preliminary ND: not detectable

A cross-sectional study of pediatric fertility values performed with IMMULITE FSH (mono/poly) at a "wellness" clinic in the southwestern United States yielded the following results.

Group	Age (yr)	n	FSH, mIU/mL	
			Median	Central 95%
Females	Cord	30	ND	
	0.1 – 3	57	2.3	0.11 – 13
	4 – 9	28	0.8	0.11 – 1.6
Males	Cord	37	0.24	ND – 1.2
	0.1 – 3	72	0.6	ND – 5.5
	4 – 9	31	0.23	ND – 1.9
Combined	Cord	67	0.11	ND – 1.1
	0.1 – 3	129	1.1	ND – 10
	4 – 9	59	0.5	ND – 1.8

ND: not detectable.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

In certain cases of infertility, treatment with human gonadotropins poses a potential problem for the accurate measurement of FSH levels. The FSH that is administered can cause the patient to produce antibodies to FSH which will interfere directly with the assay.

Because of pulsatile secretion, samples obtained within the same day from the same patient may fluctuate widely within the reference range, reflecting physiological variation rather than errors in technique or methodology.

Because EDTA would have a significant effect on results, it should not be used as an anticoagulant.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in

combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in mIU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: up to 170 mIU/mL (8.9 ng/mL) WHO 2nd IRP 78/549 [interim replacement code 94/632].

Analytical Sensitivity: 0.1 mIU/mL.

High-dose Hook Effect: None up to 50,000 mIU/mL.

Precision: Samples were assayed in duplicate in 40 runs for a total of 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three FSH solutions (150, 600 and 1,800 mIU/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for FSH. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 600 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 24 volunteers into plastic Becton Dickinson plain serum, heparinized, EDTA and SST[®] vacutainer tubes. All samples were assayed by the IMMULITE 2000 FSH procedure

(Heparin) = 1.09 (Serum) – 0.04 mIU/mL

(SST) = 0.94 (Plain Tubes) + 0.07 mIU/mL

(EDTA) = 0.75 (Serum) + 0.50 mIU/mL

Means:

13.7 mIU/mL (Serum)

14.8 mIU/mL (Heparin)

12.9 mIU/mL (SST)

10.8 mIU/mL (EDTA)

EDTA plasma should not be used as a sample type.

Method Comparison 1: The IMMULITE 2000 FSH mono/mo assay was compared to DPC's IMMULITE 2000 FSH mono/poly assay on 429 patient samples. (Concentration range: approximately 1.0 to 170 mIU/mL. See Graph 1.) By linear regression:

(IML 2000 m/m) =

1.06 (IML 2000 m/p) – 1.29 mIU/mL

r = 0.992

Means:

34.9 mIU/mL (IMMULITE 2000 mono/mo)

34.1 mIU/mL (IMMULITE 2000 mono/poly)

Method Comparison 2: The IMMULITE 2000 FSH mono/mo assay was compared to DPC's IMMULITE FSH mono/poly assay on 289 samples. (Concentration range: approximately 1.0 to 170 mIU/mL. See graph 2.) By linear regression:

(IML 2000 m/m) = 1.04 (IML m/p) –

0.06 mIU/mL

r = 0.991

Means:

36.1 mIU/mL (IMMULITE 2000 mono/mo)

34.9 mIU/mL (IMMULITE mono/poly)

References

- 1) Santner S, Santen R, Kulin H, Demers L. A model for validation of radioimmunoassay kit reagents: measurement of follitropin and lutropin in blood and urine. Clin Chem 1981;27:1892-5.
- 2) Odell W, et al. Radioimmunoassay for luteinizing hormone in human plasma or serum. J Clin Invest 1967;46:248-55.
- 3) Davidsohn I, Henry J, editors. Clinical diagnosis by laboratory methods. 15th ed. Philadelphia: W.B. Saunders, 1974: 704.
- 4) Nankin H, et al. Repetitive luteinizing hormone elevations in serum of normal men. J Clin Endo Metab 1972;33:558-60.
- 5) Beitens I, et al. Gonadotropin determinations in timed 3-hour urine collections during the menstrual cycle and LHRH testing. J Clin Endo Metab 1976;43:46-55.
- 6) Chipman J, et al. Interrelationship of plasma and urinary gonadotropins: correlations for 24 hours, for sleep/wake periods, and for 3 hours after luteinizing hormone-releasing hormone stimulation. Clin Endo Metab 1981;52:225-30.
- 7) Kulin H, et al. Integration of pulsatile gonadotropin secretion by timed urinary measurements: an accurate and sensitive 3-hour test. J Clin Endo Metab 1975;40:783-9.
- 8) Kulin H, Santner S. Timed urinary gonadotropin

measurements in normal infants, children, and adults, and in patients with disorders of sexual maturation. *J Pediatrics* 1977;90:760-5. 9) Urban M, et al. Comparison of estimates of gonadotropin levels by isolated blood samples, integrated blood concentrations, and timed urinary fractions. *J Clin Endo Metab* 1979;48:732-5. 10) Rebar R, Yen S. In: Dorothy Krieger, editor. *Endocrine rhythms*. New York: Raven Press, 1979. 11) Odell WD, et al. Radioimmunoassay for human follicle-stimulating hormone: physiological studies. *J Clin Invest* 1968;47:2551-62. 12) Sutaria UE, Crooke AC. Selection of patients for treatment with human gonadotrophins. *Int J Fertil* 1971;16:42-6.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (mIU/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	6.8	0.20	2.9%	0.28	4.1%
2	22.9	0.79	3.4%	1.11	4.8%
3	44.9	1.44	3.2%	1.86	4.1%
4	68.1	2.89	4.2%	4.27	6.3%
5	103	3.17	3.1%	8.09	7.9%

Linearity (mIU/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	10.0	—	—
	4 in 8	5.08	5.0	102%
	2 in 8	2.46	2.5	98%
	1 in 8	1.19	1.25	95%
2	8 in 8	15.0	—	—
	4 in 8	8.15	7.5	109%
	2 in 8	4.04	3.75	108%
	1 in 8	1.96	1.88	104%
3	8 in 8	29.2	—	—
	4 in 8	16.1	14.6	110%
	2 in 8	7.75	7.30	106%
	1 in 8	4.01	3.65	110%
4	8 in 8	44.1	—	—
	4 in 8	23.5	22.1	106%
	2 in 8	11.9	11.0	108%
	1 in 8	5.97	5.51	108%
5	8 in 8	63.0	—	—
	4 in 8	33.7	31.5	107%
	2 in 8	17.4	15.8	110%
	1 in 8	8.6	7.88	109%
6	8 in 8	109	—	—
	4 in 8	59.5	54.3	110%
	2 in 8	29.1	27.2	107%
	1 in 8	14.6	13.6	107%

Specificity




Compound ¹	Added ²		Apparent FSH ng/mL	% Cross reactivity ³
	ng/mL	ng/mL		
HCG mIU/mL	2,000	215	ND	ND
	200,000	21,538	ND	ND
TSH uIU/mL	10	2.0	ND	ND
	200	40.5	ND	ND
LH mIU/mL	75	7.6	ND	ND
	300	30.5	ND	ND
Prolactin ng/mL	25	25	ND	ND
	200	200	ND	ND
HGH ng/mL	5.0	5.0	ND	ND
	50	50	ND	ND
HPL ng/mL	10,000	10,000	ND	ND
	50,000	50,000	ND	ND

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH