 IMMULITE[®]
2000

LH

DPC[®]

IMMULITE[®] 2000 LH

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of LH in serum, as an aid in clinical diagnosis.

Catalog Number: **L2KLH2** (200 tests), **L2KLH6** (600 tests)

Test Code: **LH** Color: **Red**

CDC Analyte Identifier Code: 3713

CDC Test System Identifier Code: 10418

CLIA Complexity Category: Moderate

Summary and Explanation

Luteinizing hormone (lutropin, LH), a glycoprotein of 28,000 daltons, is secreted by the β -cells of the anterior pituitary under the control of the hypothalamic gonadotropin releasing hormone (GnRH). LH consists of two polypeptide chains, alpha and beta. The alpha chains of LH, FSH, TSH and HCG are biochemically identical, whereas the beta chains are biochemically unique, conferring bioactivity and biological and immunological specificity. In females, LH causes ovulation and steroid (estrogen and progesterone) production by the corpus luteum. Small quantities of LH are also necessary to promote estrogen production by the maturing follicle. In males, it stimulates interstitial (Leydig) cells to produce androgens and estrogens. Circulating levels of LH are controlled by a negative feedback effect on the hypothalamus by the steroid hormones. LH secretion, different for the two sexes and required for normal sexual function, occurs in pulses with rapid fluctuations over the entire reference range. Values for samples obtained in a single day from the same patient may therefore vary widely.

LH measurements are used to define the hypothalamic-pituitary-gonadal axis. Serum gonadotropin determinations permit distinguishing between primary gonadal failure and deficient gonadal stimulation. If LH and FSH levels are elevated, primary gonadal failure is present, whereas if gonadotropin levels

are low, deficient gonadal stimulation has resulted in the hypogonadal state. LH measurement is also clinically important because LH and growth hormone are frequently the first hormones to be affected by pituitary disease.

Serum LH determinations have been very useful in the diagnosis and treatment of infertility in women. A midcycle rise is a good indication that ovulation will occur approximately 24 hours later. Subfertile couples, and women being treated with gonadotropins for infertility, can be informed that ovulation is imminent.

Principle of the Procedure

IMMULITE 2000 LH is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 LH has not been tested with all possible variations of tube types.

Volume Required: 75 μ L serum.

Storage: 2 weeks at 2–8°C, or 2 months at –20°C.¹¹

Since LH is known to exhibit a small circadian rhythm, the time of collection should be noted

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

LH Bead Pack (L2LH12)

With barcode. 200 beads, coated with monoclonal murine anti-LH. Stable at 2–8°C until expiration date.

L2KLH2: 1 pack. **L2KLH6:** 3 packs.

LH Reagent Wedge (L2LHA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-LH in buffer. Stable at 2–8°C until expiration date.

L2KLH2: 1 wedge. **L2KLH6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

LH Adjustors (LLHL, LLHH)

Two vials (Low and High) of lyophilized LH in a nonhuman serum matrix. Reconstitute each vial with **4.0 mL** distilled or deionized

water. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KLH2: 1 set. **L2KLH6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

CON6: Tri-level, multi-constituent control

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 4 weeks.

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of LH.

Expected Values

Based on its relationship to DPC's IMMULITE LH (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

Reference ranges were generated using IMMULITE LH in a multi-national study involving women in apparent good health (age: 16 – 44 years), who volunteered to have blood samples drawn, on a daily basis, throughout one complete ovulatory cycle.

Ovulatory Cycles	n*	LH, mIU/mL	
		Median	Central 95%
Follicular Phase	54 (762)	4.6	1.1 – 11.6
Midcycle	54 (54)	39	17 – 77
Luteal Phase	54 (658)	4.3	ND – 14.7
Perimenstrual, ± 8 days	54 (959)	3.9	ND – 12.0

*Number of subjects (total number of results)

Group	n	LH, mIU/mL	
		Median	Central 95%
Adult Males	135	2.4	0.8 – 7.6
Adult Females:			
Postmenopausal*	75	24.9	11.3 – 39.8
Oral Contraceptives	104	3.1	ND – 8.0

*Preliminary
ND: Not detectable

See Menstrual Cycle Graph (Tables and Graphs section).

A cross-sectional study of pediatric fertility values performed with IMMULITE LH at a "wellness" clinic in the southwestern United States yielded the following results.

Group	Age (yr)	n	LH, mIU/mL	
			Median	Central 95%
Females	Cord	31	ND	
	0.1 – 1.5	46	0.7	ND – 2.3
	1.6 – 9	38	ND	ND – 1.3
Males	Cord	36	ND	ND – 3.6
	0.1 – 1.5	54	1.0	ND – 4.1
	1.6 – 9	46	ND	ND – 3.8
Combined	Cord	67	ND	ND – 3.5
	0.1 – 1.5	100	0.7	ND – 3.7
	1.6 – 9	84	ND	ND – 3.2

ND: Not detectable.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Because of pulsatile secretion, samples obtained within the same day from the same patient may fluctuate widely within the reference range, reflecting physiological variation rather than errors in technique or methodology.

The assay's crossreactivity to HCG is low and will not have impact under normal circumstances. However, for samples with very high levels of HCG, such as pregnancy samples, samples from trophoblastic disease or testicular cancer patients, LH may be falsely elevated due to crossreactivity to HCG.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in mIU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: up to 200 mIU/mL (WHO 1st IRP 68/40 and 2nd IS 80/552)

Analytical Sensitivity: 0.05 mIU/mL.

High-dose Hook Effect:
None up to 85,000 mIU/mL.

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with four LH solutions (200, 400, 1,000 and 2,000 mIU/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for LH. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 188 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to DPC's IMMULITE LH on 105 samples. (Concentration range: approximately 1 to 140 mIU/mL.) By linear regression:

$$(IML\ 2000) = 1.04 (IML) - 0.25\ mIU/mL$$

$$r = 0.972$$

Means:

26 mIU/mL (IMMULITE 2000)
26 mIU/mL (IMMULITE)

References

- 1) Beitens I, et al. Gonadotropin determinations in timed 3-hour urine collections during the menstrual cycle and lhrh testing. *J Clin Endo Metab* 1976;43:46-55.
- 2) Chipman J, et al. Interrelationship of plasma and urinary gonadotropins: correlations for 24 hours, for sleep/wake periods, and for 3 hours after luteinizing hormone-releasing hormone stimulation. *Clin Endo Metab* 1981;52:225-30.
- 3) Davidsohn I, Henry J, editors. *Clinical diagnosis by laboratory methods*. 15th ed. Philadelphia: W.B. Saunders, 1974:704.
- 4) Kulin H, et al. Integration of pulsatile gonadotropin secretion by timed urinary measurements: an accurate and sensitive 3-hour test. *J Clin Endo Metab* 1975;40:783-9.
- 5) Kulin H, Santner S. Timed urinary gonadotropin measurements in normal infants, children, and adults, and in patients with disorders of sexual maturation. *J Pediatrics* 1977;90:760-5.
- 6) Nankin H, et al. Repetitive luteinizing hormone elevations in serum of normal men. *J Clin Endo Metab* 1972;33:558-60.
- 7) Odell W, et al. Radioimmunoassay for luteinizing hormone in human plasma or serum. *J Clin Invest* 1967;46:248.
- 8) Rebar R, Yen S. In: Dorothy Krieger, editor. *Endocrine rhythms*. New York: Raven Press, 1979.
- 9) Santner S, Santen R, Kulin H, Demers L. A model for validation of radioimmunoassay kit reagents: measurement of follitropin and lutropin in blood and urine. *Clin Chem* 1981;27:1892-5.
- 10) Urban M, et al. Comparison of estimates of gonadotropin levels by isolated blood samples, integrated blood concentrations, and timed urinary fractions. *J Clin Endo Metab* 1979;48:732-5.
- 11) Burtis CA,

Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: W.B. Saunders, 1994: 920

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Linearity (mIU/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	1.9	—	—
	4 in 8	0.81	0.95	85%
	2 in 8	0.62	0.48	88%
	1 in 8	0.18	0.24	75%
2	8 in 8	14	—	—
	4 in 8	6.8	7.0	97%
	2 in 8	3.3	3.5	94%
	1 in 8	1.5	1.8	83%
3	8 in 8	24	—	—
	4 in 8	12	12	100%
	2 in 8	6.3	6.0	105%
	1 in 8	2.9	3.0	97%
4	8 in 8	54	—	—
	4 in 8	26	27	96%
	2 in 8	14	14	100%
	1 in 8	6.7	6.8	99%
5	8 in 8	105	—	—
	4 in 8	55	53	104%
	2 in 8	28	26	108%
	1 in 8	12	13	92%
6	8 in 8	154	—	—
	4 in 8	77	77	100%
	2 in 8	37	39	95%
	1 in 8	19	19	100%

Precision (mIU/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	0.15	0.02	13.1	0.036	23.9%
2	0.29	0.015	5.11	0.075	26.3%
3	1.04	0.032	3.04	0.069	6.6%
4	1.89	0.07	3.71	0.118	6.2%
5	8.7	0.31	3.6%	0.58	6.7%
6	20	0.68	3.4%	1.4	7.0%
7	28	1.1	3.9%	1.7	6.1%
8	96	3.7	3.9%	6.4	6.7%
9	170	6.0	3.5%	12	7.1%

Recovery (mIU/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	1.9	—	—
	A	13	12	108%
	B	23	22	105%
	C	55	52	106%
	D	108	102	106%
2	—	14	—	—
	A	23	23	100%
	B	34	33	103%
	C	69	63	110%
	D	128	113	113%
3	—	26	—	—
	A	37	35	106%
	B	51	45	113%
	C	84	75	112%
	D	144	125	115%
4	—	54	—	—
	A	62	61	102%
	B	79	71	111%
	C	105	101	104%
	D	166	151	110%
5	—	105	—	—
	A	108	110	98%
	B	128	120	107%
	C	157	150	105%
	D	220	200	110%

Specificity


Compound ¹	Amount Added ²	Apparent LH ³	% Cross reactivity ⁴
TSH	32 ng/mL (158 µIU/mL)	ND	ND
FSH	18 ng/mL (345 mIU/mL)	ND	ND
HCG	104 ng/mL (970 mIU/mL)	0.213 ng/mL (2.44 mIU/mL)*	0.20%

* The crossreactivity to multiple levels of HCG is shown below.⁵


Compound	mIU/mL Added	Apparent LH mIU/mL	% Cross reactivity ⁴
	230	0.76	0.33%
HCG	970	2.44	0.25%
	2631	6.1	0.23%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH