 IMMULITE<sup>®</sup>  
2000

## **Prolactin**

**DPC<sup>®</sup>**

# IMMULITE® 2000 Prolactin

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of prolactin in serum, as an aid in the diagnosis and treatment of pituitary disorders.

Catalog Number: **L2KPR2** (200 tests),  
**L2KPR6** (600 tests)  
Test Code: **PRL** Color: **Dark Blue**

CDC Analyte Identifier Code: 4915  
CDC Test System Identifier Code: 10418  
CLIA Complexity Category: Moderate

### Summary and Explanation

Human prolactin is a polypeptide hormone of the anterior pituitary with a molecular mass of about 22,800. It plays an essential role in the secretion of milk and has the ability to suppress gonadal function. The very existence of prolactin as a substance distinct from human growth hormone was established only in 1970. Since that time, determination of prolactin has become an important tool in the investigation of amenorrhea, galactorrhea and hypothalamic-pituitary disorders.

As a reference range for circulating prolactin, the literature suggests concentrations up to approximately 20 ng/mL. Values are distinctly elevated at birth but decline to adult levels in less than three months.

Women are reported to have slightly higher mean levels than men, with a slight rise at puberty — apparently estrogen related — and a corresponding fall at menopause. During pregnancy, the prolactin level climbs steadily to ten or twenty times its former value, then drops back down to normal after delivery — within three weeks in nonnursing mothers. In those who breast-feed, the decline to normal is more gradual because of the prompt and dramatic surges in prolactin release induced by suckling. Women taking oral contraceptives or under estrogen treatment may have prolactin levels higher than normal.

In assessing the significance of moderate elevations, it is important to keep in mind that prolactin is a stress hormone. Not only surgery, but events no more distressing than venipuncture or a clinical interview have been reported to occasion a transient rise. Moreover, the release of prolactin is inherently episodic, and day-to-day fluctuations with CVs as high as 30% have been encountered. Finally, there is a sleep-related diurnal variation: prolactin levels increase during sleep and reach their lowest a few hours after waking. The advice sometimes given to draw samples "between nine and noon" is based on the assumption that subjects observe reasonably normal waking hours.

### Principle of the Procedure

IMMULITE 2000 Prolactin is a solid-phase, two-site chemiluminescent immunometric assay.

**Incubation Cycles:** 1 × 30 minutes.

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Prolactin has not been tested with all possible variations of tube types.

**Volume Required:** 25 µL serum.

**Storage:** 7 days at 2–8°C, or 3 months at –20°C.

## Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

### Prolactin Bead Pack (L2PR12)

With barcode. 200 beads, coated with monoclonal murine anti-prolactin. Stable at 2–8°C until expiration date.

**L2KPR2:** 1 pack. **L2KPR6:** 3 packs.

### Prolactin Reagent Wedge (L2PRA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-prolactin, with preservative. Stable at 2–8°C until expiration date.

**L2KPR2:** 1 wedge. **L2KPR6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### Prolactin Adjustors (LPRL, LPRH)

Two vials (Low and High), of lyophilized human prolactin in a serum/buffer matrix. Reconstitute each vial with **2.0 mL** distilled or deionized water, and mix by gentle

inversion. Stable after reconstitution for 60 days (aliquotted) at –20°C.

**L2KPR2:** 1 set. **L2KPR6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

## Kit Components Supplied Separately

### Prolactin Sample Diluent (L2PRZ)

For on-board dilution of high samples. 25 mL of concentrated (ready-to-use) processed, prolactin-free serum/buffer matrix. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2PRZ:** 3 labels

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**L2ZT:** 250 Sample Diluent Test Tubes (16 × 100 mm)

**L2ZC:** 250 Sample Diluent Tube Caps

**CON6:** Tri-level, multi-constituent control

Also Required

Distilled or deionized water; test tubes; controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

**Recommended Adjustment Interval:**

4 weeks.

**Quality Control Samples:** Use controls or serum pools with at least two levels (low and high) of prolactin.

## Expected Values

Based on its relationship to DPC's IMMULITE Prolactin (see Method

Comparison), the assay can be expected to have essentially the same reference range.

Group	n	Median	96% Range
Adult Males	19	6.2 ng/mL 131 mIU/L	2.5 – 17 53 – 360

A study<sup>17</sup> performed with IMMULITE Prolactin yielded the following results.

Group	n	Median	95% Range
Adult Females	115	9.4 ng/mL 199 mIU/L	1.9 – 25 40 – 530

As summarized in DPC's technical report ZB157,<sup>17</sup> in a study which followed normal ovulating women on a daily basis throughout one complete cycle, somewhat higher values were obtained – an outcome consistent with the increasingly stressful sample collection process entailed by the study design.

A cross-sectional study of pediatric fertility values performed with IMMULITE Prolactin at a "wellness" clinic in the southwestern United States yielded the following results.

Group	Age (yr)	n	Prolactin, ng/mL	
			Median	Central 95%
Females	Cord	28	380	200 – 675
	0.1 – 0.5	28	15	1 – 140
	0.6 – 9	55	11	2 – 43
Males	Cord	27	295	150 – 565
	0.1 – 0.5	36	19	4 – 65
	0.6 – 9	55	8	0.6 – 29
Combined	Cord	55	340	160 – 665
	0.1 – 0.5	64	117	2 – 125
	0.6 – 9	110	9	1 – 40

Group	Age (yr)	n	Prolactin, mIU/L	
			Median	Central 95%
Females	Cord	28	8,056	4,240 – 14,310
	0.1 – 0.5	28	318	21 – 2,968
	0.6 – 9	55	233	42 – 912
Males	Cord	27	6,254	3,180 – 11,978
	0.1 – 0.5	36	403	85 – 1,378
	0.6 – 9	55	170	13 – 615
Combined	Cord	55	7,208	3,392 – 14,098
	0.1 – 0.5	64	2,480	42 – 2,650
	0.6 – 9	110	191	21 – 848

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

### Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

### Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

#### Conversion Factor:

ng/mL × 21.2 → mIU/L [3rd IS 84/500]

**Calibration Range:** Up to 150 ng/mL (Up to 3,180 mIU/L).

**Analytical Sensitivity:** 0.16 ng/mL (3.4 mIU/L).

**High-dose Hook Effect:** None up to 20,500 ng/mL (434,600 mIU/L).

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with four prolactin solutions (190, 356, 928 and 1,740 ng/mL), were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for prolactin. (See "Specificity" table.)

**Alternate Sample Type:** A *limited* study of matched samples from laboratory volunteers yielded comparable results for serum and heparinized plasma.

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Presence of hemoglobin in concentrations up to 375 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

**Method Comparison:** The assay was compared to DPC's IMMULITE Prolactin on 148 samples (Concentration range: approximately 1 to 140 ng/mL. See graph.) By linear regression:

$(\text{IML } 2000) = 1.00 (\text{IML}) - 0.96 \text{ ng/mL}$   
 $r = 0.988$

Means:  
44 ng/mL (IMMULITE 2000)  
45 ng/mL (IMMULITE)

## References

- 1) Cowden E, et al. Laboratory assessment of prolactin status. *Ann Clin Biochem* 1979;16:113–21.
- 2) Cowden E, et al. Tests of prolactin secretion. *Lancet* 1979;1155–8.
- 3) Djursing H. Short- and long-term fluctuations in plasma prolactin concentration in normal subjects. *Acta Endocrinol* 1981; 97:1–6.
- 4) Greer ME, et al. Prevalence of hyperprolactinemia in anovulatory women. *Obstet Gynecol* 1980;56:65–9.
- 5) Harrington RA, et al. Metoclopramide: an updated review. *Drugs* 1983;25:451–94.
- 6) Healy DL, et al.

Pituitary autonomy in hyperprolactinemia secondary amenorrhea: results of hypothalamic-pituitary testing. *J Clin Endocrinol Metab* 1977;44:809–18.

- 7) Kleinberg DL, Noel GL, Frantz AG. Galactorrhea: a study of 235 cases. *N Engl J Med* 1977;296:589–600.
- 8) Knuth UA, Friesen HG. Prolactin and pregnancy. *Curr Top Exp Endocrinol* 1983;4:69–96.
- 9) Maxson WS, Hammond CB. Hyperprolactinemia. *Contemp Ob/Gyn* 1982 Jan;19(1):49–63, 19:67–83.
- 10) Pepperell RJ. Prolactin and reproduction. *Fertil Steril* 1981;35:267–74.
- 11) Babson, AL. The IMMULITE Automated Immunoassay System. *J Clin Immunoassay* 1991;14:83–8.
- 12) Thorner MO. Prolactin: clinical physiology and management of hyperprolactinemia. In: Martini L, Besser GM, et al, editors. *Clinical endocrinology*. New York: Academic Press, 1977;319–61.
- 13) Tolis G, Franks S. Physiology and pathology of prolactin secretion. In: Tolis G, et al, editors. *Clinical endocrinology: a pathophysiological approach*. New York: Raven Press, 1979;291–317.
- 14) Tolis G, et al. Aspects of prolactin pathophysiology. *Periodicum Biologorum* 1983;85 Suppl 1:29–36.
- 15) Yen SSC. Neuroendocrine regulation of gonadotropin and prolactin secretion in women: disorders in reproduction. In: Vaitukaitis J, editor. *Clinical reproductive endocrinology*. New York: Elsevier Biomedical, 1982;137–76.
- 16) Zacur HA. Use of the human prolactin immunoassay. *J Clin Immunoassay* 1983 Spring;6(1):63–67.
- 17) Vankrieken L. IMMULITE reproductive hormone assays: multicenter reference range data. Los Angeles: Diagnostic Products Corporation, 2000. Document No. ZB157-D.

## Technical Assistance

In the United States, contact DPC's Technical Services department.  
Tel: 800.372.1782 or 973.927.2828  
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

## Tables and Graphs

### Precision (ng/mL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	4.9	0.11	2.2%	0.34	6.9%
2	8.8	0.25	2.8%	0.72	8.2%
3	19	0.68	3.6%	1.4	7.4%
4	22	0.51	2.3%	1.3	5.9%
5	48	1.2	2.5%	3.3	6.9%
6	164	3.8	2.3%	13	7.9%

### Linearity (ng/mL)

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	34	—	—
	4 in 8	17	17	100%
	2 in 8	8.2	8.5	96%
	1 in 8	4.0	4.3	93%
2	8 in 8	73	—	—
	4 in 8	36	37	97%
	2 in 8	17	18	94%
	1 in 8	8.1	9.1	89%
3	8 in 8	84	—	—
	4 in 8	41	42	98%
	2 in 8	20	21	95%
	1 in 8	9.6	11	87%
4	8 in 8	120	—	—
	4 in 8	59	60	98%
	2 in 8	28	30	93%
	1 in 8	13	15	87%

### Specificity


Compound <sup>1</sup>	Amount Added <sup>2</sup> (ng/mL)	% Cross-reactivity <sup>3</sup>
FSH	1,000	ND
	200	ND
HCG	1,000	ND
	100	ND
LH	1,000	ND
	500	ND
TSH	1,000	0.14%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
**EURO/DPC** LTD
EC REP DPC Biermann GmbH